

Package ‘AffyCompatible’

October 18, 2022

Type Package

Title Affymetrix GeneChip software compatibility

Version 1.56.0

Author Martin Morgan, Robert Gentleman

Maintainer Martin Morgan <martin.morgan@roswellpark.org>

Depends R (>= 2.7.0), XML (>= 2.8-1), RCurl (>= 0.8-1), methods

Imports Biostrings

Description This package provides an interface to Affymetrix chip annotation and sample attribute files. The package allows an easy way for users to download and manage local data bases of Affymetrix NetAffx annotation files. The package also provides access to GeneChip Operating System (GCOS) and GeneChip Command Console (AGCC)-compatible sample annotation files.

License Artistic-2.0

biocViews Infrastructure, Microarray, OneChannel

git_url <https://git.bioconductor.org/packages/AffyCompatible>

git_branch RELEASE_3_15

git_last_commit 37ea4bb

git_last_commit_date 2022-04-26

Date/Publication 2022-10-18

R topics documented:

| | |
|----------------------------------|---|
| AffyCompatible-package | 2 |
| accessors | 3 |
| AffyCompatible-class | 3 |
| Internal-class | 4 |
| NetAffxResource | 5 |
| NetAffxResource-class | 6 |
| readArr | 8 |
| readMage | 9 |

| | |
|--|-----------|
| readNetAffx | 10 |
| readXml | 11 |
| xclass | 12 |
| xmlValue.XMLAttributeValue | 13 |
| zzz-ARRCompatibility-class | 13 |
| zzz-DTTCompatibility-class | 14 |
| zzz-NetAffxCompatibility-class | 16 |
| Index | 18 |

AffyCompatible-package

Affymetrix chip annotation and sample attribute management

Description

This package manages Affymetrix NetAffx annotation files, and parses MAGE and ARR sample attribute files into R data structures.

NetAffx provides a way to manage the diverse annotation files provided by Affymetrix for their arrays.

ARR files provide attribute information about samples used in an experiment. These files are created, for instance, by the GeneChip Command Console or compliant applications, or

MAGE files also provide information about samples. Support in this package is for MAGE produced by Affymetrix GCOS compliant applications, and emphasizes sample description and processing.

Details

Package details are available in the DESCRIPTION file of the package, or with the R command `packageDescription("AffyCompatible")`

Overviews of how to use the functionality of this package are available in the package vignettes. See `vignette(package="AffyCompatible")`.

Author(s)

Written by Martin Morgan.

Maintainer: Martin Morgan <martin.t.morgan@gmail.com>

accessors

Access NetAffxResource-class content

Description

These functions exist to access the information represented in objects of `linkS4class{NetAffxResource}` and related classes.

Usage

```
annotationsFile(x)
directory(x)
password(x)
user(x)
```

Arguments

x An object for which the accessor is defined.

Value

See, e.g., `help('directory,NetAffxResource-method')` for values returned by specific methods.

Author(s)

Martin Morgan

See Also

[NetAffxResource](#)

AffyCompatible-class

Virtual classes defining the "AffyCompatibility" hierarchy

Description

The classes `AffyCompatible`, `DTTCompatibility`, `ARRCompatibility`, and `NetAffxCompatibility` define the overall class hierarchy in this package. Information about these classes is not usually relevant to the user.

Objects from the Class

These are virtual classes: no objects may be created from them.

Methods

show signature(object = "AffyCompatible"): for each slot, provide the slot name, and then invoke show for the slot content.

Author(s)

Martin Morgan

Examples

```
showClass(getClass("NetAffxCompatibility",
  where=getNamespace("AffyCompatible")))
```

Internal-class

(Internal) Represent 'typed associations' from XML documents

Description

This internal class represents 'typed associates' (vectors of identically-typed objects) that appear as vectors to the R user.

Objects from the Class

Objects are created internally, during XML document parsing.

Slots

.Data: Object of class "list" containing the association objects
.extends: Object of class "character" and length 1 describing the class each element in .Data must extend

Extends

Class "[list](#)", from data part. Class "[vector](#)", by class "list", distance 2.

Methods

[signature(x = ".TypedAssociation", i = "numeric", j = "missing"): Subset the typed association, returning a typed association with elements i.

[signature(x = ".TypedAssociation", i = "ANY", j = "ANY"): This method is present to stop access to subsetting methods for class list; it has no user-visible functionality.

show signature(object = ".TypedAssociation"): Compactly display the typed associate, e.g., for .extnds=Foo and length 5, display Foo(5).

Author(s)

Martin Morgan

NetAffxResource *Create and access a NetAffxResource database*

Description

Use this function to create and access a NetAffxResource database.

Usage

```
NetAffxResource(user = character(0),  
                password = character(0),  
                affxLicence = "FHCRC0607",  
                directory = tempdir(), ...)
```

Arguments

| | |
|-------------|---|
| user | Length 1 character vector containing the Affymetrix user name authorized to access the NetAffx resource. |
| password | Length 1 character vector containing the Affymetrix password of the user. |
| affxLicence | Length 1 character vector containing the Affymetrix licence (typically application specific, i.e., the specified licence covers AffyCompatible). |
| directory | (Optional) length 1 character vector providing the full path to the local directory where the NetAffx database and files will be managed. Defaults to a session-specific temporary directory. |
| ... | Additional arguments, not currently used. |

Details

NetAffxResource creates a [NetAffxResource](#) object that manages a data base of NetAffx resources. The data base is located on the user file system as indicated by the `directory` argument. It is populated when the user requests specific files.

Users should provide a `directory` argument, so that the data base can be reused across sessions. This means that large files are only downloaded once.

Additional details are provided in the vignette accompanying this package.

Value

[NetAffxResource](#).

Author(s)

Martin Morgan

See Also

[NetAffxResource](#). The vignette ‘Annotations with NetAffx’.

NetAffxResource-class *Manage NetAffx annotations*

Description

This class is for accessing and retrieving NetAffx annotation files for Affymetrix chips. It requires that the user has established an account with Affymetrix. Detailed usage is described in the 'Annotations with NetAffx' vignette.

Objects from the Class

Objects are created by calls of the form `NetAffxResource(user="user@research.org", password="terces")`, as illustrated below and in the 'Annotations with NetAffx vignette'.

Slots

user: Object of class "character" containing the full user name (typically an email address) obtained when the user registers with Affymetrix.

password: Object of class "character" containing the password associated with the Affymetrix user name. The password is entered and stored in plain text, and should not be considered secure within Bioconductor.

affxLicence: Object of class "character" containing the licence, obtained from Affymetrix, for retrieving NetAffx annotations. The licence is usually software specific (e.g., for the Affy-Compatible package of Bioconductor). Users will not usually change the default value.

directory: Object of class "character", containing the local directory path where the NetAffx data base and annotation files will reside.

annotationsFile: Object of class "character" containing the file name of the NetAffx data base. Users will not usually change the default value.

affxUrl: Object of class "character" containing the URL used to download the NetAffx annotations file. Users will not usually change the default value.

Extends

Class "[NetAffxCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "NetAffxCompatibility", distance 2.

Methods

[[signature(x = "NetAffxResource", i = "numeric", j = "missing")

[[signature(x = "NetAffxResource", i = "character", j = "missing"): Retrieve an object representing all NetAffx annotations available for an array. The argument *i* may be a numeric index or character string corresponding to an array name returned by `names(x)`.

[[signature(x = "NetAffxResource", i = "numeric", j = "numeric")

[[signature(x = "NetAffxResource", i = "character", j = "numeric")

[[signature(x = "NetAffxResource", i = "numeric", j = "character")

- [[signature(x = "NetAffxResource", i = "character", j = "character")**: Retrieve an object representing a specific annotation for an array. The argument *i* may be a numeric index or character string corresponding to an array name returned by `names(x)`. The argument *j* may be a numeric index or character string corresponding to an annotation returned by `description(x[[i]])`.
- affxLicence** signature(x = "NetAffxResource"): the character string representation of the licence used to access the NetAffx site.
- affxUrl** signature(x = "NetAffxResource"): the url used to access the NetAffx site.
- annotationsFile** signature(x = "NetAffxResource"): the file name of NetAffx data base, as stored on the local disk of the user.
- directory** signature(x = "NetAffxResource"): the directory in which the NetAffx data base and annotation files are stored on the local disk of the user.
- names** signature(x = "NetAffxResource"): the names of the arrays for which annotations are available.
- password** signature(x = "NetAffxResource"): The password associated with the user name obtained when the user registered with Affymetrix. The password is entered and stored as plain text, and is not secure.
- readAnnotation** signature(netAffxResource = "NetAffxResource", array = "numeric", annotation = "numeric", ...)
- readAnnotation** signature(netAffxResource = "NetAffxResource", array = "numeric", annotation = "character", ...)
- readAnnotation** signature(netAffxResource = "NetAffxResource", array = "character", annotation = "numeric", ...)
- readAnnotation** signature(netAffxResource = "NetAffxResource", array = "character", annotation = "character", ...)
- readAnnotation** signature(netAffxResource = "NetAffxResource", array = "missing", annotation = "AffxAnnotation", ...): Read a specific annotation file as an R object. For the first four methods, array, if present, indexes or names the array while annotation indexes or names the annotation. All methods also accept `content = "logical"` to indicate whether the content of the file should be returned if possible (`content=TRUE`, default) or only the path to the local file (`content=FALSE`) and `update = "logical"` indicating whether the file should be read from disk if available (`update=FALSE`, default) or retrieved from the NetAffx web site (`update=TRUE`). `readAnnotation` parses CSV, Tabular, and PSI files to data frames, and reads FASTA files to a list of description / sequence pairs (see `Biostrings::readFASTA`).
The ... argument is passed to corresponding read functions. CSV files are parsed using `read.csv(conn, ...)`, Tabular files with `read.delim(conn, ...)`, PSI files with `read.delim(conn, header=FALSE, skip=1, sep="\t", ...)`, and FASTA with `readFASTA(conn, ...)`.
- show** signature(object = "NetAffxResource"): compactly show the NetAffxResource object.
- user** signature(x = "NetAffxResource"): The user name (typically email address) obtained when the user registers with the Affymetrix web site.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette ‘Annotations with NetAffx’. Affymetrix user login site <https://www.affymetrix.com/site/login/login.affx>.

[.Rprofile](#) for ways to automatically set variables (such as password, and directory) for each session.

Examples

```
## password must be obtained by the user
rsrc <- NetAffxResource(user="mtmorgan@fhcrc.org", password=password)
head(names(rsrc))
affxDescription(rsrc[["Bovine"]])
## paradigm to display all elements of 'annotation'
sapply(affxAnnotation(rsrc[["Bovine"]]), force)
## obtain and read a single annotation
anno <- rsrc[["Porcine", "Annotations, CSV format"]]
df <- readAnnotation(rsrc, annotation=anno)
```

readArr

Read ARR sample attribute files

Description

This functions reads the content of Affymetrix ARR sample attribute files created by GeneChip(R) Command Console. Objects are represented as R classes. Navigation and extraction of data are described in the vignette ‘Retrieving MAGE and ARR sample attributes’.

Usage

```
readArr(fls, ...)
```

Arguments

| | |
|-----|--|
| fls | A character vector containing the paths to files to be read. |
| ... | Additional arguments, currently unused. |

Value

An object or list of objects of class `ArrayFileSet`, depending on whether one or several file paths are provided.

Author(s)

Martin Morgan

See Also

The vignette ‘Retrieving MAGE and ARR sample attributes’ for detailed illustration of use. [xclass](#) and [readXml](#) for advanced access methods.

Examples

```
dir <- system.file("extdata", "ARR", package="AffyCompatible")
fls <- list.files(dir, full=TRUE)
readArr(fls[[1]])
readArr(fls)
```

readMage

Read DTT (MAGE) sample attribute files

Description

This functions reads the content of Affymetrix DTT (MAGE) sample attribute files created by GCOS. Objects are represented as R classes. Navigation and extraction of data are described in the vignette ‘Retrieving MAGE and ARR sample attributes’.

Usage

```
readMage(fls, ...)
```

Arguments

`fls` A character vector containing the paths to files to be read.
`...` Additional arguments, currently unused.

Value

An object or list of objects of class `ArrayFileSet`, depending on whether one or several file paths are provided.

Author(s)

Martin Morgan

See Also

The vignette ‘Retrieving MAGE and ARR sample attributes’ for detailed illustration of use. [xclass](#) and [readXml](#) for advanced access methods.

Examples

```
dir <- system.file("extdata", "DTT", package="AffyCompatible")
fls <- list.files(dir, full=TRUE)
basename(fls)
readArr(fls[[1]])
```

| | |
|-------------|---|
| readNetAffx | <i>Update or read a NetAffx data base</i> |
|-------------|---|

Description

This function provides a way to read or update a NetAffx data base. The usual use for this is to force the data base description to be updated (by accessing the Affymetrix web site) or to retrieve the entire list of NetAffx resource data base. A more normal use case is to

Usage

```
readNetAffx(netAffxResource, update = FALSE, ...)
```

Arguments

| | |
|-----------------|--|
| netAffxResource | An object of class NetAffxResource containing the path to the local data base, and permissions for access the Affymetrix NetAffx resource. |
| update | An object of logical(1) indicating whether the data base is to be updated (TRUE) or not. |
| ... | Additional arguments, not currently used. |

Value

An object representing the annotations found; see [NetAffxResource](#) for additional details.

Author(s)

Martin Morgan

Examples

```
## retrieve all annotations
## Not run:
example("NetAffxResource-class")
readNetAffx(rsrc)

## End(Not run)
```

| | |
|---------|--|
| readXml | <i>(Advanced) Read XML files of Affymetrix sample attributes or annotation</i> |
|---------|--|

Description

This functions reads the content of Affymetrix DTT (MAGE) or ARR sample files, or AffxNetAffx-AnnotFile array annotation files, return an internal representation of the underlying XML. The return value can be used with [xclass](#) to rapidly create R objects representing deeply nested nodes of the document. The return value can also be used with `xpathApply` in the XML package for access to character-based representation of the data.

Usage

```
readXml(f1, ...)
```

Arguments

| | |
|------------------|--|
| <code>f1</code> | A length-1 character vector giving the full path to the file to be parsed. |
| <code>...</code> | Additional arguments, passed to <code>xmlTreeParse</code> . |

Value

An object of class `XMLInternalDocument`.

Author(s)

Martin Morgan

See Also

The vignettes ‘Retrieving MAGE and ARR sample attributes’ and ‘Annotation retrieval with NetAffxResource’ for detailed illustration of use.

[xclass](#) for instantiating objects from xpath queries.

Examples

```
mageDir <- system.file("extdata", "DTT", package="AffyCompatible")
mageFiles <- list.files(mageDir, full=TRUE)
xml <- readXml(mageFiles[[1]])
xclass(xml, "//MeasuredBioAssay")[[1]]
```

`xclass`*(Advanced) Construct R objects from XML document queries*

Description

This function takes an xml object, created from [readXml](#), and an xpath query referencing an element of the document, and returns an R object representing the element.

The xpath query must retrieve an element, not an attribute.

Usage

```
xclass(doc, xpathq, prefix = character(0))
```

Arguments

| | |
|---------------------|--|
| <code>doc</code> | An XMLInternalDocument retrieved using readXml . |
| <code>xpathq</code> | A length-1 character string representing an xpath query into the document. |
| <code>prefix</code> | A prefix used to map, when necessary, between the XML class name and the corresponding R class. For NetAffx-based classes, the prefix is 'Affx'; no prefix is required for DTT or ARR documents. |

Value

An R object corresponding to the class(es) referenced by the xpath query.

Author(s)

Martin Morgan

References

<http://www.w3.org/TR/xpath> provides a very useful description of xpath. Additional references are in the vignette 'Retrieving MAGE and ARR sample attributes'.

See Also

The vignettes 'Retrieving MAGE and ARR sample attributes' and 'Annotation retrieval with NetAffxResource' for detailed illustration of use.

[readXml](#) for appropriately reading XML documents.

Examples

```
example(readXml)
```

xmlValue.XMLAttributeValue

Additional functions defined in AffyCompatible

Description

This page documents additional, miscellaneous, functions defined by AffyCompatible

Usage

```
## S3 method for class 'XMLAttributeValue'
xmlValue(x, ignoreComments = FALSE, recursive = TRUE,
         encoding=XML:::CE_NATIVE, trim=FALSE)
```

Arguments

| | |
|----------------|---|
| x | An object of class XMLAttributeValue. |
| ignoreComments | Logical; ignored by the xmlValue.XMLAttributeValue method |
| recursive | ignored |
| encoding | File encoding |
| trim | Remove trailing whitespace? Ignored by this method |

Value

A character vector (usually length 1) containing the character representation of the XML attribute.

zzz-ARRCompatibility-class

(Auto-generated) classes and methods for navigating ARR sample attributes.

Description

These classes and methods facilitate navigation of the Affymetrix ARR sample attribute files. ARR files are generated by the Affymetrix AGCC (GeneChip(R) Command Console) application. Specifically, the classes are R representations of the Affymetrix classes defined in ArraySetAndTemplate-File.dtd. The 'show' methods for these classes indicate how objects are to be accessed; details are in the 'Retrieving MAGe and ARR sample attributes' vignette.

Objects from the Class

Objects are created internally, by parsing XML retrieved from ARR, and are not created directly by the user.

Slots

Objects have a variety of slots. To discover slot names and content type, see the `slotNames` portion of the example, below.

Extends

Class "[ARRCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "ARRCompatibility", distance 2.

Methods

'Accessor' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form `accessor: content`. To retrieve the content of an instance object, invoke `accessor(object)`.

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using `[, [[]`. A useful paradigm for displaying all instances is illustrated in the vignette.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Retrieving MAGE and ARR sample attributes'.

Examples

```
## defined and auto-generated classes
names(slot(getClass("ARRCompatibility",
                  where=getNamespace("AffyCompatible")),
       "subclasses"))
```

zzz-*DTTCompatibility-class*

(Auto-generated) classes and methods for navigating DTT sample attributes.

Description

These classes and methods facilitate navigation of the Affymetrix DTT sample attribute files. DTT files are generated by the Affymetrix GCOS application. Specifically, the classes are R representations of the Affymetrix classes defined in MAGE-ML.dtd. The 'show' methods for these class indicate how objects are to be accessed; details are in the 'Retrieving MAGE and ARR sample attributes' vignette.

Objects from the Class

Objects are created internally, by parsing XML retrieved from ARR, and are not created directly by the user.

Slots

Objects have a variety of slots. To discover slot names and content type, see the `slotNames` portion of the example, below.

Extends

Class "[DTTCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "[DTTCompatibility](#)", distance 2.

Methods

'Accessor' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form `accessor: content`. To retrieve the content of an instance object, invoke `accessor(object)`.

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using `[, []`. A useful paradigm for displaying all instances is illustrated in the vignette.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Retrieving MAGe and ARR sample attributes'.

Examples

```
## defined and auto-generated classes
names(slot(getClass("DTTCompatibility",
                  where=getNamespace("AffyCompatible")),
       "subclasses"))
```

zzz-NetAffxCompatibility-class

(Auto-generated) classes and methods for navigating NetAffx resources

Description

These classes and methods facilitate navigation of the NetAffx resource data base. They are R representations of the Affymetrix classes defined in NetAffxAnnotFileList.dtd. The 'show' methods for these classes indicate how objects are to be accessed; details are in the 'Annotations with NetAffx' vignette.

The auto-generated classes are discoverable as indicated in the example below.

Objects from the Class

Objects are created internally, by parsing XML retrieved from the local NetAffx data base or from the Affymetrix web site.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "[NetAffxCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "[NetAffxCompatibility](#)", distance 2.

Methods

'Read' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke accessor(object).

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using [, [[]. A useful paradigm for displaying all instances is illustrated on the [NetAffxResource](#) page.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Annotations with NetAffx'. Affymetrix user login site <https://www.affymetrix.com/site/login/login.affx>.

[NetAffxResource](#) for creating and managing the NetAffx data base.

[NetAffxResource](#) for navigating the data base.

Examples

```
## defined and auto-generated classes
getClass("NetAffxCompatibility", where=getNamespace("AffyCompatible"))

## slots of an auto-generated class
getSlots(getClass("AffxNetAffxAnnotFileList",
  where=getNamespace("AffyCompatible")))

## Not run:
example("NetAffxrResource-class")

## End(Not run)
```

Index

- * **classes**
 - AffyCompatible-class, 3
 - Internal-class, 4
 - NetAffxResource-class, 6
 - readArr, 8
 - readMage, 9
 - zzz-ARRCompatibility-class, 13
 - zzz-DTTCompatibility-class, 14
 - zzz-NetAffxCompatibility-class, 16
- * **database**
 - NetAffxResource, 5
 - readNetAffx, 10
- * **manip**
 - readXml, 11
 - xclass, 12
 - xmlValue.XMLAttributeValue, 13
- * **methods**
 - accessors, 3
- * **package**
 - AffyCompatible-package, 2
 - .Rprofile, 8
 - .TypedAssociation, 14–16
 - .TypedAssociation-class
 - (Internal-class), 4
 - [, .TypedAssociation, ANY, ANY-method
 - (Internal-class), 4
 - [, .TypedAssociation, numeric, missing, ANY-method
 - (Internal-class), 4
 - [, .TypedAssociation, numeric, missing-method
 - (Internal-class), 4
 - [[, NetAffxResource, character, character-method
 - (NetAffxResource-class), 6
 - [[, NetAffxResource, character, missing-method
 - (NetAffxResource-class), 6
 - [[, NetAffxResource, character, numeric-method
 - (NetAffxResource-class), 6
 - [[, NetAffxResource, numeric, character-method
 - (NetAffxResource-class), 6
 - [[, NetAffxResource, numeric, missing-method
 - (NetAffxResource-class), 6
 - [[, NetAffxResource, numeric, numeric-method
 - (NetAffxResource-class), 6
 - accession (zzz-DTTCompatibility-class), 14
 - accession, DatabaseEntry-method
 - (zzz-DTTCompatibility-class), 14
 - accessions_assnlist
 - (zzz-DTTCompatibility-class), 14
 - accessions_assnlist, BibliographicReference-method
 - (zzz-DTTCompatibility-class), 14
 - Accessions_assnlist-class
 - (zzz-DTTCompatibility-class), 14
 - accessionVersion
 - (zzz-DTTCompatibility-class), 14
 - accessionVersion, DatabaseEntry-method
 - (zzz-DTTCompatibility-class), 14
 - accessors, 3
 - acpassword (NetAffxResource), 5
 - action (zzz-DTTCompatibility-class), 14
 - action, Audit-method
 - (zzz-DTTCompatibility-class), 14
 - action_assn
 - (zzz-DTTCompatibility-class), 14
 - action_assn, Treatment-method
 - (zzz-DTTCompatibility-class), 14
 - Action_assn-class
 - (zzz-DTTCompatibility-class), 14
 - actionMeasurement_assn

- (zzz-DTTCompatibility-class),
14
- actionMeasurement_assn,Treatment-method
(zzz-DTTCompatibility-class),
14
- ActionMeasurement_assn-class
(zzz-DTTCompatibility-class),
14
- activityDate
(zzz-DTTCompatibility-class),
14
- activityDate,ProtocolApplication-method
(zzz-DTTCompatibility-class),
14
- address(zzz-DTTCompatibility-class), 14
- address,Organization-method
(zzz-DTTCompatibility-class),
14
- address,Person-method
(zzz-DTTCompatibility-class),
14
- adjustments_assnlist
(zzz-DTTCompatibility-class),
14
- adjustments_assnlist,ArrayManufactureDeviationFileRecord
(zzz-DTTCompatibility-class),
14
- Adjustments_assnlist-class
(zzz-DTTCompatibility-class),
14
- affiliation_assnref
(zzz-DTTCompatibility-class),
14
- affiliation_assnref,Person-method
(zzz-DTTCompatibility-class),
14
- Affiliation_assnref-class
(zzz-DTTCompatibility-class),
14
- affxAnnotation
(zzz-NetAffxCompatibility-class),
16
- affxAnnotation,AffxArray-method
(zzz-NetAffxCompatibility-class),
16
- AffxAnnotation-class
(zzz-NetAffxCompatibility-class),
16
- affxArray
(zzz-NetAffxCompatibility-class),
16
- affxArray,AffxNetAffxAnnotFileList-method
(zzz-NetAffxCompatibility-class),
16
- affxArray,Array_assnlist-method
(zzz-DTTCompatibility-class),
14
- AffxArray-class
(zzz-NetAffxCompatibility-class),
16
- affxCatgory
(zzz-DTTCompatibility-class),
14
- affxCatgory,OntologyEntry-method
(zzz-DTTCompatibility-class),
14
- affxComment
(zzz-NetAffxCompatibility-class),
16
- affxComment,AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
16
- affxComment,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- AffxComment-class
(zzz-NetAffxCompatibility-class),
16
- affxCompression
(zzz-NetAffxCompatibility-class),
16
- affxCompression,AffxUrl-method
(zzz-NetAffxCompatibility-class),
16
- affxCrc32
(zzz-NetAffxCompatibility-class),
16
- affxCrc32,AffxFile-method
(zzz-NetAffxCompatibility-class),
16
- affxCrc32,AffxUrl-method
(zzz-NetAffxCompatibility-class),
16
- affxDate
(zzz-NetAffxCompatibility-class),
16

- affxDate, AffxFile-method
(zzz-NetAffxCompatibility-class),
16
- affxDate, Audit-method
(zzz-DTTCompatibility-class),
14
- affxDescription
(zzz-NetAffxCompatibility-class),
16
- affxDescription, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
16
- affxDescription, AffxArray-method
(zzz-NetAffxCompatibility-class),
16
- affxEnd (zzz-DTTCompatibility-class), 14
- affxEnd, CompositePosition-method
(zzz-DTTCompatibility-class),
14
- affxEnd, ReporterPosition-method
(zzz-DTTCompatibility-class),
14
- affxEnd, SequencePosition-method
(zzz-DTTCompatibility-class),
14
- affxFile
(zzz-NetAffxCompatibility-class),
16
- affxFile, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
16
- AffxFile-class
(zzz-NetAffxCompatibility-class),
16
- affxImage (zzz-DTTCompatibility-class),
14
- affxImage, PhysicalBioAssayData_assnlist-method
(zzz-DTTCompatibility-class),
14
- affxLicence
(zzz-NetAffxCompatibility-class),
16
- affxLicence, NetAffxResource-method
(NetAffxResource-class), 6
- affxName
(zzz-NetAffxCompatibility-class),
16
- affxName, AffxArray-method
(zzz-NetAffxCompatibility-class),
16
- affxName, AffxFile-method
(zzz-NetAffxCompatibility-class),
16
- AffxNetAffxAnnotFileList-class
(zzz-NetAffxCompatibility-class),
16
- affxOrder (zzz-DTTCompatibility-class),
14
- affxOrder, BioDataCube-method
(zzz-DTTCompatibility-class),
14
- affxOrder, Treatment-method
(zzz-DTTCompatibility-class),
14
- affxPerson
(zzz-DTTCompatibility-class),
14
- affxPerson, Contact_assnlist-method
(zzz-DTTCompatibility-class),
14
- affxPerson, Person_assnlist-method
(zzz-DTTCompatibility-class),
14
- affxRow (zzz-DTTCompatibility-class), 14
- affxRow, FeatureLocation-method
(zzz-DTTCompatibility-class),
14
- affxRow, Zone-method
(zzz-DTTCompatibility-class),
14
- affxSequence
(zzz-DTTCompatibility-class),
14
- affxSequence, BioSequence-method
(zzz-DTTCompatibility-class),
14
- affxSize
(zzz-NetAffxCompatibility-class),
16
- affxSize, AffxFile-method
(zzz-NetAffxCompatibility-class),
16
- affxSize, AffxUrl-method
(zzz-NetAffxCompatibility-class),
16
- affxStart (zzz-DTTCompatibility-class),

- 14
- affxStart, CompositePosition-method
(zzz-DTTCompatibility-class),
14
- affxStart, ReporterPosition-method
(zzz-DTTCompatibility-class),
14
- affxStart, SequencePosition-method
(zzz-DTTCompatibility-class),
14
- affxStatus
(zzz-NetAffxCompatibility-class),
16
- affxStatus, AffxNetAffxAnnotFileList-method
(zzz-NetAffxCompatibility-class),
16
- AffxStatus-class
(zzz-NetAffxCompatibility-class),
16
- affxText (zzz-DTTCompatibility-class),
14
- affxText, Description-method
(zzz-DTTCompatibility-class),
14
- affxText, Protocol-method
(zzz-DTTCompatibility-class),
14
- affxTitle (zzz-DTTCompatibility-class),
14
- affxTitle, BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- affxTitle, Protocol-method
(zzz-DTTCompatibility-class),
14
- affxType
(zzz-NetAffxCompatibility-class),
16
- affxType, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
16
- affxType, AffxFile-method
(zzz-NetAffxCompatibility-class),
16
- affxUrl
(zzz-NetAffxCompatibility-class),
16
- affxUrl, AffxFile-method
(zzz-NetAffxCompatibility-class),
16
- affxUrl, NetAffxResource-method
(NetAffxResource-class), 6
- AffxUrl-class
(zzz-NetAffxCompatibility-class),
16
- affyBarcode
(zzz-ARRCompatibility-class),
13
- affyBarcode, PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- AffyCompatible, 6, 14–16
- AffyCompatible
(AffyCompatible-package), 2
- AffyCompatible-class, 3
- AffyCompatible-package, 2
- analysisResults_assnreflist
(zzz-DTTCompatibility-class),
14
- analysisResults_assnreflist, Experiment-method
(zzz-DTTCompatibility-class),
14
- AnalysisResults_assnreflist-class
(zzz-DTTCompatibility-class),
14
- annotations_assnlist
(zzz-DTTCompatibility-class),
14
- annotations_assnlist, Description-method
(zzz-DTTCompatibility-class),
14
- annotations_assnlist, ExperimentalFactor-method
(zzz-DTTCompatibility-class),
14
- Annotations_assnlist-class
(zzz-DTTCompatibility-class),
14
- annotationsFile (accessors), 3
- annotationsFile, NetAffxResource-method
(NetAffxResource-class), 6
- Array-class
(zzz-DTTCompatibility-class),
14
- array_assnlist
(zzz-DTTCompatibility-class),
14

array_assnlist,Array_package-method
 (zzz-DTTCompatibility-class),
 14
 Array_assnlist-class
 (zzz-DTTCompatibility-class),
 14
 array_assnref
 (zzz-DTTCompatibility-class),
 14
 array_assnref,BioAssayCreation-method
 (zzz-DTTCompatibility-class),
 14
 array_assnref,Hybridization-method
 (zzz-DTTCompatibility-class),
 14
 Array_assnref-class
 (zzz-DTTCompatibility-class),
 14
 array_package
 (zzz-DTTCompatibility-class),
 14
 array_package,MAGE-ML-method
 (zzz-DTTCompatibility-class),
 14
 Array_package-class
 (zzz-DTTCompatibility-class),
 14
 array_ref (zzz-DTTCompatibility-class),
 14
 array_ref,Array_assnref-method
 (zzz-DTTCompatibility-class),
 14
 array_ref,Arrays_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
 Array_ref-class
 (zzz-DTTCompatibility-class),
 14
 arrayAttribute
 (zzz-ARRCompatibility-class),
 13
 arrayAttribute,PhysicalArray-method
 (zzz-ARRCompatibility-class),
 13
 ArrayAttribute-class
 (zzz-ARRCompatibility-class),
 13
 arrayDesign
 (zzz-DTTCompatibility-class),
 14
 arrayDesign,ArrayDesign_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 ArrayDesign-class
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_assnlist
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_assnlist,ArrayDesign_package-method
 (zzz-DTTCompatibility-class),
 14
 ArrayDesign_assnlist-class
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_assnref
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_assnref,Array-method
 (zzz-DTTCompatibility-class),
 14
 ArrayDesign_assnref-class
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_package
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_package,MAGE-ML-method
 (zzz-DTTCompatibility-class),
 14
 ArrayDesign_package-class
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_ref
 (zzz-DTTCompatibility-class),
 14
 arrayDesign_ref,ArrayDesign_assnref-method
 (zzz-DTTCompatibility-class),
 14
 ArrayDesign_ref-class
 (zzz-DTTCompatibility-class),
 14
 arrayGroup
 (zzz-DTTCompatibility-class),
 14
 arrayGroup,ArrayGroup_assnlist-method

- (zzz-DTTCompatibility-class),
14
- ArrayGroup-class
(zzz-DTTCompatibility-class),
14
- arrayGroup_assnlist
(zzz-DTTCompatibility-class),
14
- arrayGroup_assnlist,Array_package-method
(zzz-DTTCompatibility-class),
14
- ArrayGroup_assnlist-class
(zzz-DTTCompatibility-class),
14
- arrayGroup_assnref
(zzz-DTTCompatibility-class),
14
- arrayGroup_assnref,Array-method
(zzz-DTTCompatibility-class),
14
- ArrayGroup_assnref-class
(zzz-DTTCompatibility-class),
14
- arrayGroup_ref
(zzz-DTTCompatibility-class),
14
- arrayGroup_ref,ArrayGroup_assnref-method
(zzz-DTTCompatibility-class),
14
- ArrayGroup_ref-class
(zzz-DTTCompatibility-class),
14
- arrayIdentifier
(zzz-DTTCompatibility-class),
14
- arrayIdentifier,Array-method
(zzz-DTTCompatibility-class),
14
- arrayManufacture
(zzz-DTTCompatibility-class),
14
- arrayManufacture,ArrayManufacture_assnlist-method
(zzz-DTTCompatibility-class),
14
- ArrayManufacture-class
(zzz-DTTCompatibility-class),
14
- arrayManufacture_assnlist
(zzz-DTTCompatibility-class),
14
- arrayManufacture_assnlist,Array_package-method
(zzz-DTTCompatibility-class),
14
- ArrayManufacture_assnlist-class
(zzz-DTTCompatibility-class),
14
- arrayManufacture_ref
(zzz-DTTCompatibility-class),
14
- arrayManufacture_ref,Information_assnref-method
(zzz-DTTCompatibility-class),
14
- ArrayManufacture_ref-class
(zzz-DTTCompatibility-class),
14
- arrayManufactureDeviation
(zzz-DTTCompatibility-class),
14
- arrayManufactureDeviation,ArrayManufactureDeviations_assnlist
(zzz-DTTCompatibility-class),
14
- ArrayManufactureDeviation-class
(zzz-DTTCompatibility-class),
14
- arrayManufactureDeviations_assnlist
(zzz-DTTCompatibility-class),
14
- arrayManufactureDeviations_assnlist,Array-method
(zzz-DTTCompatibility-class),
14
- ArrayManufactureDeviations_assnlist-class
(zzz-DTTCompatibility-class),
14
- arrayManufacturers_assnreflist
(zzz-DTTCompatibility-class),
14
- arrayManufacturers_assnreflist,ArrayManufacture-method
(zzz-DTTCompatibility-class),
14
- ArrayManufacturers_assnreflist-class
(zzz-DTTCompatibility-class),
14
- arrayName (zzz-ARRCompatibility-class),
13
- arrayName,PhysicalArray-method
(zzz-ARRCompatibility-class),

- 13
- arrays_assnreflist
(zzz-DTTCompatibility-class),
14
- arrays_assnreflist,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- arrays_assnreflist,ArrayManufacture-method
(zzz-DTTCompatibility-class),
14
- Arrays_assnreflist-class
(zzz-DTTCompatibility-class),
14
- ArraySetFile-class
(zzz-ARRCompatibility-class),
13
- arraySpacingX
(zzz-DTTCompatibility-class),
14
- arraySpacingX,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- arraySpacingY
(zzz-DTTCompatibility-class),
14
- arraySpacingY,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- arrayXOrigin
(zzz-DTTCompatibility-class),
14
- arrayXOrigin,Array-method
(zzz-DTTCompatibility-class),
14
- arrayYOrigin
(zzz-DTTCompatibility-class),
14
- arrayYOrigin,Array-method
(zzz-DTTCompatibility-class),
14
- ARRCompatibility, 14
- ARRCompatibility-class
(AffyCompatible-class), 3
- associations_assnlist
(zzz-DTTCompatibility-class),
14
- associations_assnlist,OntologyEntry-method
(zzz-DTTCompatibility-class),
14
- Associations_assnlist-class
(zzz-DTTCompatibility-class),
14
- audit (zzz-DTTCompatibility-class), 14
- audit,AuditTrail_assnlist-method
(zzz-DTTCompatibility-class),
14
- Audit-class
(zzz-DTTCompatibility-class),
14
- auditAndSecurity_package
(zzz-DTTCompatibility-class),
14
- auditAndSecurity_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- AuditAndSecurity_package-class
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Array-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ArrayManufacture-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Audit-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioAssayCreation-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioAssayDimension-method

- (zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioAssayMap-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioSample-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioSequence-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,BioSource-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Channel-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,CompositeGroup-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,CompositeSequence-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,CompositeSequenceDimensionMethod
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Compound-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Database-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,DerivedBioAssay-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,DerivedSignal-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Description-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Error-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ExpectedValue-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Experiment-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ExperimentDesign-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,FactorValue-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Failed-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Feature-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,FeatureDimension-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,FeatureExtraction-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Fiducial-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Hardware-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,HardwareApplication-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Hybridization-method
(zzz-DTTCompatibility-class),
14

- (zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Image-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ImageAcquisition-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ManufactureLIMSBiomaterials-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,MeasuredSignal-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Node-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,NodeContents-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Organization-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Parameter-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Person-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,PresentAbsent-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Protocol-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,PValue-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Ratio-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Reporter-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ReporterDimension-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Security-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,Software-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist,SoftwareApplication-method

- (zzz-DTTCompatibility-class),
14
- auditTrail_assnlist, SpecializedQuantitationTypeAssayAssnlist-class, BioAssay_package-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist, Transformation-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist, Treatment-method
(zzz-DTTCompatibility-class),
14
- auditTrail_assnlist, Zone-method
(zzz-DTTCompatibility-class),
14
- AuditTrail_assnlist-class
(zzz-DTTCompatibility-class),
14
- authors (zzz-DTTCompatibility-class), 14
- authors, BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- barcode (zzz-DTTCompatibility-class), 14
- barcode, ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- basis (zzz-DTTCompatibility-class), 14
- basis, SeqFeature-method
(zzz-DTTCompatibility-class),
14
- bibliographicReference
(zzz-DTTCompatibility-class),
14
- bibliographicReference, BibliographicReferences_assnlist-method
(zzz-DTTCompatibility-class),
14
- BibliographicReference-class
(zzz-DTTCompatibility-class),
14
- bibliographicReferences_assnlist
(zzz-DTTCompatibility-class),
14
- bibliographicReferences_assnlist, DescriptionBioAssayCreation-assn-method
(zzz-DTTCompatibility-class),
14
- BibliographicReferences_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssay_assnlist
(zzz-DTTCompatibility-class),
14
- bioAssay_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssay_assnref
(zzz-DTTCompatibility-class),
14
- bioAssay_assnref, BioAssayTuple-method
(zzz-DTTCompatibility-class),
14
- BioAssay_assnref-class
(zzz-DTTCompatibility-class),
14
- bioAssay_package
(zzz-DTTCompatibility-class),
14
- bioAssay_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- BioAssay_package-class
(zzz-DTTCompatibility-class),
14
- BioAssay_ref-class
(zzz-DTTCompatibility-class),
14
- bioAssayCreation
(zzz-DTTCompatibility-class),
14
- bioAssayCreation, BioAssayCreation_assn-method
(zzz-DTTCompatibility-class),
14
- BioAssayCreation-class
(zzz-DTTCompatibility-class),
14
- bioAssayCreation_assn
(zzz-DTTCompatibility-class),
14
- bioAssayCreation_assn, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
14
- BioAssayCreation_assn-class
(zzz-DTTCompatibility-class),
14
- bioAssayData_assnlist

- (zzz-DTTCompatibility-class),
14
- bioAssayData_assnlist, BioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- BioAssayData_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssayData_assnreflist
(zzz-DTTCompatibility-class),
14
- bioAssayData_assnreflist, Experiment-method
(zzz-DTTCompatibility-class),
14
- BioAssayData_assnreflist-class
(zzz-DTTCompatibility-class),
14
- bioAssayData_package
(zzz-DTTCompatibility-class),
14
- bioAssayData_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- BioAssayData_package-class
(zzz-DTTCompatibility-class),
14
- BioAssayData_ref-class
(zzz-DTTCompatibility-class),
14
- bioAssayDataCluster
(zzz-DTTCompatibility-class),
14
- bioAssayDataCluster, BioAssayDataCluster_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioAssayDataCluster-class
(zzz-DTTCompatibility-class),
14
- bioAssayDataCluster_assnlist
(zzz-DTTCompatibility-class),
14
- bioAssayDataCluster_assnlist, HigherLevelAnalysis-method
(zzz-DTTCompatibility-class),
14
- BioAssayDataCluster_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssayDataCluster_ref
(zzz-DTTCompatibility-class),
14
- bioAssayDataCluster_ref, AnalysisResults_assnreflist-method
(zzz-DTTCompatibility-class),
14
- BioAssayDataCluster_ref-class
(zzz-DTTCompatibility-class),
14
- bioAssayDataSources_assnreflist
(zzz-DTTCompatibility-class),
14
- bioAssayDataSources_assnreflist, Transformation-method
(zzz-DTTCompatibility-class),
14
- BioAssayDataSources_assnreflist-class
(zzz-DTTCompatibility-class),
14
- bioAssayDimension
(zzz-DTTCompatibility-class),
14
- bioAssayDimension, BioAssayDimension_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioAssayDimension-class
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_assnlist
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_assnlist, BioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- BioAssayDimension_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_assnref
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_assnref, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_assnref, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_assnref, NodeContents-method
(zzz-DTTCompatibility-class),
14
- BioAssayDimension_assnref-class

- (zzz-DTTCompatibility-class),
14
- bioAssayDimension_ref
(zzz-DTTCompatibility-class),
14
- bioAssayDimension_ref,BioAssayDimension_assnreflist-method
(zzz-DTTCompatibility-class),
14
- BioAssayDimension_ref-class
(zzz-DTTCompatibility-class),
14
- bioAssayFactorValues_assnreflist
(zzz-DTTCompatibility-class),
14
- bioAssayFactorValues_assnreflist,DerivedBioAssayFactorValues-assnreflist-method
(zzz-DTTCompatibility-class),
14
- bioAssayFactorValues_assnreflist,MeasuredBioAssayFactorValues-assnreflist-method
(zzz-DTTCompatibility-class),
14
- bioAssayFactorValues_assnreflist,PhysicalBioAssayFactorValues-assnreflist-method
(zzz-DTTCompatibility-class),
14
- BioAssayFactorValues_assnreflist-class
(zzz-DTTCompatibility-class),
14
- bioAssayMap
(zzz-DTTCompatibility-class),
14
- bioAssayMap,BioAssayMap_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioAssayMap-class
(zzz-DTTCompatibility-class),
14
- bioAssayMap_assnlist
(zzz-DTTCompatibility-class),
14
- bioAssayMap_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- BioAssayMap_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssayMap_ref
(zzz-DTTCompatibility-class),
14
- bioAssayMap_ref,BioAssayMaps_assnreflist-method
(zzz-DTTCompatibility-class),
14
- (zzz-DTTCompatibility-class),
14
- bioAssayMap_ref,DerivedBioAssayMap_assnreflist-method
(zzz-DTTCompatibility-class),
14
- BioAssayMap_ref-class
(zzz-DTTCompatibility-class),
14
- bioAssayMapping
(zzz-DTTCompatibility-class),
14
- bioAssayMapping,BioAssayMapping_assn-method
(zzz-DTTCompatibility-class),
14
- BioAssayMapping-class
(zzz-DTTCompatibility-class),
14
- BioAssayMapping_assn
(zzz-DTTCompatibility-class),
14
- BioAssayMapping_assn,Transformation-method
(zzz-DTTCompatibility-class),
14
- BioAssayMapping_assn-class
(zzz-DTTCompatibility-class),
14
- bioAssayMaps_assnreflist
(zzz-DTTCompatibility-class),
14
- bioAssayMaps_assnreflist,BioAssayMapping-method
(zzz-DTTCompatibility-class),
14
- BioAssayMaps_assnreflist-class
(zzz-DTTCompatibility-class),
14
- bioAssayMapTarget_assnref
(zzz-DTTCompatibility-class),
14
- bioAssayMapTarget_assnref,BioAssayMap-method
(zzz-DTTCompatibility-class),
14
- BioAssayMapTarget_assnref-class
(zzz-DTTCompatibility-class),
14
- bioAssays_assnreflist
(zzz-DTTCompatibility-class),
14
- bioAssays_assnreflist,BioAssayDimension-method
(zzz-DTTCompatibility-class),
14

- (zzz-DTTCompatibility-class),
14
- bioAssays_assnreflist, Experiment-method
(zzz-DTTCompatibility-class),
14
- BioAssays_assnreflist-class
(zzz-DTTCompatibility-class),
14
- bioAssayTreatment
(zzz-DTTCompatibility-class),
14
- bioAssayTreatment, BioAssayTreatments_assnlist
(zzz-DTTCompatibility-class),
14
- BioAssayTreatment-class
(zzz-DTTCompatibility-class),
14
- bioAssayTreatments_assnlist
(zzz-DTTCompatibility-class),
14
- bioAssayTreatments_assnlist, PhysicalBioAssayTreatment
(zzz-DTTCompatibility-class),
14
- BioAssayTreatments_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioAssayTuple
(zzz-DTTCompatibility-class),
14
- bioAssayTuple, BioAssayTuples_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioAssayTuple-class
(zzz-DTTCompatibility-class),
14
- bioAssayTuples_assnlist
(zzz-DTTCompatibility-class),
14
- bioAssayTuples_assnlist, BioDataTuples-method
(zzz-DTTCompatibility-class),
14
- BioAssayTuples_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioDataCube
(zzz-DTTCompatibility-class),
14
- bioDataCube, BioDataValues_assn-method
(zzz-DTTCompatibility-class),
14
- BioDataCube-class
(zzz-DTTCompatibility-class),
14
- bioDataTuples
(zzz-DTTCompatibility-class),
14
- bioDataTuples, BioDataValues_assn-method
(zzz-DTTCompatibility-class),
14
- BioDataTuples-class
(zzz-DTTCompatibility-class),
14
- bioDataValues_assn
(zzz-DTTCompatibility-class),
14
- bioDataValues_assn, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
14
- bioDataValues_assn, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
14
- BioDataValues_assn-class
(zzz-DTTCompatibility-class),
14
- bioEvent_package
(zzz-DTTCompatibility-class),
14
- bioEvent_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- BioEvent_package-class
(zzz-DTTCompatibility-class),
14
- biologicalCharacteristics_assnreflist
(zzz-DTTCompatibility-class),
14
- biologicalCharacteristics_assnreflist, CompositeSequence-method
(zzz-DTTCompatibility-class),
14
- BiologicalCharacteristics_assnreflist-class
(zzz-DTTCompatibility-class),
14
- bioMaterial_assnlist
(zzz-DTTCompatibility-class),
14
- bioMaterial_assnlist, BioMaterial_package-method
(zzz-DTTCompatibility-class),
14

- (zzz-DTTCompatibility-class),
14
- BioMaterial_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioMaterial_assnref
(zzz-DTTCompatibility-class),
14
- bioMaterial_assnref,BioMaterialMeasurement-method
(zzz-DTTCompatibility-class),
14
- bioMaterial_assnref,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
14
- bioMaterial_assnref,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
14
- BioMaterial_assnref-class
(zzz-DTTCompatibility-class),
14
- bioMaterial_package
(zzz-DTTCompatibility-class),
14
- bioMaterial_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- BioMaterial_package-class
(zzz-DTTCompatibility-class),
14
- BioMaterial_ref-class
(zzz-DTTCompatibility-class),
14
- bioMaterialMeasurement
(zzz-DTTCompatibility-class),
14
- bioMaterialMeasurement,SourceBioMaterialMeasurements-assnlist-method
(zzz-DTTCompatibility-class),
14
- BioMaterialMeasurement-class
(zzz-DTTCompatibility-class),
14
- bioMaterialPlateCol
(zzz-DTTCompatibility-class),
14
- bioMaterialPlateCol,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
14
- bioMaterialPlateIdentifier
(zzz-DTTCompatibility-class),
14
- bioMaterialPlateIdentifier,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
14
- bioMaterialPlateRow
(zzz-DTTCompatibility-class),
14
- bioMaterialPlateRow,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
14
- bioSample
(zzz-DTTCompatibility-class),
14
- bioSample,BioMaterial_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioSample-class
(zzz-DTTCompatibility-class),
14
- bioSample_ref
(zzz-DTTCompatibility-class),
14
- bioSample_ref,BioMaterial_assnref-method
(zzz-DTTCompatibility-class),
14
- BioSample_ref-class
(zzz-DTTCompatibility-class),
14
- bioSequence
(zzz-DTTCompatibility-class),
14
- bioSequence,BioSequence_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioSequence-class
(zzz-DTTCompatibility-class),
14
- bioSequence_assnlist
(zzz-DTTCompatibility-class),
14
- bioSequence_assnlist,BioSequence_package-method
(zzz-DTTCompatibility-class),
14
- BioSequence_assnlist-class
(zzz-DTTCompatibility-class),
14
- bioSequence_package
(zzz-DTTCompatibility-class),

- 14
- bioSequence_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- BioSequence_package-class
(zzz-DTTCompatibility-class),
14
- bioSequence_ref
(zzz-DTTCompatibility-class),
14
- bioSequence_ref, BiologicalCharacteristics_assnreflist-method
(zzz-DTTCompatibility-class),
14
- bioSequence_ref, ImmobilizedCharacteristics_assnreflist-method
(zzz-DTTCompatibility-class),
14
- BioSequence_ref-class
(zzz-DTTCompatibility-class),
14
- bioSource (zzz-DTTCompatibility-class),
14
- bioSource, BioMaterial_assnlist-method
(zzz-DTTCompatibility-class),
14
- BioSource-class
(zzz-DTTCompatibility-class),
14
- bioSource_ref
(zzz-DTTCompatibility-class),
14
- bioSource_ref, BioMaterial_assnref-method
(zzz-DTTCompatibility-class),
14
- BioSource_ref-class
(zzz-DTTCompatibility-class),
14
- bqs_package
(zzz-DTTCompatibility-class),
14
- bqs_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- BQS_package-class
(zzz-DTTCompatibility-class),
14
- category_assn
(zzz-DTTCompatibility-class),
14
- category_assn, ExperimentalFactor-method
(zzz-DTTCompatibility-class),
14
- Category_assn-class
(zzz-DTTCompatibility-class),
14
- channel (zzz-DTTCompatibility-class), 14
- channel, Channel_assnlist-method
(zzz-DTTCompatibility-class),
14
- channel, Channel_reflist-method
(zzz-DTTCompatibility-class),
14
- channel, Channel_assnlist
(zzz-DTTCompatibility-class),
14
- channel_assnlist, BioAssay_package-method
(zzz-DTTCompatibility-class),
14
- Channel_assnlist-class
(zzz-DTTCompatibility-class),
14
- channel_assnref
(zzz-DTTCompatibility-class),
14
- channel_assnref, DerivedSignal-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, Error-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, ExpectedValue-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, Failed-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, MeasuredSignal-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, PresentAbsent-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, PValue-method
(zzz-DTTCompatibility-class),
14
- channel_assnref, Ratio-method
(zzz-DTTCompatibility-class),

- 14
channel_assnref, SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
14
- Channel_assnref-class
(zzz-DTTCompatibility-class),
14
- channel_ref
(zzz-DTTCompatibility-class),
14
- channel_ref, Channel_assnref-method
(zzz-DTTCompatibility-class),
14
- channel_ref, Channels_assnreflist-method
(zzz-DTTCompatibility-class),
14
- Channel_ref-class
(zzz-DTTCompatibility-class),
14
- channels_assnreflist
(zzz-DTTCompatibility-class),
14
- channels_assnreflist, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
14
- channels_assnreflist, Image-method
(zzz-DTTCompatibility-class),
14
- channels_assnreflist, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
14
- channels_assnreflist, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
14
- Channels_assnreflist-class
(zzz-DTTCompatibility-class),
14
- characteristics_assnlist
(zzz-DTTCompatibility-class),
14
- characteristics_assnlist, BioSample-method
(zzz-DTTCompatibility-class),
14
- characteristics_assnlist, BioSource-method
(zzz-DTTCompatibility-class),
14
- characteristics_assnlist, LabeledExtract-method
(zzz-DTTCompatibility-class),
14
- 14
characteristics_assnlist-class
(zzz-DTTCompatibility-class),
14
- clusterBioAssayData_assnref
(zzz-DTTCompatibility-class),
14
- clusterBioAssayData_assnref, BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
14
- ClusterBioAssayData_assnref-class
(zzz-DTTCompatibility-class),
14
- column (zzz-DTTCompatibility-class), 14
- column, FeatureLocation-method
(zzz-DTTCompatibility-class),
14
- column, Zone-method
(zzz-DTTCompatibility-class),
14
- componentCompounds_assnlist
(zzz-DTTCompatibility-class),
14
- componentCompounds_assnlist, Compound-method
(zzz-DTTCompatibility-class),
14
- ComponentCompounds_assnlist-class
(zzz-DTTCompatibility-class),
14
- composite_assnref
(zzz-DTTCompatibility-class),
14
- composite_assnref, CompositePosition-method
(zzz-DTTCompatibility-class),
14
- Composite_assnref-class
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap, CompositeCompositeMap_assnlist-method
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap, DesignElementMap_assnlist-method
(zzz-DTTCompatibility-class),
14
- CompositeCompositeMap-class

- (zzz-DTTCompatibility-class),
14
- compositeCompositeMap_assnlist
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap_assnlist,DesignElement_package_ref,CompositeGroups_assnreflist-method
(zzz-DTTCompatibility-class),
14
- CompositeCompositeMap_assnlist-class
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap_ref
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap_ref,CompositeCompositeMap_assnreflist-method,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- compositeCompositeMap_ref,DesignElementMaps_assnreflist-method,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- CompositeCompositeMap_ref-class
(zzz-DTTCompatibility-class),
14
- compositeCompositeMaps_assnreflist
(zzz-DTTCompatibility-class),
14
- compositeCompositeMaps_assnreflist,CompositeSequenceRef-method,CompositePositionSources_assnlist-method
(zzz-DTTCompatibility-class),
14
- CompositeCompositeMaps_assnreflist-class
(zzz-DTTCompatibility-class),
14
- compositeGroup
(zzz-DTTCompatibility-class),
14
- compositeGroup,CompositeGroup_assnlist-method
(zzz-DTTCompatibility-class),
14
- CompositeGroup-class
(zzz-DTTCompatibility-class),
14
- compositeGroup_assnlist
(zzz-DTTCompatibility-class),
14
- compositeGroup_assnlist,ArrayDesign_package_methods
(zzz-DTTCompatibility-class),
14
- CompositeGroup_assnlist-class
- (zzz-DTTCompatibility-class),
14
- compositeGroup_ref
(zzz-DTTCompatibility-class),
14
- CompositeGroup_ref-class
(zzz-DTTCompatibility-class),
14
- compositeGroups_assnreflist
(zzz-DTTCompatibility-class),
14
- compositeGroups_assnreflist,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- CompositeGroups_assnreflist-class
(zzz-DTTCompatibility-class),
14
- compositePosition
(zzz-DTTCompatibility-class),
14
- compositePosition,CompositePositionSources_assnlist-method
(zzz-DTTCompatibility-class),
14
- compositePosition,Coordinate_assn-method
(zzz-DTTCompatibility-class),
14
- CompositePosition-class
(zzz-DTTCompatibility-class),
14
- compositePositionSources_assnlist
(zzz-DTTCompatibility-class),
14
- compositePositionSources_assnlist,CompositeCompositeMap-me
(zzz-DTTCompatibility-class),
14
- CompositePositionSources_assnlist-class
(zzz-DTTCompatibility-class),
14
- compositeSequence
(zzz-DTTCompatibility-class),
14
- compositeSequence,CompositeSequence_assnlist-method

| | |
|--|---|
| (zzz-DTTCompatibility-class), 14 | (zzz-DTTCompatibility-class), 14 |
| compositeSequence,DesignElement_assnlist-method (zzz-DTTCompatibility-class), 14 | compositeSequenceDimension,DesignElementDimension_assnlist (zzz-DTTCompatibility-class), 14 |
| CompositeSequence-class (zzz-DTTCompatibility-class), 14 | CompositeSequenceDimension-class (zzz-DTTCompatibility-class), 14 |
| compositeSequence_assnlist (zzz-DTTCompatibility-class), 14 | compositeSequenceDimension_ref (zzz-DTTCompatibility-class), 14 |
| compositeSequence_assnlist,DesignElement_package-method (zzz-DTTCompatibility-class), 14 | compositeSequenceDimension_ref,DesignElementDimension_assnlist (zzz-DTTCompatibility-class), 14 |
| CompositeSequence_assnlist-class (zzz-DTTCompatibility-class), 14 | CompositeSequenceDimension_ref-class (zzz-DTTCompatibility-class), 14 |
| compositeSequence_assnref (zzz-DTTCompatibility-class), 14 | compositeSequences_assnreflist (zzz-DTTCompatibility-class), 14 |
| compositeSequence_assnref,CompositeMap-method (zzz-DTTCompatibility-class), 14 | compositeSequences_assnreflist,CompositeGroup-method (zzz-DTTCompatibility-class), 14 |
| compositeSequence_assnref,ReporterCompositeMap-method (zzz-DTTCompatibility-class), 14 | compositeSequences_assnreflist,CompositeSequenceDimension (zzz-DTTCompatibility-class), 14 |
| CompositeSequence_assnref-class (zzz-DTTCompatibility-class), 14 | CompositeSequences_assnreflist-class (zzz-DTTCompatibility-class), 14 |
| compositeSequence_ref (zzz-DTTCompatibility-class), 14 | compound (zzz-DTTCompatibility-class), 14 |
| compositeSequence_ref,Composite_assnref-method (zzz-DTTCompatibility-class), 14 | compound,Compound_assnlist-method (zzz-DTTCompatibility-class), 14 |
| compositeSequence_ref,CompositeSequence_assnref-method (zzz-DTTCompatibility-class), 14 | Compound-class (zzz-DTTCompatibility-class), 14 |
| compositeSequence_ref,CompositeSequences_assnreflist-method (zzz-DTTCompatibility-class), 14 | compound_assnlist (zzz-DTTCompatibility-class), 14 |
| compositeSequence_ref,DesignElement_assnref-method (zzz-DTTCompatibility-class), 14 | compound_assnlist,BioMaterial_package-method (zzz-DTTCompatibility-class), 14 |
| CompositeSequence_ref-class (zzz-DTTCompatibility-class), 14 | Compound_assnlist-class (zzz-DTTCompatibility-class), 14 |
| compositeSequenceDimension | compound_assnref (zzz-DTTCompatibility-class), |

- 14
 compound_assnref,CompoundMeasurement-method
 (zzz-DTTCompatibility-class),
 14
- Compound_assnref-class
 (zzz-DTTCompatibility-class),
 14
- compound_ref
 (zzz-DTTCompatibility-class),
 14
- compound_ref,Compound_assnref-method
 (zzz-DTTCompatibility-class),
 14
- compound_ref,Labels_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- Compound_ref-class
 (zzz-DTTCompatibility-class),
 14
- compoundIndices_assnlist
 (zzz-DTTCompatibility-class),
 14
- compoundIndices_assnlist,Compound-method
 (zzz-DTTCompatibility-class),
 14
- CompoundIndices_assnlist-class
 (zzz-DTTCompatibility-class),
 14
- compoundMeasurement
 (zzz-DTTCompatibility-class),
 14
- compoundMeasurement,ComponentCompounds_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- compoundMeasurement,CompoundMeasurements_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- CompoundMeasurement-class
 (zzz-DTTCompatibility-class),
 14
- compoundMeasurements_assnlist
 (zzz-DTTCompatibility-class),
 14
- compoundMeasurements_assnlist,Treatment-method
 (zzz-DTTCompatibility-class),
 14
- CompoundMeasurements_assnlist-class
 (zzz-DTTCompatibility-class),
- 14
 concentrationUnit
 (zzz-DTTCompatibility-class),
 14
- concentrationUnit,Unit_assn-method
 (zzz-DTTCompatibility-class),
 14
- ConcentrationUnit-class
 (zzz-DTTCompatibility-class),
 14
- ConfidenceIndicator_assnlist-class
 (zzz-DTTCompatibility-class),
 14
- ConfidenceIndicator_ref-class
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,DerivedSignal-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,Error-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,ExpectedValue-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,Failed-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,MeasuredSignal-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,PresentAbsent-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,PValue-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,Ratio-method
 (zzz-DTTCompatibility-class),
 14
- confidenceIndicators_assnreflist,SpecializedQuantitationType
 (zzz-DTTCompatibility-class),
 14
- ConfidenceIndicators_assnreflist-class
 (zzz-DTTCompatibility-class),

- 14
- contact_assnlist
 - (zzz-DTTCompatibility-class), 14
- contact_assnlist, AuditAndSecurity_package-method
 - (zzz-DTTCompatibility-class), 14
- Contact_assnlist-class
 - (zzz-DTTCompatibility-class), 14
- Contact_ref-class
 - (zzz-DTTCompatibility-class), 14
- contacts_assnreflist
 - (zzz-DTTCompatibility-class), 14
- contacts_assnreflist, Database-method
 - (zzz-DTTCompatibility-class), 14
- Contacts_assnreflist-class
 - (zzz-DTTCompatibility-class), 14
- containedFeatures_assnreflist
 - (zzz-DTTCompatibility-class), 14
- containedFeatures_assnreflist, FeatureDimension-method
 - (zzz-DTTCompatibility-class), 14
- ContainedFeatures_assnreflist-class
 - (zzz-DTTCompatibility-class), 14
- control (zzz-ARRCompatibility-class), 13
- control, UserAttribute-method
 - (zzz-ARRCompatibility-class), 13
- Control-class
 - (zzz-ARRCompatibility-class), 13
- controlFeatures_assnreflist
 - (zzz-DTTCompatibility-class), 14
- controlFeatures_assnreflist, Feature-method
 - (zzz-DTTCompatibility-class), 14
- ControlFeatures_assnreflist-class
 - (zzz-DTTCompatibility-class), 14
- controlledFeatures_assnreflist
 - (zzz-DTTCompatibility-class), 14
- controlledFeatures_assnreflist, Feature-method
 - (zzz-DTTCompatibility-class), 14
- ControlledFeatures_assnreflist-class
 - (zzz-DTTCompatibility-class), 14
- controlType_assn
 - (zzz-DTTCompatibility-class), 14
- controlType_assn, CompositeSequence-method
 - (zzz-DTTCompatibility-class), 14
- controlType_assn, Feature-method
 - (zzz-DTTCompatibility-class), 14
- controlType_assn, Reporter-method
 - (zzz-DTTCompatibility-class), 14
- ControlType_assn-class
 - (zzz-DTTCompatibility-class), 14
- coordinate_assn
 - (zzz-DTTCompatibility-class), 14
- coordinate_assn, SeqFeatureLocation-method
 - (zzz-DTTCompatibility-class), 14
- Coordinate_assn-class
 - (zzz-DTTCompatibility-class), 14
- createdBy (zzz-ARRCompatibility-class), 13
- createdBy, ArraySetFile-method
 - (zzz-ARRCompatibility-class), 13
- createdBy, PhysicalArray-method
 - (zzz-ARRCompatibility-class), 13
- createdBy, TemplateFile-method
 - (zzz-ARRCompatibility-class), 13
- createdDateTime
 - (zzz-ARRCompatibility-class), 13
- createdDateTime, ArraySetFile-method
 - (zzz-ARRCompatibility-class), 13

- 13
- createdDateTime,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- createdDateTime,TemplateFile-method
(zzz-ARRCompatibility-class),
13
- createdStep
(zzz-ARRCompatibility-class),
13
- createdStep,ArraySetFile-method
(zzz-ARRCompatibility-class),
13
- createdStep,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- database (zzz-DTTCompatibility-class),
14
- database,Database_assnlist-method
(zzz-DTTCompatibility-class),
14
- Database-class
(zzz-DTTCompatibility-class),
14
- database_assnlist
(zzz-DTTCompatibility-class),
14
- database_assnlist,Description_package-method
(zzz-DTTCompatibility-class),
14
- Database_assnlist-class
(zzz-DTTCompatibility-class),
14
- database_assnref
(zzz-DTTCompatibility-class),
14
- database_assnref,DatabaseEntry-method
(zzz-DTTCompatibility-class),
14
- Database_assnref-class
(zzz-DTTCompatibility-class),
14
- database_ref
(zzz-DTTCompatibility-class),
14
- database_ref,Database_assnref-method
(zzz-DTTCompatibility-class),
14
- Database_ref-class
(zzz-DTTCompatibility-class),
14
- databaseEntry
(zzz-DTTCompatibility-class),
14
- databaseEntry,Accessions_assnlist-method
(zzz-DTTCompatibility-class),
14
- databaseEntry,DatabaseReferences_assnlist-method
(zzz-DTTCompatibility-class),
14
- databaseEntry,ExternalLIMS_assn-method
(zzz-DTTCompatibility-class),
14
- databaseEntry,IdentifierLIMS_assn-method
(zzz-DTTCompatibility-class),
14
- databaseEntry,OntologyReference_assn-method
(zzz-DTTCompatibility-class),
14
- databaseEntry,SequenceDatabases_assnlist-method
(zzz-DTTCompatibility-class),
14
- DatabaseEntry-class
(zzz-DTTCompatibility-class),
14
- databaseReferences_assnlist
(zzz-DTTCompatibility-class),
14
- databaseReferences_assnlist,Description-method
(zzz-DTTCompatibility-class),
14
- DatabaseReferences_assnlist-class
(zzz-DTTCompatibility-class),
14
- dataExternal
(zzz-DTTCompatibility-class),
14
- dataExternal,DataExternal_assn-method
(zzz-DTTCompatibility-class),
14
- DataExternal-class
(zzz-DTTCompatibility-class),
14
- dataExternal_assn
(zzz-DTTCompatibility-class),
14

- dataExternal_assn, BioDataCube-method
(zzz-DTTCompatibility-class),
[14](#)
- DataExternal_assn-class
(zzz-DTTCompatibility-class),
[14](#)
- dataFormat
(zzz-DTTCompatibility-class),
[14](#)
- dataFormat, DataExternal-method
(zzz-DTTCompatibility-class),
[14](#)
- dataFormatInfoURI
(zzz-DTTCompatibility-class),
[14](#)
- dataFormatInfoURI, DataExternal-method
(zzz-DTTCompatibility-class),
[14](#)
- dataInternal
(zzz-DTTCompatibility-class),
[14](#)
- dataInternal, DataInternal_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- DataInternal-class
(zzz-DTTCompatibility-class),
[14](#)
- dataInternal_assn
(zzz-DTTCompatibility-class),
[14](#)
- dataInternal_assn, BioDataCube-method
(zzz-DTTCompatibility-class),
[14](#)
- DataInternal_assn-class
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, DerivedSignal-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, Error-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, ExpectedValue-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, Failed-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, MeasuredSignal-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, NodeValue-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, Parameter-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, PresentAbsent-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, PValue-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, Ratio-method
(zzz-DTTCompatibility-class),
[14](#)
- dataType_assn, SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
[14](#)
- DataType_assn-class
(zzz-DTTCompatibility-class),
[14](#)
- datum (zzz-DTTCompatibility-class), [14](#)
- datum, Datum_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- Datum-class
(zzz-DTTCompatibility-class),
[14](#)
- datum_assn
(zzz-DTTCompatibility-class),
[14](#)
- datum_assn, QuantitationTypeTuple-method
(zzz-DTTCompatibility-class),
[14](#)
- Datum_assn-class
(zzz-DTTCompatibility-class),
[14](#)
- defaultValue
(zzz-ARRCompatibility-class),
[13](#)
- defaultValue, UserAttribute-method
(zzz-ARRCompatibility-class),

- [13](#)
 defaultValue_assn
 (zzz-DTTCompatibility-class),
 [14](#)
 defaultValue_assn,Parameter-method
 (zzz-DTTCompatibility-class),
 [14](#)
 DefaultValue_assn-class
 (zzz-DTTCompatibility-class),
 [14](#)
 defectType_assn
 (zzz-DTTCompatibility-class),
 [14](#)
 defectType_assn,FeatureDefect-method
 (zzz-DTTCompatibility-class),
 [14](#)
 defectType_assn,ZoneDefect-method
 (zzz-DTTCompatibility-class),
 [14](#)
 DefectType_assn-class
 (zzz-DTTCompatibility-class),
 [14](#)
 deltaX(zzz-DTTCompatibility-class), [14](#)
 deltaX,PositionDelta-method
 (zzz-DTTCompatibility-class),
 [14](#)
 deltaY(zzz-DTTCompatibility-class), [14](#)
 deltaY,PositionDelta-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay,BioAssay_assnlist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay,DerivedBioAssay_assnlist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 DerivedBioAssay-class
 (zzz-DTTCompatibility-class),
 [14](#)
 DerivedBioAssay_assnlist-class
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay_ref
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay_ref,BioAssay_assnref-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay_ref,BioAssayMapTarget_assnref-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay_ref,BioAssays_assnreflist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay_ref,SourceBioAssays_assnreflist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssay_ref,TopLevelBioAssays_assnreflist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 DerivedBioAssay_ref-class
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData,BioAssayData_assnlist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData,DerivedBioAssayData_assnlist-method
 (zzz-DTTCompatibility-class),
 [14](#)
 DerivedBioAssayData-class
 (zzz-DTTCompatibility-class),
 [14](#)
 DerivedBioAssayData_assnlist-class
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData_assnreflist
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData_assnreflist,DerivedBioAssay-method
 (zzz-DTTCompatibility-class),
 [14](#)
 DerivedBioAssayData_assnreflist-class
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData_ref
 (zzz-DTTCompatibility-class),
 [14](#)
 derivedBioAssayData_ref,BioAssayData_assnreflist-method
 (zzz-DTTCompatibility-class),
 [14](#)

| | |
|--|--|
| derivedBioAssayData_ref,BioAssayDataSources_assnref-list-method (zzz-DTTCompatibility-class), 14 | derivedSignal_ref,TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class), 14 |
| derivedBioAssayData_ref,ClusterBioAssayData_assnref-class (zzz-DTTCompatibility-class), 14 | derivedSignal_ref-class (zzz-DTTCompatibility-class), 14 |
| derivedBioAssayData_ref,DerivedBioAssayData_assnref-list-method (zzz-DTTCompatibility-class), 14 | description (zzz-DTTCompatibility-class), 14 |
| derivedBioAssayData_ref,DerivedBioAssayDataTarget_assnref-list-method (zzz-DTTCompatibility-class), 14 | description,DescriptionOptions_assnlist-method (zzz-DTTCompatibility-class), 14 |
| DerivedBioAssayData_ref-class (zzz-DTTCompatibility-class), 14 | description,NormalizationDescription_assn-method (zzz-DTTCompatibility-class), 14 |
| DerivedBioAssayDataTarget_assnref-class (zzz-DTTCompatibility-class), 14 | description,OntologyEntry-method (zzz-DTTCompatibility-class), 14 |
| derivedBioAssayMap_assnreflist (zzz-DTTCompatibility-class), 14 | description,QualityControlDescription_assn-method (zzz-DTTCompatibility-class), 14 |
| derivedBioAssayMap_assnreflist,DerivedBioAssayMap_assnref-list-method (zzz-DTTCompatibility-class), 14 | description,ReplicateDescription_assn-method (zzz-DTTCompatibility-class), 14 |
| DerivedBioAssayMap_assnreflist-class (zzz-DTTCompatibility-class), 14 | Description-class (zzz-DTTCompatibility-class), 14 |
| derivedSignal (zzz-DTTCompatibility-class), 14 | description_package (zzz-DTTCompatibility-class), 14 |
| derivedSignal,QuantitationType_assnlist-method (zzz-DTTCompatibility-class), 14 | description_package,MAGE-ML-method (zzz-DTTCompatibility-class), 14 |
| DerivedSignal-class (zzz-DTTCompatibility-class), 14 | Description_package-class (zzz-DTTCompatibility-class), 14 |
| derivedSignal_ref (zzz-DTTCompatibility-class), 14 | descriptions_assnlist (zzz-DTTCompatibility-class), 14 |
| derivedSignal_ref,QuantitationType_assnref-method (zzz-DTTCompatibility-class), 14 | descriptions_assnlist,Array-method (zzz-DTTCompatibility-class), 14 |
| derivedSignal_ref,QuantitationTypes_assnref-list-method (zzz-DTTCompatibility-class), 14 | descriptions_assnlist,ArrayDesign-method (zzz-DTTCompatibility-class), 14 |
| derivedSignal_ref,SourcesQuantitationType_assnref-list-method (zzz-DTTCompatibility-class), 14 | descriptions_assnlist,ArrayGroup-method (zzz-DTTCompatibility-class), 14 |

descriptions_assnlist,ArrayManufacture-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Audit-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BibliographicReference-description
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioAssayCreation-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioAssayDataCluster-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioAssayDimension-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioAssayMap-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioAssayTreatment-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioSample-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioSequence-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,BioSource-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Channel-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,CompositeCompositeMap-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,CompositeGroup-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,CompositeSequence-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,CompositeSequenceDimension-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Compound-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Database-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,DerivedBioAssay-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,DerivedBioAssayData-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,DerivedSignal-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Description-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Error-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,ExpectedValue-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Experiment-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,ExperimentalFactor-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,ExperimentDesign-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,FactorValue-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Failed-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,Feature-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,FeatureDimension-method
 (zzz-DTTCompatibility-class), 14

descriptions_assnlist,FeatureExtraction-method
 (zzz-DTTCompatibility-class), 14

- descriptions_assnlist,FeatureGroup-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,FeatureReporterMap-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Fiducial-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Hardware-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,HardwareApplication-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Hybridization-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Image-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,ImageAcquisition-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,LabeledExtract-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,MAGE-ML-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,ManufactureLIMS-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,ManufactureLIMSBiomaterials-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,MeasuredBioAssay-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,MeasuredBioAssayData-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,MeasuredSignal-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Node-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,NodeContents-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Organization-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Parameter-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Person-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,PhysicalArrayDesign-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,PhysicalBioAssay-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,PresentAbsent-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Protocol-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,ProtocolApplication-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,PValue-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,QuantitationTypeDimension-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,QuantitationTypeMap-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Ratio-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Reporter-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,ReporterCompositeMap-method
 (zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,ReporterDimension-method
 (zzz-DTTCompatibility-class),
[14](#)

- descriptions_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Security-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Software-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Transformation-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Treatment-method
(zzz-DTTCompatibility-class),
[14](#)
- descriptions_assnlist,Zone-method
(zzz-DTTCompatibility-class),
[14](#)
- Descriptions_assnlist-class
(zzz-DTTCompatibility-class),
[14](#)
- DesignElement_assnlist-class
(zzz-DTTCompatibility-class),
[14](#)
- designElement_assnref
(zzz-DTTCompatibility-class),
[14](#)
- designElement_assnref,DesignElementTuple-method
(zzz-DTTCompatibility-class),
[14](#)
- DesignElement_assnref-class
(zzz-DTTCompatibility-class),
[14](#)
- designElement_package
(zzz-DTTCompatibility-class),
[14](#)
- designElement_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
[14](#)
- DesignElement_package-class
(zzz-DTTCompatibility-class),
[14](#)
- DesignElement_ref-class
(zzz-DTTCompatibility-class),
[14](#)
- designElementDimension_assnlist
(zzz-DTTCompatibility-class),
[14](#)
- designElementDimension_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class),
[14](#)
- DesignElementDimension_assnlist-class
(zzz-DTTCompatibility-class),
[14](#)
- DesignElementDimension_assnref
(zzz-DTTCompatibility-class),
[14](#)
- designElementDimension_assnref,DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
[14](#)
- designElementDimension_assnref,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
[14](#)
- designElementDimension_assnref,NodeContents-method
(zzz-DTTCompatibility-class),
[14](#)
- DesignElementDimension_assnref-class
(zzz-DTTCompatibility-class),
[14](#)
- DesignElementDimension_ref-class
(zzz-DTTCompatibility-class),
[14](#)
- DesignElementMap_assnlist-class
(zzz-DTTCompatibility-class),
[14](#)
- DesignElementMap_ref-class
(zzz-DTTCompatibility-class),
[14](#)
- designElementMapping
(zzz-DTTCompatibility-class),
[14](#)
- designElementMapping,DesignElementMapping_assn-method
(zzz-DTTCompatibility-class),
[14](#)

- DesignElementMapping-class
(zzz-DTTCompatibility-class),
14
- designElementMapping_assn
(zzz-DTTCompatibility-class),
14
- designElementMapping_assn,Transformation-method
(zzz-DTTCompatibility-class),
14
- DesignElementMapping_assn-class
(zzz-DTTCompatibility-class),
14
- designElementMaps_assnreflist
(zzz-DTTCompatibility-class),
14
- designElementMaps_assnreflist,DesignElementMapping-method
(zzz-DTTCompatibility-class),
14
- DesignElementMaps_assnreflist-class
(zzz-DTTCompatibility-class),
14
- designElementTuple
(zzz-DTTCompatibility-class),
14
- designElementTuple,DesignElementTuples_assnlist-method
(zzz-DTTCompatibility-class),
14
- DesignElementTuple-class
(zzz-DTTCompatibility-class),
14
- designElementTuples_assnlist
(zzz-DTTCompatibility-class),
14
- designElementTuples_assnlist,BioAssayTuple-method
(zzz-DTTCompatibility-class),
14
- DesignElementTuples_assnlist-class
(zzz-DTTCompatibility-class),
14
- designProviders_assnreflist
(zzz-DTTCompatibility-class),
14
- designProviders_assnreflist,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- designProviders_assnreflist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- DesignProviders_assnreflist-class
(zzz-DTTCompatibility-class),
14
- directory (accessors), 3
- directory,NetAffxResource-method
(NetAffxResource-class), 6
- DistanceUnit
(zzz-DTTCompatibility-class),
14
- distanceUnit,DistanceUnit_assn-method
(zzz-DTTCompatibility-class),
14
- distanceUnit,Unit_assn-method
(zzz-DTTCompatibility-class),
14
- DistanceUnit-class
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,FeatureGroup-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,Fiducial-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,Position-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,PositionDelta-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,Zone-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,ZoneGroup-method
(zzz-DTTCompatibility-class),
14
- distanceUnit_assn,ZoneLayout-method
(zzz-DTTCompatibility-class),
14
- DistanceUnit_assn-class
(zzz-DTTCompatibility-class),
14

- DTTCompatibility, 15
- DTTCompatibility-class
(AffyCompatible-class), 3
- editor (zzz-DTTCompatibility-class), 14
- editor, BibliographicReference-method
(zzz-DTTCompatibility-class), 14
- email (zzz-DTTCompatibility-class), 14
- email, Organization-method
(zzz-DTTCompatibility-class), 14
- email, Person-method
(zzz-DTTCompatibility-class), 14
- error (zzz-DTTCompatibility-class), 14
- error, ConfidenceIndicator_assnlist-method
(zzz-DTTCompatibility-class), 14
- error, QuantitationType_assnlist-method
(zzz-DTTCompatibility-class), 14
- Error-class
(zzz-DTTCompatibility-class), 14
- error_ref (zzz-DTTCompatibility-class), 14
- error_ref, ConfidenceIndicators_assnreflist-method
(zzz-DTTCompatibility-class), 14
- error_ref, QuantitationType_assnref-method
(zzz-DTTCompatibility-class), 14
- error_ref, QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class), 14
- error_ref, SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class), 14
- error_ref, TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class), 14
- Error_ref-class
(zzz-DTTCompatibility-class), 14
- expectedValue
(zzz-DTTCompatibility-class), 14
- expectedValue, ConfidenceIndicator_assnlist-method
(zzz-DTTCompatibility-class), 14
- expectedValue, QuantitationType_assnlist-method
(zzz-DTTCompatibility-class), 14
- ExpectedValue-class
(zzz-DTTCompatibility-class), 14
- expectedValue_ref
(zzz-DTTCompatibility-class), 14
- expectedValue_ref, ConfidenceIndicators_assnreflist-method
(zzz-DTTCompatibility-class), 14
- expectedValue_ref, QuantitationType_assnref-method
(zzz-DTTCompatibility-class), 14
- expectedValue_ref, QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class), 14
- expectedValue_ref, SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class), 14
- expectedValue_ref, TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class), 14
- ExpectedValue_ref-class
(zzz-DTTCompatibility-class), 14
- experiment
(zzz-DTTCompatibility-class), 14
- experiment, Experiment_assnlist-method
(zzz-DTTCompatibility-class), 14
- Experiment-class
(zzz-DTTCompatibility-class), 14
- experiment_assnlist
(zzz-DTTCompatibility-class), 14
- experiment_assnlist, Experiment_package-method
(zzz-DTTCompatibility-class), 14
- Experiment_assnlist-class
(zzz-DTTCompatibility-class), 14
- experiment_package

- (zzz-DTTCompatibility-class),
14
- experiment_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- Experiment_package-class
(zzz-DTTCompatibility-class),
14
- experimentalFactor
(zzz-DTTCompatibility-class),
14
- experimentalFactor, ExperimentalFactors_assnlist-method
(zzz-DTTCompatibility-class),
14
- ExperimentalFactor-class
(zzz-DTTCompatibility-class),
14
- ExperimentalFactor_assnref-class
(zzz-DTTCompatibility-class),
14
- experimentalFactor_ref
(zzz-DTTCompatibility-class),
14
- experimentalFactor_ref, ExperimentalFactor_assnref-method
(zzz-DTTCompatibility-class),
14
- ExperimentalFactor_ref-class
(zzz-DTTCompatibility-class),
14
- experimentalFactors_assnlist
(zzz-DTTCompatibility-class),
14
- experimentalFactors_assnlist, ExperimentDesign-method
(zzz-DTTCompatibility-class),
14
- ExperimentalFactors_assnlist-class
(zzz-DTTCompatibility-class),
14
- experimentDesign
(zzz-DTTCompatibility-class),
14
- experimentDesign, ExperimentDesigns_assnlist-method
(zzz-DTTCompatibility-class),
14
- ExperimentDesign-class
(zzz-DTTCompatibility-class),
14
- experimentDesigns_assnlist
(zzz-DTTCompatibility-class),
14
- experimentDesigns_assnlist, Experiment-method
(zzz-DTTCompatibility-class),
14
- ExperimentDesigns_assnlist-class
(zzz-DTTCompatibility-class),
14
- exportedFromDB
(zzz-DTTCompatibility-class),
14
- exportedFromDB, ExternalReference-method
(zzz-DTTCompatibility-class),
14
- exportedFromServer
(zzz-DTTCompatibility-class),
14
- exportedFromServer, ExternalReference-method
(zzz-DTTCompatibility-class),
14
- exportID (zzz-DTTCompatibility-class),
14
- exportID, ExternalReference-method
(zzz-DTTCompatibility-class),
14
- exportName
(zzz-DTTCompatibility-class),
14
- exportName, ExternalReference-method
(zzz-DTTCompatibility-class),
14
- externalLIMS_assn
(zzz-DTTCompatibility-class),
14
- externalLIMS_assn, Compound-method
(zzz-DTTCompatibility-class),
14
- ExternalLIMS_assn-class
(zzz-DTTCompatibility-class),
14
- externalReference
(zzz-DTTCompatibility-class),
14
- externalReference, ExternalReference_assn-method
(zzz-DTTCompatibility-class),
14
- externalReference_assn
(zzz-DTTCompatibility-class),

- 14
- externalReference_assn,Description-method
(zzz-DTTCompatibility-class),
14
- ExternalReference_assn-class
(zzz-DTTCompatibility-class),
14
- factorValue
(zzz-DTTCompatibility-class),
14
- factorValue,FactorValues_assnlist-method
(zzz-DTTCompatibility-class),
14
- FactorValue-class
(zzz-DTTCompatibility-class),
14
- factorValue_ref
(zzz-DTTCompatibility-class),
14
- factorValue_ref,BioAssayFactorValues_assnreflist-method
(zzz-DTTCompatibility-class),
14
- FactorValue_ref-class
(zzz-DTTCompatibility-class),
14
- factorValues_assnlist
(zzz-DTTCompatibility-class),
14
- factorValues_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class),
14
- FactorValues_assnlist-class
(zzz-DTTCompatibility-class),
14
- failed(zzz-DTTCompatibility-class), 14
- failed,QuantitationType_assnlist-method
(zzz-DTTCompatibility-class),
14
- Failed-class
(zzz-DTTCompatibility-class),
14
- failed_ref
(zzz-DTTCompatibility-class),
14
- failed_ref,QuantitationType_assnref-method
(zzz-DTTCompatibility-class),
14
- failed_ref,QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class),
14
- failed_ref,SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
14
- failed_ref,TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class),
14
- Failed_ref-class
(zzz-DTTCompatibility-class),
14
- failTypes_assnlist
(zzz-DTTCompatibility-class),
14
- failTypes_assnlist,Reporter-method
(zzz-DTTCompatibility-class),
14
- FailTypes_assnlist-class
(zzz-DTTCompatibility-class),
14
- fax(zzz-DTTCompatibility-class), 14
- fax,Organization-method
(zzz-DTTCompatibility-class),
14
- fax,Person-method
(zzz-DTTCompatibility-class),
14
- feature(zzz-DTTCompatibility-class), 14
- feature,DesignElement_assnlist-method
(zzz-DTTCompatibility-class),
14
- feature,Features_assnlist-method
(zzz-DTTCompatibility-class),
14
- Feature-class
(zzz-DTTCompatibility-class),
14
- feature_assnref
(zzz-DTTCompatibility-class),
14
- feature_assnref,FeatureDefect-method
(zzz-DTTCompatibility-class),
14
- feature_assnref,FeatureInformation-method
(zzz-DTTCompatibility-class),
14
- feature_assnref,ManufactureLIMS-method
(zzz-DTTCompatibility-class),

- 14
 feature_assnref,ManufactureLIMSBiomaterial-method
 (zzz-DTTCompatibility-class),
 14
- Feature_assnref-class
 (zzz-DTTCompatibility-class),
 14
- feature_ref
 (zzz-DTTCompatibility-class),
 14
- feature_ref,ContainedFeatures_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- feature_ref,ControlFeatures_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- feature_ref,ControlledFeatures_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- feature_ref,DesignElement_assnref-method
 (zzz-DTTCompatibility-class),
 14
- feature_ref,Feature_assnref-method
 (zzz-DTTCompatibility-class),
 14
- Feature_ref-class
 (zzz-DTTCompatibility-class),
 14
- featureDefect
 (zzz-DTTCompatibility-class),
 14
- featureDefect,FeatureDefects_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- FeatureDefect-class
 (zzz-DTTCompatibility-class),
 14
- featureDefects_assnlist
 (zzz-DTTCompatibility-class),
 14
- featureDefects_assnlist,ArrayManufactureDevialt
 (zzz-DTTCompatibility-class),
 14
- FeatureDefects_assnlist-class
 (zzz-DTTCompatibility-class),
 14
- featureDimension
 (zzz-DTTCompatibility-class),
 14
- 14
 featureDimension,DesignElementDimension_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- FeatureDimension-class
 (zzz-DTTCompatibility-class),
 14
- featureDimension_ref
 (zzz-DTTCompatibility-class),
 14
- featureDimension_ref,DesignElementDimension_assnref-method
 (zzz-DTTCompatibility-class),
 14
- FeatureDimension_ref-class
 (zzz-DTTCompatibility-class),
 14
- featureExtraction
 (zzz-DTTCompatibility-class),
 14
- featureExtraction,FeatureExtraction_assn-method
 (zzz-DTTCompatibility-class),
 14
- FeatureExtraction-class
 (zzz-DTTCompatibility-class),
 14
- featureExtraction_assn
 (zzz-DTTCompatibility-class),
 14
- featureExtraction_assn,MeasuredBioAssay-method
 (zzz-DTTCompatibility-class),
 14
- FeatureExtraction_assn-class
 (zzz-DTTCompatibility-class),
 14
- featureGroup
 (zzz-DTTCompatibility-class),
 14
- featureGroup,FeatureGroups_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- featureGroup-class
 (zzz-DTTCompatibility-class),
 14
- FeatureGroup_assnref-class
 (zzz-DTTCompatibility-class),
 14
- featureGroup_ref
 (zzz-DTTCompatibility-class),
 14

- 14
featureGroup_ref,FeatureGroup_assnref-method (zzz-DTTCompatibility-class), 14
- FeatureGroup_ref-class (zzz-DTTCompatibility-class), 14
- featureGroups_assnlist (zzz-DTTCompatibility-class), 14
- featureGroups_assnlist,ArrayDesign-method (zzz-DTTCompatibility-class), 14
- featureGroups_assnlist,PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 14
- FeatureGroups_assnlist-class (zzz-DTTCompatibility-class), 14
- featureHeight (zzz-DTTCompatibility-class), 14
- featureHeight,FeatureGroup-method (zzz-DTTCompatibility-class), 14
- featureInformation (zzz-DTTCompatibility-class), 14
- featureInformation,FeatureInformationSources_assnlist-method (zzz-DTTCompatibility-class), 14
- FeatureInformation-class (zzz-DTTCompatibility-class), 14
- featureInformationSources_assnlist (zzz-DTTCompatibility-class), 14
- featureInformationSources_assnlist,FeatureReporterMap-method (zzz-DTTCompatibility-class), 14
- FeatureInformationSources_assnlist-class (zzz-DTTCompatibility-class), 14
- featureLength (zzz-DTTCompatibility-class), 14
- featureLength,FeatureGroup-method (zzz-DTTCompatibility-class), 14
- featureLIMs_assnlist (zzz-DTTCompatibility-class), 14
- featureLIMs_assnlist,ArrayManufacture-method (zzz-DTTCompatibility-class), 14
- FeatureLIMs_assnlist-class (zzz-DTTCompatibility-class), 14
- featureLocation (zzz-DTTCompatibility-class), 14
- featureLocation,FeatureLocation_assn-method (zzz-DTTCompatibility-class), 14
- FeatureLocation-class (zzz-DTTCompatibility-class), 14
- featureLocation_assn (zzz-DTTCompatibility-class), 14
- featureLocation_assn,Feature-method (zzz-DTTCompatibility-class), 14
- FeatureLocation_assn-class (zzz-DTTCompatibility-class), 14
- FeatureReporterMap (zzz-DTTCompatibility-class), 14
- featureReporterMap,DesignElementMap_assnlist-method (zzz-DTTCompatibility-class), 14
- featureReporterMap,FeatureReporterMap_assnlist-method (zzz-DTTCompatibility-class), 14
- FeatureReporterMap-class (zzz-DTTCompatibility-class), 14
- featureReporterMap_assnlist (zzz-DTTCompatibility-class), 14
- featureReporterMap_assnlist,DesignElement_package-method (zzz-DTTCompatibility-class), 14
- FeatureReporterMap_assnlist-class (zzz-DTTCompatibility-class), 14

- 14
- featureReporterMap_ref
(zzz-DTTCompatibility-class),
14
- featureReporterMap_ref,DesignElementMaps_assnreflist-method
(zzz-DTTCompatibility-class),
14
- featureReporterMap_ref,FeatureReporterMaps_assnreflist-method
(zzz-DTTCompatibility-class),
14
- FeatureReporterMap_ref-class
(zzz-DTTCompatibility-class),
14
- featureReporterMaps_assnreflist
(zzz-DTTCompatibility-class),
14
- featureReporterMaps_assnreflist,Reporter-method
(zzz-DTTCompatibility-class),
14
- FeatureReporterMaps_assnreflist-class
(zzz-DTTCompatibility-class),
14
- features_assnlist
(zzz-DTTCompatibility-class),
14
- features_assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class),
14
- Features_assnlist-class
(zzz-DTTCompatibility-class),
14
- featureShape_assn
(zzz-DTTCompatibility-class),
14
- featureShape_assn,FeatureGroup-method
(zzz-DTTCompatibility-class),
14
- FeatureShape_assn-class
(zzz-DTTCompatibility-class),
14
- featureWidth
(zzz-DTTCompatibility-class),
14
- featureWidth,FeatureGroup-method
(zzz-DTTCompatibility-class),
14
- fiducial (zzz-DTTCompatibility-class),
14
- fiducial,Fiducials_assnlist-method
(zzz-DTTCompatibility-class),
14
- Fiducial-class
(zzz-DTTCompatibility-class),
14
- fiducials_assnlist
(zzz-DTTCompatibility-class),
14
- fiducials_assnlist,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- Fiducials_assnlist-class
(zzz-DTTCompatibility-class),
14
- fiducialType_assn
(zzz-DTTCompatibility-class),
14
- fiducialType_assn,Fiducial-method
(zzz-DTTCompatibility-class),
14
- FiducialType_assn-class
(zzz-DTTCompatibility-class),
14
- filenameURI
(zzz-DTTCompatibility-class),
14
- filenameURI,DataExternal-method
(zzz-DTTCompatibility-class),
14
- firstName (zzz-DTTCompatibility-class),
14
- firstName,Person-method
(zzz-DTTCompatibility-class),
14
- format_assn
(zzz-DTTCompatibility-class),
14
- format_assn,Image-method
(zzz-DTTCompatibility-class),
14
- Format_assn-class
(zzz-DTTCompatibility-class),
14
- guid (zzz-ARRCompatibility-class), 13
- guid,ArraySetFile-method
(zzz-ARRCompatibility-class),
13

- guid,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- guid,TemplateFile-method
(zzz-ARRCompatibility-class),
13
- hardware (zzz-DTTCompatibility-class),
14
- hardware,Hardware_assnlist-method
(zzz-DTTCompatibility-class),
14
- Hardware-class
(zzz-DTTCompatibility-class),
14
- hardware_assnlist
(zzz-DTTCompatibility-class),
14
- hardware_assnlist,Protocol_package-method
(zzz-DTTCompatibility-class),
14
- Hardware_assnlist-class
(zzz-DTTCompatibility-class),
14
- hardware_assnref
(zzz-DTTCompatibility-class),
14
- hardware_assnref,HardwareApplication-method
(zzz-DTTCompatibility-class),
14
- hardware_assnref,Software-method
(zzz-DTTCompatibility-class),
14
- Hardware_assnref-class
(zzz-DTTCompatibility-class),
14
- hardware_ref
(zzz-DTTCompatibility-class),
14
- hardware_ref,Hardware_assnref-method
(zzz-DTTCompatibility-class),
14
- hardware_ref,Hardwares_assnreflist-method
(zzz-DTTCompatibility-class),
14
- Hardware_ref-class
(zzz-DTTCompatibility-class),
14
- hardwareApplication
(zzz-DTTCompatibility-class),
14
- hardwareApplication,HardwareApplications_assnlist-method
(zzz-DTTCompatibility-class),
14
- HardwareApplication-class
(zzz-DTTCompatibility-class),
14
- hardwareApplications_assnlist
(zzz-DTTCompatibility-class),
14
- hardwareApplications_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
14
- HardwareApplications_assnlist-class
(zzz-DTTCompatibility-class),
14
- hardwareManufacturers_assnreflist
(zzz-DTTCompatibility-class),
14
- hardwareManufacturers_assnreflist,Hardware-method
(zzz-DTTCompatibility-class),
14
- HardwareManufacturers_assnreflist-class
(zzz-DTTCompatibility-class),
14
- hardwares_assnreflist
(zzz-DTTCompatibility-class),
14
- hardwares_assnreflist,Protocol-method
(zzz-DTTCompatibility-class),
14
- Hardwares_assnreflist-class
(zzz-DTTCompatibility-class),
14
- higherLevelAnalysis_package
(zzz-DTTCompatibility-class),
14
- higherLevelAnalysis_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- HigherLevelAnalysis_package-class
(zzz-DTTCompatibility-class),
14
- hybridization
(zzz-DTTCompatibility-class),
14
- hybridization,BioAssayCreation_assn-method

- (zzz-DTTCompatibility-class),
14
- Hybridization-class
(zzz-DTTCompatibility-class),
14
- identifier
(zzz-DTTCompatibility-class),
14
- identifier,Array-method
(zzz-DTTCompatibility-class),
14
- identifier,Array_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- identifier,ArrayDesign_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- identifier,ArrayGroup_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,ArrayManufacture-method
(zzz-DTTCompatibility-class),
14
- identifier,ArrayManufacture_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssay_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayCreation-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayData_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayDataCluster_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayDimension-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayDimension_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayMap-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayMap_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
14
- identifier,BioMaterial_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioSample-method
(zzz-DTTCompatibility-class),
14
- identifier,BioSample_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioSequence-method
(zzz-DTTCompatibility-class),
14
- identifier,BioSequence_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,BioSource-method
(zzz-DTTCompatibility-class),
14
- identifier,BioSource_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Channel-method
(zzz-DTTCompatibility-class),
14
- identifier,Channel_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
14
- identifier,CompositeCompositeMap_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,CompositeGroup-method

- (zzz-DTTCompatibility-class),
14
- identifier, CompositeGroup_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, CompositeSequence-method
(zzz-DTTCompatibility-class),
14
- identifier, CompositeSequence_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
14
- identifier, CompositeSequenceDimension_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Compound-method
(zzz-DTTCompatibility-class),
14
- identifier, Compound_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, ConfidenceIndicator_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Contact_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Database-method
(zzz-DTTCompatibility-class),
14
- identifier, Database_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
14
- identifier, DerivedBioAssay_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
14
- identifier, DerivedBioAssayData_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, DerivedSignal-method
(zzz-DTTCompatibility-class),
14
- identifier, DerivedSignal_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, DesignElement_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, DesignElementDimension_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, DesignElementMap_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Error-method
(zzz-DTTCompatibility-class),
14
- identifier, Error_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, ExpectedValue-method
(zzz-DTTCompatibility-class),
14
- identifier, ExpectedValue_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Experiment-method
(zzz-DTTCompatibility-class),
14
- identifier, ExperimentalFactor-method
(zzz-DTTCompatibility-class),
14
- identifier, ExperimentalFactor_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, FactorValue-method
(zzz-DTTCompatibility-class),
14
- identifier, FactorValue_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Failed-method
(zzz-DTTCompatibility-class),
14
- identifier, Failed_ref-method
(zzz-DTTCompatibility-class),
14
- identifier, Feature-method

- (zzz-DTTCompatibility-class),
[14](#)
- identifier, Feature_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureDimension-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureDimension_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureExtraction-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureGroup-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureGroup_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, FeatureReporterMap_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Hardware-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Hardware_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Hybridization-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Image-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Image_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, ImageAcquisition-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, LabeledExtract-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, LabeledExtract_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MAGE-ML-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MeasuredBioAssay_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MeasuredBioAssayData_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MeasuredSignal-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, MeasuredSignal_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Organization-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Organization_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Parameter-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Parameter_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Person-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, Person_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, PhysicalArrayDesign_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- identifier, PhysicalBioAssay-method

- (zzz-DTTCompatibility-class),
14
- identifier,PhysicalBioAssay_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,PresentAbsent-method
(zzz-DTTCompatibility-class),
14
- identifier,PresentAbsent_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Protocol-method
(zzz-DTTCompatibility-class),
14
- identifier,Protocol_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,PValue-method
(zzz-DTTCompatibility-class),
14
- identifier,PValue_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,QuantitationType_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
14
- identifier,QuantitationTypeDimension_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
14
- identifier,QuantitationTypeMap_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Ratio-method
(zzz-DTTCompatibility-class),
14
- identifier,Ratio_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Reporter-method
(zzz-DTTCompatibility-class),
14
- identifier,Reporter_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
14
- identifier,ReporterCompositeMap_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,ReporterDimension-method
(zzz-DTTCompatibility-class),
14
- identifier,ReporterDimension_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,ReporterGroup-method
(zzz-DTTCompatibility-class),
14
- identifier,ReporterGroup_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Security-method
(zzz-DTTCompatibility-class),
14
- identifier,Security_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,SecurityGroup-method
(zzz-DTTCompatibility-class),
14
- identifier,SecurityGroup_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Software-method
(zzz-DTTCompatibility-class),
14
- identifier,Software_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
14
- identifier,SpecializedQuantitationType_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,StandardQuantitationType_ref-method
(zzz-DTTCompatibility-class),
14
- identifier,Transformation-method

- (zzz-DTTCompatibility-class),
14
- identifier,Treatment-method
(zzz-DTTCompatibility-class),
14
- identifier,Zone-method
(zzz-DTTCompatibility-class),
14
- identifier,Zone_ref-method
(zzz-DTTCompatibility-class),
14
- identifierLIMS_assn
(zzz-DTTCompatibility-class),
14
- identifierLIMS_assn,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
14
- identifierLIMS_assn,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
14
- IdentifierLIMS_assn-class
(zzz-DTTCompatibility-class),
14
- Image-class
(zzz-DTTCompatibility-class),
14
- image_ref (zzz-DTTCompatibility-class),
14
- image_ref,Images_assnreflist-method
(zzz-DTTCompatibility-class),
14
- Image_ref-class
(zzz-DTTCompatibility-class),
14
- imageAcquisition
(zzz-DTTCompatibility-class),
14
- imageAcquisition,BioAssayTreatments_assnlist-method
(zzz-DTTCompatibility-class),
14
- ImageAcquisition-class
(zzz-DTTCompatibility-class),
14
- images_assnreflist
(zzz-DTTCompatibility-class),
14
- images_assnreflist,ImageAcquisition-method
(zzz-DTTCompatibility-class),
14
- 14
- Images_assnreflist-class
(zzz-DTTCompatibility-class),
14
- immobilizedCharacteristics_assnreflist
(zzz-DTTCompatibility-class),
14
- immobilizedCharacteristics_assnreflist,Reporter-method
(zzz-DTTCompatibility-class),
14
- ImmobilizedCharacteristics_assnreflist-class
(zzz-DTTCompatibility-class),
14
- information_assnref
(zzz-DTTCompatibility-class),
14
- information_assnref,Array-method
(zzz-DTTCompatibility-class),
14
- Information_assnref-class
(zzz-DTTCompatibility-class),
14
- Internal-class,4
- isApproximateLength
(zzz-DTTCompatibility-class),
14
- isApproximateLength,BioSequence-method
(zzz-DTTCompatibility-class),
14
- isBackground
(zzz-DTTCompatibility-class),
14
- isBackground,DerivedSignal-method
(zzz-DTTCompatibility-class),
14
- isBackground,Error-method
(zzz-DTTCompatibility-class),
14
- isBackground,ExpectedValue-method
(zzz-DTTCompatibility-class),
14
- isBackground,Failed-method
(zzz-DTTCompatibility-class),
14
- isBackground,MeasuredSignal-method
(zzz-DTTCompatibility-class),
14
- isBackground,PresentAbsent-method

- (zzz-DTTCompatibility-class),
14
- isBackground,PValue-method
(zzz-DTTCompatibility-class),
14
- isBackground,Ratio-method
(zzz-DTTCompatibility-class),
14
- isBackground,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
14
- isCircular
(zzz-DTTCompatibility-class),
14
- isCircular,BioSequence-method
(zzz-DTTCompatibility-class),
14
- isSolvent(zzz-DTTCompatibility-class),
14
- isSolvent,Compound-method
(zzz-DTTCompatibility-class),
14
- issue(zzz-DTTCompatibility-class), 14
- issue,BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- kindCV(zzz-DTTCompatibility-class), 14
- kindCV,Measurement-method
(zzz-DTTCompatibility-class),
14
- labeledExtract
(zzz-DTTCompatibility-class),
14
- labeledExtract,BioMaterial_assnlist-method
(zzz-DTTCompatibility-class),
14
- LabeledExtract-class
(zzz-DTTCompatibility-class),
14
- labeledExtract_ref
(zzz-DTTCompatibility-class),
14
- labeledExtract_ref,BioMaterial_assnref-method
(zzz-DTTCompatibility-class),
14
- LabeledExtract_ref-class
(zzz-DTTCompatibility-class),
14
- labels_assnreflist
(zzz-DTTCompatibility-class),
14
- labels_assnreflist,Channel-method
(zzz-DTTCompatibility-class),
14
- labels_assnreflist,LabeledExtract-method
(zzz-DTTCompatibility-class),
14
- Labels_assnreflist-class
(zzz-DTTCompatibility-class),
14
- lastName(zzz-DTTCompatibility-class),
14
- lastName,Person-method
(zzz-DTTCompatibility-class),
14
- length,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- length,BioSequence-method
(zzz-DTTCompatibility-class),
14
- libraryPackageName
(zzz-ARRCompatibility-class),
13
- libraryPackageName,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- list, 4
- lowerRightX
(zzz-DTTCompatibility-class),
14
- lowerRightX,Zone-method
(zzz-DTTCompatibility-class),
14
- lowerRightY
(zzz-DTTCompatibility-class),
14
- lowerRightY,Zone-method
(zzz-DTTCompatibility-class),
14
- MAGE-ML-class
(zzz-DTTCompatibility-class),
14
- make(zzz-DTTCompatibility-class), 14
- make,Hardware-method

- (zzz-DTTCompatibility-class),
14
- manufactureLIMS
(zzz-DTTCompatibility-class),
14
- manufactureLIMS, FeatureLIMSS_assnlist-method
(zzz-DTTCompatibility-class),
14
- ManufactureLIMS-class
(zzz-DTTCompatibility-class),
14
- manufactureLIMSBiomaterial
(zzz-DTTCompatibility-class),
14
- manufactureLIMSBiomaterial, FeatureLIMSS_assnlist-method
(zzz-DTTCompatibility-class),
14
- ManufactureLIMSBiomaterial-class
(zzz-DTTCompatibility-class),
14
- manufacturingDate
(zzz-DTTCompatibility-class),
14
- manufacturingDate, ArrayManufacture-method
(zzz-DTTCompatibility-class),
14
- massUnit (zzz-DTTCompatibility-class),
14
- massUnit, Unit_assn-method
(zzz-DTTCompatibility-class),
14
- MassUnit-class
(zzz-DTTCompatibility-class),
14
- masterFileGUID
(zzz-ARRCompatibility-class),
13
- masterFileGUID, PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- masterFileName
(zzz-ARRCompatibility-class),
13
- masterFileName, PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- materialType_assn
(zzz-DTTCompatibility-class),
14
- materialType_assn, BioSample-method
(zzz-DTTCompatibility-class),
14
- materialType_assn, BioSource-method
(zzz-DTTCompatibility-class),
14
- materialType_assn, LabeledExtract-method
(zzz-DTTCompatibility-class),
14
- MaterialType_assn-class
(zzz-DTTCompatibility-class),
14
- measuredBioAssay
(zzz-DTTCompatibility-class),
14
- measuredBioAssay, BioAssay_assnlist-method
(zzz-DTTCompatibility-class),
14
- measuredBioAssay, MeasuredBioAssay_assnlist-method
(zzz-DTTCompatibility-class),
14
- MeasuredBioAssay-class
(zzz-DTTCompatibility-class),
14
- MeasuredBioAssay_assnlist-class
(zzz-DTTCompatibility-class),
14
- measuredBioAssay_ref
(zzz-DTTCompatibility-class),
14
- measuredBioAssay_ref, BioAssay_assnref-method
(zzz-DTTCompatibility-class),
14
- measuredBioAssay_ref, BioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
14
- measuredBioAssay_ref, MeasuredBioAssayTarget_assnref-method
(zzz-DTTCompatibility-class),
14
- measuredBioAssay_ref, SourceBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
14
- measuredBioAssay_ref, TopLevelBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
14
- MeasuredBioAssay_ref-class
(zzz-DTTCompatibility-class),

[14](#)
 measuredBioAssayData
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData,BioAssayData_assnlist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData,MeasuredBioAssayData_assnlist-method
 (zzz-DTTCompatibility-class),
[14](#)
 MeasuredBioAssayData-class
 (zzz-DTTCompatibility-class),
[14](#)
 MeasuredBioAssayData_assnlist-class
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_assnreflist
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_assnreflist,MeasuredBioAssayData_assnreflist
 (zzz-DTTCompatibility-class),
[14](#)
 MeasuredBioAssayData_assnreflist-class
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_ref
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_ref,BioAssayData_assnreflist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_ref,BioAssayDataSources_assnreflist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_ref,ClusterBioAssayData_assnreflist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredBioAssayData_ref,MeasuredBioAssayData_assnreflist-method
 (zzz-DTTCompatibility-class),
[14](#)
 MeasuredBioAssayData_ref-class
 (zzz-DTTCompatibility-class),
[14](#)
 MeasuredBioAssayTarget_assnref-class
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal,QuantitationType_assnlist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal-class
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal_ref
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal_ref,QuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal_ref,QuantitationTypes_assnreflist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal_ref,SourcesQuantitationType_assnreflist-method
 (zzz-DTTCompatibility-class),
[14](#)
 measuredSignal_ref,TargetQuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
[14](#)
 MeasuredSignal_ref-class
 (zzz-DTTCompatibility-class),
[14](#)
 measurement
 (zzz-DTTCompatibility-class),
[14](#)
 measurement,ActionMeasurement_assn-method
 (zzz-DTTCompatibility-class),
[14](#)
 measurement,DefaultValue_assn-method
 (zzz-DTTCompatibility-class),
[14](#)
 measurement,Measurement_assn-method
 (zzz-DTTCompatibility-class),
[14](#)
 Measurement-class
 (zzz-DTTCompatibility-class),
[14](#)
 measurement_assn
 (zzz-DTTCompatibility-class),
[14](#)
 measurement_assn,BioMaterialMeasurement-method
 (zzz-DTTCompatibility-class),
[14](#)
 measurement_assn,CompoundMeasurement-method
 (zzz-DTTCompatibility-class),
[14](#)

- 14
- measurement_assn,FactorValue-method
(zzz-DTTCompatibility-class),
14
- Measurement_assn-class
(zzz-DTTCompatibility-class),
14
- measurement_package
(zzz-DTTCompatibility-class),
14
- measurement_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- Measurement_package-class
(zzz-DTTCompatibility-class),
14
- mediaCol (zzz-ARRCompatibility-class),
13
- mediaCol,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- mediaFileGUID
(zzz-ARRCompatibility-class),
13
- mediaFileGUID,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- mediaFileName
(zzz-ARRCompatibility-class),
13
- mediaFileName,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- mediaRow (zzz-ARRCompatibility-class),
13
- mediaRow,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- mediaType (zzz-ARRCompatibility-class),
13
- mediaType,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- members_assnreflist
(zzz-DTTCompatibility-class),
14
- members_assnreflist,SecurityGroup-method
(zzz-DTTCompatibility-class),
14
- Members_assnreflist-class
(zzz-DTTCompatibility-class),
14
- midInitials
(zzz-DTTCompatibility-class),
14
- midInitials,Person-method
(zzz-DTTCompatibility-class),
14
- mismatchInformation
(zzz-DTTCompatibility-class),
14
- mismatchInformation,MismatchInformation_assnlist-method
(zzz-DTTCompatibility-class),
14
- MismatchInformation-class
(zzz-DTTCompatibility-class),
14
- mismatchInformation_assnlist
(zzz-DTTCompatibility-class),
14
- mismatchInformation_assnlist,CompositePosition-method
(zzz-DTTCompatibility-class),
14
- mismatchInformation_assnlist,FeatureInformation-method
(zzz-DTTCompatibility-class),
14
- mismatchInformation_assnlist,ReporterPosition-method
(zzz-DTTCompatibility-class),
14
- MismatchInformation_assnlist-class
(zzz-DTTCompatibility-class),
14
- model (zzz-DTTCompatibility-class), 14
- model,Hardware-method
(zzz-DTTCompatibility-class),
14
- name (zzz-DTTCompatibility-class), 14
- name,Array-method
(zzz-DTTCompatibility-class),
14
- name,Array_ref-method
(zzz-DTTCompatibility-class),
14
- name,ArrayAttribute-method
(zzz-ARRCompatibility-class),
13

name,ArrayDesign-method
 (zzz-DTTCompatibility-class),
 14

name,ArrayDesign_ref-method
 (zzz-DTTCompatibility-class),
 14

name,ArrayGroup-method
 (zzz-DTTCompatibility-class),
 14

name,ArrayGroup_ref-method
 (zzz-DTTCompatibility-class),
 14

name,ArrayManufacture-method
 (zzz-DTTCompatibility-class),
 14

name,ArrayManufacture_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssay_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayCreation-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayData_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayDataCluster-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayDataCluster_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayDimension-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayDimension_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayMap-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayMap_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioAssayTreatment-method
 (zzz-DTTCompatibility-class),
 14

name,BioMaterial_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioSample-method
 (zzz-DTTCompatibility-class),
 14

name,BioSample_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioSequence-method
 (zzz-DTTCompatibility-class),
 14

name,BioSequence_ref-method
 (zzz-DTTCompatibility-class),
 14

name,BioSource-method
 (zzz-DTTCompatibility-class),
 14

name,BioSource_ref-method
 (zzz-DTTCompatibility-class),
 14

name,Channel-method
 (zzz-DTTCompatibility-class),
 14

name,Channel_ref-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeCompositeMap-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeCompositeMap_ref-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeGroup-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeGroup_ref-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeSequence-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeSequence_ref-method
 (zzz-DTTCompatibility-class),
 14

name,CompositeSequenceDimension-method
 (zzz-DTTCompatibility-class),
 14

name, CompositeSequenceDimension_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Compound-method
(zzz-DTTCompatibility-class),
[14](#)

name, Compound_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, ConfidenceIndicator_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Contact_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Database-method
(zzz-DTTCompatibility-class),
[14](#)

name, Database_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
[14](#)

name, DerivedBioAssay_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
[14](#)

name, DerivedBioAssayData_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, DerivedSignal-method
(zzz-DTTCompatibility-class),
[14](#)

name, DerivedSignal_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, DesignElement_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, DesignElementDimension_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, DesignElementMap_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Error-method
(zzz-DTTCompatibility-class),
[14](#)

name, Error_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, ExpectedValue-method
(zzz-DTTCompatibility-class),
[14](#)

name, ExpectedValue_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Experiment-method
(zzz-DTTCompatibility-class),
[14](#)

name, ExperimentalFactor-method
(zzz-DTTCompatibility-class),
[14](#)

name, ExperimentalFactor_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, FactorValue-method
(zzz-DTTCompatibility-class),
[14](#)

name, FactorValue_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Failed-method
(zzz-DTTCompatibility-class),
[14](#)

name, Failed_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Feature-method
(zzz-DTTCompatibility-class),
[14](#)

name, Feature_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureDimension-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureDimension_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureExtraction-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureGroup-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureGroup_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
[14](#)

name, FeatureReporterMap_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Hardware-method
(zzz-DTTCompatibility-class),
[14](#)

name, Hardware_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Hybridization-method
(zzz-DTTCompatibility-class),
[14](#)

name, Image-method
(zzz-DTTCompatibility-class),
[14](#)

name, Image_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, ImageAcquisition-method
(zzz-DTTCompatibility-class),
[14](#)

name, LabeledExtract-method
(zzz-DTTCompatibility-class),
[14](#)

name, LabeledExtract_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, MAGE-ML-method
(zzz-DTTCompatibility-class),
[14](#)

name, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
[14](#)

name, MeasuredBioAssay_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
[14](#)

name, MeasuredBioAssayData_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, MeasuredSignal-method
(zzz-DTTCompatibility-class),
[14](#)

name, MeasuredSignal_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, NameValueType-method
(zzz-DTTCompatibility-class),
[14](#)

name, NodeValue-method
(zzz-DTTCompatibility-class),
[14](#)

name, Organization-method
(zzz-DTTCompatibility-class),
[14](#)

name, Organization_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Parameter-method
(zzz-DTTCompatibility-class),
[14](#)

name, Parameter_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, Person-method
(zzz-DTTCompatibility-class),
[14](#)

name, Person_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
[14](#)

name, PhysicalArrayDesign_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
[14](#)

name, PhysicalBioAssay_ref-method
(zzz-DTTCompatibility-class),
[14](#)

name, PresentAbsent-method
(zzz-DTTCompatibility-class),
[14](#)

- name, PresentAbsent_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Protocol-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Protocol_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, PValue-method
(zzz-DTTCompatibility-class),
[14](#)
- name, PValue_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, QuantitationType_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
[14](#)
- name, QuantitationTypeDimension_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
[14](#)
- name, QuantitationTypeMap_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Ratio-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Ratio_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Reporter-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Reporter_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
[14](#)
- name, ReporterCompositeMap_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, ReporterDimension-method
(zzz-DTTCompatibility-class),
[14](#)
- name, ReporterDimension_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, ReporterGroup-method
(zzz-DTTCompatibility-class),
[14](#)
- name, ReporterGroup_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Security-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Security_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, SecurityGroup-method
(zzz-DTTCompatibility-class),
[14](#)
- name, SecurityGroup_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Software-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Software_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
[14](#)
- name, SpecializedQuantitationType_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, StandardQuantitationType_ref-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Transformation-method
(zzz-DTTCompatibility-class),
[14](#)
- name, Treatment-method
(zzz-DTTCompatibility-class),
[14](#)
- name, UserAttribute-method
(zzz-ARRCompatibility-class),
[13](#)

name, Zone-method
 (zzz-DTTCompatibility-class),
 14
 name, Zone_ref-method
 (zzz-DTTCompatibility-class),
 14
 names, AffxArray-method
 (zzz-NetAffxCompatibility-class),
 16
 names, NetAffxResource-method
 (NetAffxResource-class), 6
 namespace (zzz-ARRCompatibility-class),
 13
 namespace, UserAttribute-method
 (zzz-ARRCompatibility-class),
 13
 nameValueType
 (zzz-DTTCompatibility-class),
 14
 nameValueType, PropertySets_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 nameValueType, QualityControlStatistics_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 nameValueType, SummaryStatistics_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 NameValueType-class
 (zzz-DTTCompatibility-class),
 14
 NetAffxCompatibility, 6, 16
 NetAffxCompatibility-class
 (AffyCompatible-class), 3
 NetAffxResource, 3, 5, 5, 10, 16
 NetAffxResource-accessors (accessors), 3
 NetAffxResource-class, 6
 newSequence
 (zzz-DTTCompatibility-class),
 14
 newSequence, MismatchInformation-method
 (zzz-DTTCompatibility-class),
 14
 node (zzz-DTTCompatibility-class), 14
 node, Nodes_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 Node-class
 (zzz-DTTCompatibility-class),
 14
 nodeContents
 (zzz-DTTCompatibility-class),
 14
 nodeContents, NodeContents_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 NodeContents-class
 (zzz-DTTCompatibility-class),
 14
 nodeContents_assnlist
 (zzz-DTTCompatibility-class),
 14
 nodeContents_assnlist, Node-method
 (zzz-DTTCompatibility-class),
 14
 NodeContents_assnlist-class
 (zzz-DTTCompatibility-class),
 14
 nodes_assnlist
 (zzz-DTTCompatibility-class),
 14
 nodes_assnlist, BioAssayDataCluster-method
 (zzz-DTTCompatibility-class),
 14
 nodes_assnlist, Node-method
 (zzz-DTTCompatibility-class),
 14
 Nodes_assnlist-class
 (zzz-DTTCompatibility-class),
 14
 nodeValue (zzz-DTTCompatibility-class),
 14
 nodeValue, NodeValue_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 NodeValue-class
 (zzz-DTTCompatibility-class),
 14
 nodeValue_assnlist
 (zzz-DTTCompatibility-class),
 14
 nodeValue_assnlist, Node-method
 (zzz-DTTCompatibility-class),
 14
 NodeValue_assnlist-class
 (zzz-DTTCompatibility-class),
 14

- 14
- normalizationDescription_assn
(zzz-DTTCompatibility-class),
14
- normalizationDescription_assn,ExperimentDesign-method(zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class),
14
- NormalizationDescription_assn-class
(zzz-DTTCompatibility-class),
14
- numArrays (zzz-DTTCompatibility-class),
14
- numArrays,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- numberOfFeatures
(zzz-DTTCompatibility-class),
14
- numberOfFeatures,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- numberOfFeatures,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- numFeaturesPerCol
(zzz-DTTCompatibility-class),
14
- numFeaturesPerCol,ZoneLayout-method
(zzz-DTTCompatibility-class),
14
- numFeaturesPerRow
(zzz-DTTCompatibility-class),
14
- numFeaturesPerRow,ZoneLayout-method
(zzz-DTTCompatibility-class),
14
- ontologyEntries_assnlist
(zzz-DTTCompatibility-class),
14
- ontologyEntries_assnlist,BioSequence-method
(zzz-DTTCompatibility-class),
14
- OntologyEntries_assnlist-class
(zzz-DTTCompatibility-class),
14
- ontologyEntry
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Action_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Annotations_assnlist-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Associations_assnlist-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Category_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Characteristics_assnlist-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,CompoundIndices_assnlist-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,ControlType_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,DataType_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,DefectType_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,FailTypes_assnlist-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,FeatureShape_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,FiducialType_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Format_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,MaterialType_assn-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,OntologyEntries_assnlist-method
(zzz-DTTCompatibility-class),
14
- ontologyEntry,Parameters_assnlist-method
(zzz-DTTCompatibility-class),
14

- ontologyEntry, PolymerType_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, Roles_assnlist-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, Scale_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, Species_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, SubstrateType_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, SurfaceType_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, TechnologyType_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, Type_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, Types_assnlist-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, Value_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- ontologyEntry, WarningType_assn-method
(zzz-DTTCompatibility-class),
[14](#)
- OntologyEntry-class
(zzz-DTTCompatibility-class),
[14](#)
- ontologyReference_assn
(zzz-DTTCompatibility-class),
[14](#)
- ontologyReference_assn, OntologyEntry-method
(zzz-DTTCompatibility-class),
[14](#)
- OntologyReference_assn-class
(zzz-DTTCompatibility-class),
[14](#)
- organization
(zzz-DTTCompatibility-class),
[14](#)
- organization, Contact_assnlist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization, Organization_assnlist-method
(zzz-DTTCompatibility-class),
[14](#)
- Organization-class
(zzz-DTTCompatibility-class),
[14](#)
- Organization_assnlist-class
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Affiliation_assnref-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, ArrayManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Contacts_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, DesignProviders_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, HardwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Members_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Owner_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Parent_assnref-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Performer_assnref-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, Providers_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)
- organization_ref, SoftwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
[14](#)

- organization_ref,SourceContact_assnreflist-method
(zzz-DTTCompatibility-class),
14
- Organization_ref-class
(zzz-DTTCompatibility-class),
14
- orientationMark
(zzz-DTTCompatibility-class),
14
- orientationMark,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- orientationMarkPosition
(zzz-DTTCompatibility-class),
14
- orientationMarkPosition,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- originalProjectName
(zzz-ARRCompatibility-class),
13
- originalProjectName,ArraySetFile-method
(zzz-ARRCompatibility-class),
13
- originRelativeTo
(zzz-DTTCompatibility-class),
14
- originRelativeTo,Array-method
(zzz-DTTCompatibility-class),
14
- otherKind (zzz-DTTCompatibility-class),
14
- otherKind,Measurement-method
(zzz-DTTCompatibility-class),
14
- owner_assnreflist
(zzz-DTTCompatibility-class),
14
- owner_assnreflist,Security-method
(zzz-DTTCompatibility-class),
14
- Owner_assnreflist-class
(zzz-DTTCompatibility-class),
14
- pages (zzz-DTTCompatibility-class), 14
- pages,BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- parameter (zzz-DTTCompatibility-class),
14
- parameter,ParameterTypes_assnlist-method
(zzz-DTTCompatibility-class),
14
- Parameter-class
(zzz-DTTCompatibility-class),
14
- parameter_ref
(zzz-DTTCompatibility-class),
14
- parameter_ref,ParameterType_assnref-method
(zzz-DTTCompatibility-class),
14
- Parameter_ref-class
(zzz-DTTCompatibility-class),
14
- parameters_assnlist
(zzz-DTTCompatibility-class),
14
- parameters_assnlist,BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- Parameters_assnlist-class
(zzz-DTTCompatibility-class),
14
- parameterType_assnref
(zzz-DTTCompatibility-class),
14
- parameterType_assnref,ParameterValue-method
(zzz-DTTCompatibility-class),
14
- ParameterType_assnref-class
(zzz-DTTCompatibility-class),
14
- parameterTypes_assnlist
(zzz-DTTCompatibility-class),
14
- parameterTypes_assnlist,Hardware-method
(zzz-DTTCompatibility-class),
14
- parameterTypes_assnlist,Protocol-method
(zzz-DTTCompatibility-class),
14
- parameterTypes_assnlist,Software-method
(zzz-DTTCompatibility-class),
14
- ParameterTypes_assnlist-class

- (zzz-DTTCompatibility-class),
14
- parameterValue
(zzz-DTTCompatibility-class),
14
- parameterValue,ParameterValues_assnlist-method
(zzz-DTTCompatibility-class),
14
- ParameterValue-class
(zzz-DTTCompatibility-class),
14
- parameterValues_assnlist
(zzz-DTTCompatibility-class),
14
- parameterValues_assnlist,HardwareApplication-Method
(zzz-DTTCompatibility-class),
14
- parameterValues_assnlist,ProtocolApplication-Method
(zzz-DTTCompatibility-class),
14
- parameterValues_assnlist,SoftwareApplication-Method
(zzz-DTTCompatibility-class),
14
- ParameterValues_assnlist-class
(zzz-DTTCompatibility-class),
14
- parent_assnref
(zzz-DTTCompatibility-class),
14
- parent_assnref,Organization-method
(zzz-DTTCompatibility-class),
14
- Parent_assnref-class
(zzz-DTTCompatibility-class),
14
- password (accessors), 3
- password,NetAffxResource-method
(NetAffxResource-class), 6
- patassignmentMethod
(zzz-ARRCompatibility-class),
13
- patassignmentMethod,PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- performer_assnref
(zzz-DTTCompatibility-class),
14
- performer_assnref,Audit-method
(zzz-DTTCompatibility-class),
14
- (zzz-DTTCompatibility-class),
14
- Performer_assnref-class
(zzz-DTTCompatibility-class),
14
- performers_assnreflist
(zzz-DTTCompatibility-class),
14
- performers_assnreflist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
14
- Performers_assnreflist-class
(zzz-DTTCompatibility-class),
14
- Performers_assnreflist-class
(zzz-DTTCompatibility-class),
14
- person_ref,ArrayManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,Contacts_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,DesignProviders_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,HardwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,Members_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,Owner_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,Performer_assnref-method
(zzz-DTTCompatibility-class),
14
- person_ref,Performers_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref,Providers_assnreflist-method
(zzz-DTTCompatibility-class),
14

- (zzz-DTTCompatibility-class),
14
- person_ref, SoftwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
14
- person_ref, SourceContact_assnreflist-method
(zzz-DTTCompatibility-class),
14
- Person_ref-class
(zzz-DTTCompatibility-class),
14
- phone (zzz-DTTCompatibility-class), 14
- phone, Organization-method
(zzz-DTTCompatibility-class),
14
- phone, Person-method
(zzz-DTTCompatibility-class),
14
- physicalArray
(zzz-ARRCompatibility-class),
13
- physicalArray, PhysicalArrays-method
(zzz-ARRCompatibility-class),
13
- PhysicalArray-class
(zzz-ARRCompatibility-class),
13
- physicalArrayDesign
(zzz-DTTCompatibility-class),
14
- physicalArrayDesign, ArrayDesign_assnlist-method
(zzz-DTTCompatibility-class),
14
- PhysicalArrayDesign-class
(zzz-DTTCompatibility-class),
14
- physicalArrayDesign_ref
(zzz-DTTCompatibility-class),
14
- physicalArrayDesign_ref, ArrayDesign_assnref-method
(zzz-DTTCompatibility-class),
14
- PhysicalArrayDesign_ref-class
(zzz-DTTCompatibility-class),
14
- physicalArrays
(zzz-ARRCompatibility-class),
13
- physicalArrays, ArraySetFile-method
(zzz-ARRCompatibility-class),
13
- PhysicalArrays-class
(zzz-ARRCompatibility-class),
13
- physicalBioAssay
(zzz-DTTCompatibility-class),
14
- physicalBioAssay, BioAssay_assnlist-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay, PhysicalBioAssay_assnlist-method
(zzz-DTTCompatibility-class),
14
- PhysicalBioAssay-class
(zzz-DTTCompatibility-class),
14
- PhysicalBioAssay_assnlist-class
(zzz-DTTCompatibility-class),
14
- PhysicalBioAssay_assnref-class
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, BioAssay_assnref-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, BioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, PhysicalBioAssay_assnref-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, PhysicalBioAssaySource_assnref-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, PhysicalBioAssayTarget_assnref-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, SourceBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
14
- physicalBioAssay_ref, Target_assnref-method
(zzz-DTTCompatibility-class),
14

- physicalBioAssay_ref,TopLevelBioAssays_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 14
- PhysicalBioAssay_ref-class
 (zzz-DTTCompatibility-class), 14
- physicalBioAssayData_assnlist
 (zzz-DTTCompatibility-class), 14
- physicalBioAssayData_assnlist,PhysicalBioAssay-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 14
- PhysicalBioAssayData_assnlist-class
 (zzz-DTTCompatibility-class), 14
- physicalBioAssaySource_assnref
 (zzz-DTTCompatibility-class), 14
- physicalBioAssaySource_assnref,FeatureExtraction-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 14
- PhysicalBioAssaySource_assnref-class
 (zzz-DTTCompatibility-class), 14
- PhysicalBioAssayTarget_assnref-class
 (zzz-DTTCompatibility-class), 14
- polymerType_assn
 (zzz-DTTCompatibility-class), 14
- polymerType_assn,BioSequence-method
 (zzz-DTTCompatibility-class), 14
- PolymerType_assn-class
 (zzz-DTTCompatibility-class), 14
- position (zzz-DTTCompatibility-class),
 14
- position,Position_assn-method
 (zzz-DTTCompatibility-class), 14
- Position-class
 (zzz-DTTCompatibility-class), 14
- position_assn
 (zzz-DTTCompatibility-class), 14
- position_assn,Feature-method
 (zzz-DTTCompatibility-class), 14
- position_assn,Fiducial-method
 (zzz-DTTCompatibility-class), 14
- Position_assn-class
 (zzz-DTTCompatibility-class), 14
- positionDelta
 (zzz-DTTCompatibility-class), 14
- positionDelta,PositionDelta_assn-method
 (zzz-DTTCompatibility-class), 14
- PositionDelta-class
 (zzz-DTTCompatibility-class), 14
- positionDelta_assn
 (zzz-DTTCompatibility-class), 14
- positionDelta_assn,FeatureDefect-method
 (zzz-DTTCompatibility-class), 14
- positionDelta_assn,ZoneDefect-method
 (zzz-DTTCompatibility-class), 14
- PositionDelta_assn-class
 (zzz-DTTCompatibility-class), 14
- presentAbsent
 (zzz-DTTCompatibility-class), 14
- presentAbsent,QuantitationType_assnlist-method
 (zzz-DTTCompatibility-class), 14
- PresentAbsent-class
 (zzz-DTTCompatibility-class), 14
- presentAbsent_ref
 (zzz-DTTCompatibility-class), 14
- presentAbsent_ref,QuantitationType_assnref-method
 (zzz-DTTCompatibility-class), 14
- presentAbsent_ref,QuantitationTypes_assnreflist-method
 (zzz-DTTCompatibility-class), 14
- presentAbsent_ref,SourcesQuantitationType_assnreflist-method

- (zzz-DTTCompatibility-class),
14
- presentAbsent_ref,TargetQuantitationType_assnpropertySets_assnlist,BioAssayMap-method
(zzz-DTTCompatibility-class),
14
- PresentAbsent_ref-class
(zzz-DTTCompatibility-class),
14
- producerTransformation_assn
(zzz-DTTCompatibility-class),
14
- producerTransformation_assn,DerivedBioAssayDatapropertySets_assnlist,BioDataCube-method
(zzz-DTTCompatibility-class),
14
- ProducerTransformation_assn-class
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Array-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ArrayManufacture-methodpropertySets_assnlist,Channel-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ArrayManufactureDeviationpropertySets_assnlist,CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Audit-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BibliographicReferencepropertySets_assnlist,CompositePosition-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioAssayCreation-methodpropertySets_assnlist,CompositeSequence-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioAssayDataCluster-methodpropertySets_assnlist,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioAssayDimension-methodpropertySets_assnlist,Compound-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioAssayMapping-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioDataTuples-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioMaterialMeasurement-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioSample-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioSequence-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,BioSource-method
(zzz-DTTCompatibility-class),
14

| | |
|--|---|
| (zzz-DTTCompatibility-class), 14 | (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,CompoundMeasurement-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FactorValue-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,ConcentrationUnit-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,Failed-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,Database-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,Feature-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,DatabaseEntry-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureDefect-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,DerivedBioAssay-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureDimension-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,DerivedBioAssayData-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureExtraction-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,DerivedSignal-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureGroup-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,Description-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureInformation-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,DesignElementMapping-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureLocation-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,DistanceUnit-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,FeatureReporterMap-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,Error-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,Fiducial-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,ExpectedValue-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,Hardware-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,Experiment-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,HardwareApplication-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,ExperimentalFactor-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,Hybridization-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,ExperimentDesign-method (zzz-DTTCompatibility-class), 14 | propertySets_assnlist,Image-method (zzz-DTTCompatibility-class), 14 |
| propertySets_assnlist,ExternalReference-method | propertySets_assnlist,ImageAcquisition-method |

- (zzz-DTTCompatibility-class),
14
- propertySets_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,MassUnit-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,MeasuredSignal-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Measurement-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,MismatchInformation-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,NameValueType-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Node-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,NodeContents-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,NodeValue-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,OntologyEntry-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Organization-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Parameter-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ParameterValue-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Person-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Position-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,PositionDelta-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,PresentAbsent-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Protocol-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,PValue-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,QuantitationTypeMapping-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,QuantityUnit-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Ratio-method

- (zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Reporter-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ReporterDimension-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ReporterPosition-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Security-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,SeqFeatureLocation-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,SequencePosition-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Software-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,TemperatureUnit-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,TimeUnit-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Transformation-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Treatment-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,VolumeUnit-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,Zone-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ZoneDefect-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ZoneGroup-method
(zzz-DTTCompatibility-class),
14
- propertySets_assnlist,ZoneLayout-method
(zzz-DTTCompatibility-class),
14
- PropertySets_assnlist-class
(zzz-DTTCompatibility-class),
14
- protocol (zzz-DTTCompatibility-class),
14
- protocol,Protocol_assnlist-method
(zzz-DTTCompatibility-class),
14
- Protocol-class
(zzz-DTTCompatibility-class),
14
- protocol_assnlist
(zzz-DTTCompatibility-class),
14
- protocol_assnlist,Protocol_package-method
(zzz-DTTCompatibility-class),
14
- Protocol_assnlist-class
(zzz-DTTCompatibility-class),
14
- protocol_assnref
(zzz-DTTCompatibility-class),
14
- protocol_assnref,ProtocolApplication-method
(zzz-DTTCompatibility-class),
14
- Protocol_assnref-class
(zzz-DTTCompatibility-class),

- 14
 protocol_package
 (zzz-DTTCompatibility-class),
 14
- protocol_package,MAGE-ML-method
 (zzz-DTTCompatibility-class),
 14
- Protocol_package-class
 (zzz-DTTCompatibility-class),
 14
- protocol_ref
 (zzz-DTTCompatibility-class),
 14
- protocol_ref,Protocol_assnref-method
 (zzz-DTTCompatibility-class),
 14
- Protocol_ref-class
 (zzz-DTTCompatibility-class),
 14
- protocolApplication
 (zzz-DTTCompatibility-class),
 14
- protocolApplication,ProtocolApplications_assnref-method
 (zzz-DTTCompatibility-class),
 14
- ProtocolApplication-class
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,ArrayDesign-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,ArrayManufacture-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,BioAssayCreation-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,BioAssayMap-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,BioAssayTreatment-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,CompositeComposition-method
 (zzz-DTTCompatibility-class),
 14
- 14
 protocolApplications_assnlist,FeatureExtraction-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,FeatureReporterMap-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,Hybridization-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,ImageAcquisition-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,PhysicalArrayDesign-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,QuantitationTypeMap-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,ReporterCompositeMap-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,Transformation-method
 (zzz-DTTCompatibility-class),
 14
- protocolApplications_assnlist,Treatment-method
 (zzz-DTTCompatibility-class),
 14
- ProtocolApplications_assnlist-class
 (zzz-DTTCompatibility-class),
 14
- providers_assnreflist
 (zzz-DTTCompatibility-class),
 14
- providers_assnreflist,Experiment-method
 (zzz-DTTCompatibility-class),
 14
- providers_assnreflist-class
 (zzz-DTTCompatibility-class),
 14
- publication
 (zzz-DTTCompatibility-class),
 14
- publication,BibliographicReference-method
 (zzz-DTTCompatibility-class),
 14
- publication,PublicationList-method
 (zzz-DTTCompatibility-class),
 14

- publisher, BibliographicReference-method
 (zzz-DTTCompatibility-class),
 14
- pvalue (zzz-DTTCompatibility-class), 14
- pvalue, ConfidenceIndicator_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- pvalue, QuantitationType_assnlist-method
 (zzz-DTTCompatibility-class),
 14
- PValue-class
 (zzz-DTTCompatibility-class),
 14
- pvalue_ref
 (zzz-DTTCompatibility-class),
 14
- pvalue_ref, ConfidenceIndicators_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- pvalue_ref, QuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 14
- pvalue_ref, QuantitationTypes_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- pvalue_ref, SourcesQuantitationType_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
- pvalue_ref, TargetQuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 14
- PValue_ref-class
 (zzz-DTTCompatibility-class),
 14
- quality (zzz-DTTCompatibility-class), 14
- quality, ManufactureLIMS-method
 (zzz-DTTCompatibility-class),
 14
- quality, ManufactureLIMSBiomaterial-method
 (zzz-DTTCompatibility-class),
 14
- qualityControlDescription_assn
 (zzz-DTTCompatibility-class),
 14
- qualityControlDescription_assn, ExperimentDesign-method
 (zzz-DTTCompatibility-class),
 14
- QualityControlDescription_assn-class
 (zzz-DTTCompatibility-class),
 14
- qualityControlStatistics_assnlist
 (zzz-DTTCompatibility-class),
 14
- qualityControlStatistics_assnlist, ArrayManufacture-method
 (zzz-DTTCompatibility-class),
 14
- qualityControlStatistics_assnlist, BioSample-method
 (zzz-DTTCompatibility-class),
 14
- qualityControlStatistics_assnlist, BioSource-method
 (zzz-DTTCompatibility-class),
 14
- qualityControlStatistics_assnlist, LabeledExtract-method
 (zzz-DTTCompatibility-class),
 14
- QualityControlStatistics_assnlist-class
 (zzz-DTTCompatibility-class),
 14
- quantitationDimension_assnref
 (zzz-DTTCompatibility-class),
 14
- quantitationDimension_assnref, NodeContents-method
 (zzz-DTTCompatibility-class),
 14
- QuantitationDimension_assnref-class
 (zzz-DTTCompatibility-class),
 14
- quantitationType_assnlist
 (zzz-DTTCompatibility-class),
 14
- quantitationType_assnlist, QuantitationType_package-method
 (zzz-DTTCompatibility-class),
 14
- QuantitationType_assnlist-class
 (zzz-DTTCompatibility-class),
 14
- quantitationType_assnref
 (zzz-DTTCompatibility-class),
 14
- quantitationType_assnref, QuantitationTypeTuple-method
 (zzz-DTTCompatibility-class),
 14
- QuantitationType_assnref-class
 (zzz-DTTCompatibility-class),
 14
- quantitationType_package

- (zzz-DTTCompatibility-class),
14
- quantitationType_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
14
- QuantitationType_package-class
(zzz-DTTCompatibility-class),
14
- QuantitationType_ref-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension,QuantitationTypeDimension_ref-class
(zzz-DTTCompatibility-class),
14
- QuantitationTypeDimension-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_assnlist
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeDimension_assnlist-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_assnref
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_assnref,DerivedBioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_assnref,MeasuredBioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeDimension_assnref-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_ref
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_ref,QuantitationTypeDimension_ref-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeDimension_ref,QuantitationTypeDimension_ref-class
(zzz-DTTCompatibility-class),
14
- (zzz-DTTCompatibility-class),
14
- QuantitationTypeDimension_ref-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeMap
(zzz-DTTCompatibility-class),
14
- quantitationTypeMap,QuantitationTypeMap_assnlist-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeMap-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeMap_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeMap_assnlist-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeMap_ref
(zzz-DTTCompatibility-class),
14
- quantitationTypeMap_ref,QuantitationTypeMaps_assnreflist-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeMap_ref-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeMapping
(zzz-DTTCompatibility-class),
14
- quantitationTypeMapping,QuantitationTypeMapping_assn-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeMapping-class
(zzz-DTTCompatibility-class),
14
- quantitationTypeMapping_assn
(zzz-DTTCompatibility-class),
14
- quantitationTypeMapping_assn,Transformation-method
(zzz-DTTCompatibility-class),
14
- QuantitationTypeMapping_assn-class
(zzz-DTTCompatibility-class),
14

(zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,DerivedSignal-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,Error-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,ExpectedValue-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,Failed-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,MeasuredSignal-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,PresentAbsent-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,PValue-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,QuantitationTypeMapping-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,Ratio-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeMaps_assnreflist,SpecializedQuantitationType-method
 (zzz-DTTCompatibility-class),
 14
 QuantitationTypeMaps_assnreflist-class
 (zzz-DTTCompatibility-class),
 14
 quantitationTypes_assnreflist
 (zzz-DTTCompatibility-class),
 14
 quantitationTypes_assnreflist,QuantitationTypeDimensions-method
 (zzz-DTTCompatibility-class),
 14
 QuantitationTypes_assnreflist-class
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeTuple
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeTuple,QuantitationTypeTuples_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeTuple-class
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeTuples_assnlist
 (zzz-DTTCompatibility-class),
 14
 quantitationTypeTuples_assnlist,DesignElementTuple-method
 (zzz-DTTCompatibility-class),
 14
 QuantitationTypeTuples_assnlist-class
 (zzz-DTTCompatibility-class),
 14
 quantityUnit
 (zzz-DTTCompatibility-class),
 14
 quantityUnit,Unit_assn-method
 (zzz-DTTCompatibility-class),
 14
 QuantityUnit-class
 (zzz-DTTCompatibility-class),
 14
 ratio (zzz-DTTCompatibility-class), 14
 ratio,QuantitationType_assnlist-method
 (zzz-DTTCompatibility-class),
 14
 Ratio-class
 (zzz-DTTCompatibility-class),
 14
 ratio_ref (zzz-DTTCompatibility-class),
 14
 ratio_ref,QuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 14
 ratio_ref,QuantitationTypes_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
 ratio_ref,QuantitationType_assnreflist-method
 (zzz-DTTCompatibility-class),
 14
 ratio_ref,TargetQuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 14
 Ratio_ref-class

- (zzz-DTTCompatibility-class),
14
- readAnnotation (NetAffxResource-class),
6
- readAnnotation, NetAffxResource, character, character-method
(NetAffxResource-class), 6
- readAnnotation, NetAffxResource, character, numeric-method
(NetAffxResource-class), 6
- readAnnotation, NetAffxResource, missing, AffxAnnotation-class
(NetAffxResource-class), 6
- readAnnotation, NetAffxResource, numeric, character-method
(NetAffxResource-class), 6
- readAnnotation, NetAffxResource, numeric, numeric-method
(NetAffxResource-class), 6
- readArr, 8
- readMage, 9
- readNetAffx, 10
- readXml, 9, 11, 12
- regions_assnlist
(zzz-DTTCompatibility-class),
14
- regions_assnlist, SeqFeature-method
(zzz-DTTCompatibility-class),
14
- Regions_assnlist-class
(zzz-DTTCompatibility-class),
14
- releaseDate
(zzz-DTTCompatibility-class),
14
- releaseDate, SoftwareApplication-method
(zzz-DTTCompatibility-class),
14
- replacedLength
(zzz-DTTCompatibility-class),
14
- replacedLength, MismatchInformation-method
(zzz-DTTCompatibility-class),
14
- replicateDescription_assn
(zzz-DTTCompatibility-class),
14
- replicateDescription_assn, ExperimentDesign-method
(zzz-DTTCompatibility-class),
14
- ReplicateDescription_assn-class
(zzz-DTTCompatibility-class),
14
- reporter (zzz-DTTCompatibility-class),
14
- reporter, DesignElement_assnlist-method
(zzz-DTTCompatibility-class),
14
- reporter, Reporter_assnlist-method
(zzz-DTTCompatibility-class),
14
- reporter_assnlist
(zzz-DTTCompatibility-class),
14
- reporter_assnlist, DesignElement_package-method
(zzz-DTTCompatibility-class),
14
- Reporter_assnlist-class
(zzz-DTTCompatibility-class),
14
- reporter_assnref
(zzz-DTTCompatibility-class),
14
- reporter_assnref, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
14
- reporter_assnref, ReporterPosition-method
(zzz-DTTCompatibility-class),
14
- Reporter_assnref-class
(zzz-DTTCompatibility-class),
14
- reporter_ref
(zzz-DTTCompatibility-class),
14
- reporter_ref, DesignElement_assnref-method
(zzz-DTTCompatibility-class),
14
- reporter_ref, Reporter_assnref-method
(zzz-DTTCompatibility-class),
14
- reporter_ref, Reporters_assnreflist-method
(zzz-DTTCompatibility-class),
14
- Reporter_ref-class
(zzz-DTTCompatibility-class),
14
- reporterCompositeMap

| | |
|---|--|
| (zzz-DTTCompatibility-class), 14 | (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap,DesignElementMap_assnlist-reported (zzz-DTTCompatibility-class), 14 | Dimension_ref (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap,ReporterCompositeMap_assnlist-reported (zzz-DTTCompatibility-class), 14 | Dimension_ref,DesignElementDimension_assnref-method (zzz-DTTCompatibility-class), 14 |
| ReporterCompositeMap-class (zzz-DTTCompatibility-class), 14 | ReporterDimension_ref-class (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap_assnlist (zzz-DTTCompatibility-class), 14 | reporterGroup (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap_assnlist,DesignElement_package-method (zzz-DTTCompatibility-class), 14 | ReporterGroup_assnlist-method (zzz-DTTCompatibility-class), 14 |
| ReporterCompositeMap_assnlist-class (zzz-DTTCompatibility-class), 14 | ReporterGroup-class (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap_ref (zzz-DTTCompatibility-class), 14 | reporterGroup_assnlist (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap_ref,DesignElementMaps_assnref-list-method (zzz-DTTCompatibility-class), 14 | ReporterGroup_assnlist,ArrayDesign_package-method (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMap_ref,ReporterCompositeMap_assnlist (zzz-DTTCompatibility-class), 14 | ReporterGroup_assnlist-class (zzz-DTTCompatibility-class), 14 |
| ReporterCompositeMap_ref-class (zzz-DTTCompatibility-class), 14 | reporterGroup_ref (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMaps_assnreflist (zzz-DTTCompatibility-class), 14 | reporterGroup_ref,ReporterGroups_assnreflist-method (zzz-DTTCompatibility-class), 14 |
| reporterCompositeMaps_assnreflist,CompositeReporterGroup (zzz-DTTCompatibility-class), 14 | ReporterGroup_ref-class (zzz-DTTCompatibility-class), 14 |
| ReporterCompositeMaps_assnreflist-class (zzz-DTTCompatibility-class), 14 | reporterGroups_assnreflist (zzz-DTTCompatibility-class), 14 |
| reporterDimension (zzz-DTTCompatibility-class), 14 | reporterGroups_assnreflist,ArrayDesign-method (zzz-DTTCompatibility-class), 14 |
| reporterDimension,DesignElementDimension_assnref-list-method (zzz-DTTCompatibility-class), 14 | ReporterGroups_assnreflist,PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 14 |
| ReporterDimension-class | ReporterGroups_assnreflist-class |

- (zzz-DTTCompatibility-class),
14
- reporterPosition
(zzz-DTTCompatibility-class),
14
- reporterPosition,Coordinate_assn-method
(zzz-DTTCompatibility-class),
14
- reporterPosition,ReporterPositionSources_assnlist-method
(zzz-DTTCompatibility-class),
14
- ReporterPosition-class
(zzz-DTTCompatibility-class),
14
- reporterPositionSources_assnlist
(zzz-DTTCompatibility-class),
14
- reporterPositionSources_assnlist,ReporterComposition-method
(zzz-DTTCompatibility-class),
14
- ReporterPositionSources_assnlist-class
(zzz-DTTCompatibility-class),
14
- reporters_assnreflist
(zzz-DTTCompatibility-class),
14
- reporters_assnreflist,ReporterDimension-method
(zzz-DTTCompatibility-class),
14
- reporters_assnreflist,ReporterGroup-method
(zzz-DTTCompatibility-class),
14
- Reporters_assnreflist-class
(zzz-DTTCompatibility-class),
14
- required(zzz-ARRCompatibility-class),
13
- required,UserAttribute-method
(zzz-ARRCompatibility-class),
13
- roles_assnlist
(zzz-DTTCompatibility-class),
14
- roles_assnlist,Organization-method
(zzz-DTTCompatibility-class),
14
- roles_assnlist,Person-method
(zzz-DTTCompatibility-class),
14
- Roles_assnlist-class
(zzz-DTTCompatibility-class),
14
- scale_assn
(zzz-DTTCompatibility-class),
14
- scale_assn,DerivedSignal-method
(zzz-DTTCompatibility-class),
14
- scale_assn,Error-method
(zzz-DTTCompatibility-class),
14
- scale_assn,ExpectedValue-method
(zzz-DTTCompatibility-class),
14
- scale_assn,Failed-method
(zzz-DTTCompatibility-class),
14
- scale_assn,MeasuredSignal-method
(zzz-DTTCompatibility-class),
14
- scale_assn,NodeValue-method
(zzz-DTTCompatibility-class),
14
- scale_assn,PresentAbsent-method
(zzz-DTTCompatibility-class),
14
- scale_assn,PValue-method
(zzz-DTTCompatibility-class),
14
- scale_assn,Ratio-method
(zzz-DTTCompatibility-class),
14
- scale_assn,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
14
- Scale_assn-class
(zzz-DTTCompatibility-class),
14
- security(zzz-DTTCompatibility-class),
14
- security,Security_assnlist-method
(zzz-DTTCompatibility-class),
14
- Security-class
(zzz-DTTCompatibility-class),
14

- [security_assnlist](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnlist, AuditAndSecurity_package-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [Security_assnlist-class](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Array-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, ArrayDesign-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, ArrayGroup-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, ArrayManufacture-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Audit-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BibliographicReference-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioAssayCreation-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioAssayDataCluster-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioAssayDimension-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioAssayMap-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioAssayTreatment-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioSample-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioSequence-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, BioSource-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Channel-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, CompositeCompositeMap-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, CompositeGroup-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, CompositeSequence-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, CompositeSequenceDimension-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Compound-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Database-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, DerivedBioAssay-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, DerivedBioAssayData-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, DerivedSignal-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Description-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Error-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, ExpectedValue-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref, Experiment-method](#)
 (zzz-DTTCompatibility-class),
[14](#)

- security_assnref, ExperimentalFactor-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, ExperimentDesign-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, FactorValue-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Failed-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Feature-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, FeatureDimension-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, FeatureExtraction-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, FeatureGroup-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Fiducial-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Hardware-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, HardwareApplication-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Hybridization-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Image-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, ImageAcquisition-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, LabeledExtract-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, MAGE-ML-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, ManufactureLIMS-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, MeasuredSignal-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Node-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, NodeContents-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Organization-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Parameter-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Person-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, PresentAbsent-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, Protocol-method
(zzz-DTTCompatibility-class),
[14](#)
- security_assnref, ProtocolApplication-method
(zzz-DTTCompatibility-class),
[14](#)

- [security_assnref,PValue-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,QuantitationTypeDimension-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,QuantitationTypeMap-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Ratio-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Reporter-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,ReporterCompositeMap-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,ReporterDimension-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,ReporterGroup-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Security-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,SecurityGroup-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,SeqFeature-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Software-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,SoftwareApplication-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,SpecializedQuantitationType-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Transformation-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Treatment-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [security_assnref,Zone-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [Security-assnref-class](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [Security_ref](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [Security_ref,Security_assnref-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [Security_ref-class](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup,SecurityGroup_assnlist-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup-class](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup_assnlist](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup_assnlist,AuditAndSecurity_package-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup_assnlist-class](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup_ref](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup_ref,SecurityGroups_assnreflist-method](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroup_ref-class](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroups_assnreflist](#)
 (zzz-DTTCompatibility-class),
[14](#)
- [SecurityGroups_assnreflist,Security-method](#)
 (zzz-DTTCompatibility-class),
[14](#)

- SecurityGroups_assnreflist-class
(zzz-DTTCompatibility-class),
14
- seqFeature
(zzz-DTTCompatibility-class),
14
- seqFeature, SeqFeatures_assnlist-method
(zzz-DTTCompatibility-class),
14
- SeqFeature-class
(zzz-DTTCompatibility-class),
14
- seqFeatureLocation
(zzz-DTTCompatibility-class),
14
- seqFeatureLocation, Regions_assnlist-method
(zzz-DTTCompatibility-class),
14
- seqFeatureLocation, Subregions_assnlist-method
(zzz-DTTCompatibility-class),
14
- SeqFeatureLocation-class
(zzz-DTTCompatibility-class),
14
- seqFeatures_assnlist
(zzz-DTTCompatibility-class),
14
- seqFeatures_assnlist, BioSequence-method
(zzz-DTTCompatibility-class),
14
- SeqFeatures_assnlist-class
(zzz-DTTCompatibility-class),
14
- sequenceDatabases_assnlist
(zzz-DTTCompatibility-class),
14
- sequenceDatabases_assnlist, BioSequence-method
(zzz-DTTCompatibility-class),
14
- SequenceDatabases_assnlist-class
(zzz-DTTCompatibility-class),
14
- sequencePosition
(zzz-DTTCompatibility-class),
14
- sequencePosition, Coordinate_assn-method
(zzz-DTTCompatibility-class),
14
- SequencePosition-class
(zzz-DTTCompatibility-class),
14
- serialNumber
(zzz-DTTCompatibility-class),
14
- serialNumber, HardwareApplication-method
(zzz-DTTCompatibility-class),
14
- show, .TypedAssociation-method
(Internal-class), 4
- show, AffyCompatible-method
(AffyCompatible-class), 3
- show, NetAffxResource-method
(NetAffxResource-class), 6
- software (zzz-DTTCompatibility-class),
14
- software, Software_assnlist-method
(zzz-DTTCompatibility-class),
14
- Software-class
(zzz-DTTCompatibility-class),
14
- software_assnlist
(zzz-DTTCompatibility-class),
14
- software_assnlist, Protocol_package-method
(zzz-DTTCompatibility-class),
14
- Software_assnlist-class
(zzz-DTTCompatibility-class),
14
- software_assnref
(zzz-DTTCompatibility-class),
14
- software_assnref, SoftwareApplication-method
(zzz-DTTCompatibility-class),
14
- Software_assnref-class
(zzz-DTTCompatibility-class),
14
- software_ref
(zzz-DTTCompatibility-class),
14
- software_ref, Software_assnref-method
(zzz-DTTCompatibility-class),
14
- software_ref, Softwares_assnreflist-method

- (zzz-DTTCompatibility-class),
14
- Software_ref-class
(zzz-DTTCompatibility-class),
14
- softwareApplication
(zzz-DTTCompatibility-class),
14
- softwareApplication,SoftwareApplications_assnreflist
(zzz-DTTCompatibility-class),
14
- SoftwareApplication-class
(zzz-DTTCompatibility-class),
14
- softwareApplications_assnreflist
(zzz-DTTCompatibility-class),
14
- softwareApplications_assnreflist,ProtocolApplications_assnreflist
(zzz-DTTCompatibility-class),
14
- SoftwareApplications_assnreflist-class
(zzz-DTTCompatibility-class),
14
- softwareManufacturers_assnreflist
(zzz-DTTCompatibility-class),
14
- softwareManufacturers_assnreflist,SoftwareManufacturers_assnreflist
(zzz-DTTCompatibility-class),
14
- SoftwareManufacturers_assnreflist-class
(zzz-DTTCompatibility-class),
14
- softwares_assnreflist
(zzz-DTTCompatibility-class),
14
- softwares_assnreflist,Hardware-method
(zzz-DTTCompatibility-class),
14
- softwares_assnreflist,Protocol-method
(zzz-DTTCompatibility-class),
14
- softwares_assnreflist,Software-method
(zzz-DTTCompatibility-class),
14
- Softwares_assnreflist-class
(zzz-DTTCompatibility-class),
14
- sourceBioAssays_assnreflist
(zzz-DTTCompatibility-class),
14
- sourceBioAssays_assnreflist,BioAssayMap-method
(zzz-DTTCompatibility-class),
14
- SourceBioAssays_assnreflist-class
(zzz-DTTCompatibility-class),
14
- sourceBioMaterialMeasurements_assnreflist
(zzz-DTTCompatibility-class),
14
- sourceBioMaterialMeasurements_assnreflist,BioAssayCreation-method
(zzz-DTTCompatibility-class),
14
- sourceBioMaterialMeasurements_assnreflist,Hybridization-method
(zzz-DTTCompatibility-class),
14
- sourceBioMaterialMeasurements_assnreflist,Treatment-method
(zzz-DTTCompatibility-class),
14
- SourceBioMaterialMeasurements_assnreflist-class
(zzz-DTTCompatibility-class),
14
- sourceContact_assnreflist
(zzz-DTTCompatibility-class),
14
- sourceContact_assnreflist,BioSource-method
(zzz-DTTCompatibility-class),
14
- SourceContact_assnreflist-class
(zzz-DTTCompatibility-class),
14
- sourcesQuantitationType_assnreflist
(zzz-DTTCompatibility-class),
14
- sourcesQuantitationType_assnreflist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
14
- SourcesQuantitationType_assnreflist-class
(zzz-DTTCompatibility-class),
14
- spacingBetweenCols
(zzz-DTTCompatibility-class),
14
- spacingBetweenCols,ZoneLayout-method
(zzz-DTTCompatibility-class),
14
- spacingBetweenRows

- (zzz-DTTCompatibility-class),
14
- spacingBetweenRows, ZoneLayout-method
(zzz-DTTCompatibility-class),
14
- spacingsBetweenZonesX
(zzz-DTTCompatibility-class),
14
- spacingsBetweenZonesX, ZoneGroup-method
(zzz-DTTCompatibility-class),
14
- spacingsBetweenZonesY
(zzz-DTTCompatibility-class),
14
- spacingsBetweenZonesY, ZoneGroup-method
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType, QuantitationType_assignment-method
(zzz-DTTCompatibility-class),
14
- SpecializedQuantitationType-class
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType_ref
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType_ref, QuantitationType_assignment-method
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType_ref, QuantitationType_assignment-method
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType_ref, SourcesQuantitationType_assignment-method
(zzz-DTTCompatibility-class),
14
- specializedQuantitationType_ref, TargetQuantitationType_assignment-method
(zzz-DTTCompatibility-class),
14
- SpecializedQuantitationType_ref-class
(zzz-DTTCompatibility-class),
14
- species_assn
(zzz-DTTCompatibility-class),
14
- species_assn, BioSequence-method
(zzz-DTTCompatibility-class),
14
- species_assn, CompositeGroup-method
(zzz-DTTCompatibility-class),
14
- species_assn, FeatureGroup-method
(zzz-DTTCompatibility-class),
14
- species_assn, ReporterGroup-method
(zzz-DTTCompatibility-class),
14
- Species_assn-class
(zzz-DTTCompatibility-class),
14
- StandardQuantitationType_ref-class
(zzz-DTTCompatibility-class),
14
- startCoord
(zzz-DTTCompatibility-class),
14
- startCoord, MatchInformation-method
(zzz-DTTCompatibility-class),
14
- strandType
(zzz-DTTCompatibility-class),
14
- strandType, SeqFeatureLocation-method
(zzz-DTTCompatibility-class),
14
- subregions_assignment-method
(zzz-DTTCompatibility-class),
14
- subregions_assignment-method
(zzz-DTTCompatibility-class),
14
- SubregionType_assignment-method
(zzz-DTTCompatibility-class),
14
- SubregionType_assignment-method
(zzz-DTTCompatibility-class),
14
- substrateType_assignment, ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- SubstrateType_assignment-class
(zzz-DTTCompatibility-class),
14
- summaryStatistics_assnlist

| | |
|---|--|
| (zzz-DTTCompatibility-class), 14 | (zzz-DTTCompatibility-class), 14 |
| summaryStatistics_assnlist,DerivedBioAssayDataTechnologyType_assn (zzz-DTTCompatibility-class), 14 | (zzz-DTTCompatibility-class), 14 |
| summaryStatistics_assnlist,MeasuredBioAssayDataTechnologyType_assn,FeatureGroup-method (zzz-DTTCompatibility-class), 14 | (zzz-DTTCompatibility-class), 14 |
| SummaryStatistics_assnlist-class (zzz-DTTCompatibility-class), 14 | TechnologyType_assn-class (zzz-DTTCompatibility-class), 14 |
| surfaceType_assn (zzz-DTTCompatibility-class), 14 | temperatureUnit (zzz-DTTCompatibility-class), 14 |
| surfaceType_assn,PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 14 | temperatureUnit,Unit_assn-method (zzz-DTTCompatibility-class), 14 |
| SurfaceType_assn-class (zzz-DTTCompatibility-class), 14 | TemperatureUnit-class (zzz-DTTCompatibility-class), 14 |
| target_assnref (zzz-DTTCompatibility-class), 14 | TemplateFile-class (zzz-ARRCompatibility-class), 13 |
| target_assnref,BioAssayTreatment-method (zzz-DTTCompatibility-class), 14 | timeUnit (zzz-DTTCompatibility-class), 14 |
| target_assnref,ImageAcquisition-method (zzz-DTTCompatibility-class), 14 | timeUnit,Unit_assn-method (zzz-DTTCompatibility-class), 14 |
| Target_assnref-class (zzz-DTTCompatibility-class), 14 | TimeUnit-class (zzz-DTTCompatibility-class), 14 |
| targetQuantitationType_assnref (zzz-DTTCompatibility-class), 14 | tolerance (zzz-DTTCompatibility-class), 14 |
| targetQuantitationType_assnref,Error-method (zzz-DTTCompatibility-class), 14 | tolerance,ArrayManufacture-method (zzz-DTTCompatibility-class), 14 |
| targetQuantitationType_assnref,ExpectedValue-method (zzz-DTTCompatibility-class), 14 | tollFreePhone (zzz-DTTCompatibility-class), 14 |
| targetQuantitationType_assnref,PValue-method (zzz-DTTCompatibility-class), 14 | tollFreePhone,Organization-method (zzz-DTTCompatibility-class), 14 |
| targetQuantitationType_assnref,QuantitationTypeMap-method (zzz-DTTCompatibility-class), 14 | tollFreePhone,Person-method (zzz-DTTCompatibility-class), 14 |
| TargetQuantitationType_assnref-class | topLevelBioAssays_assnreflist (zzz-DTTCompatibility-class), 14 |

- topLevelBioAssays_assnreflist, ExperimentDesign-method
(zzz-DTTCompatibility-class),
14
- topLevelBioAssays_assnreflist, PhysicalArray-method
(zzz-ARRCompatibility-class),
13
- TopLevelBioAssays_assnreflist-class
(zzz-DTTCompatibility-class),
14
- type, TemplateFile-method
(zzz-ARRCompatibility-class),
13
- transformation
(zzz-DTTCompatibility-class),
14
- type, UserAttribute-method
(zzz-ARRCompatibility-class),
13
- transformation, ProducerTransformation_assn-method
(zzz-DTTCompatibility-class),
14
- type_assn (zzz-DTTCompatibility-class),
14
- Transformation-class
(zzz-DTTCompatibility-class),
14
- type_assn, BioSample-method
(zzz-DTTCompatibility-class),
14
- treatment (zzz-DTTCompatibility-class),
14
- type_assn, BioSequence-method
(zzz-DTTCompatibility-class),
14
- treatment, Treatments_assnlist-method
(zzz-DTTCompatibility-class),
14
- type_assn, DatabaseEntry-method
(zzz-DTTCompatibility-class),
14
- Treatment-class
(zzz-DTTCompatibility-class),
14
- type_assn, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
14
- treatments_assnlist
(zzz-DTTCompatibility-class),
14
- type_assn, Hardware-method
(zzz-DTTCompatibility-class),
14
- treatments_assnlist, BioSample-method
(zzz-DTTCompatibility-class),
14
- type_assn, NodeValue-method
(zzz-DTTCompatibility-class),
14
- treatments_assnlist, BioSource-method
(zzz-DTTCompatibility-class),
14
- type_assn, Protocol-method
(zzz-DTTCompatibility-class),
14
- treatments_assnlist, LabeledExtract-method
(zzz-DTTCompatibility-class),
14
- type_assn, Software-method
(zzz-DTTCompatibility-class),
14
- Treatments_assnlist-class
(zzz-DTTCompatibility-class),
14
- Type_assn-class
(zzz-DTTCompatibility-class),
14
- type (zzz-DTTCompatibility-class), 14
- types_assnlist
(zzz-DTTCompatibility-class),
14
- type, ArraySetFile-method
(zzz-ARRCompatibility-class),
13
- types_assnlist, CompositeGroup-method
(zzz-DTTCompatibility-class),
14
- type, Measurement-method
(zzz-DTTCompatibility-class),
14
- types_assnlist, ExperimentDesign-method
(zzz-DTTCompatibility-class),
14
- type, NameValueType-method
(zzz-DTTCompatibility-class),
14
- types_assnlist, FeatureGroup-method

- (zzz-DTTCompatibility-class),
14
- types_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class),
14
- Types_assnlist-class
(zzz-DTTCompatibility-class),
14
- unit_assn (zzz-DTTCompatibility-class),
14
- unit_assn,Measurement-method
(zzz-DTTCompatibility-class),
14
- Unit_assn-class
(zzz-DTTCompatibility-class),
14
- unitName (zzz-DTTCompatibility-class),
14
- unitName,ConcentrationUnit-method
(zzz-DTTCompatibility-class),
14
- unitName,DistanceUnit-method
(zzz-DTTCompatibility-class),
14
- unitName,MassUnit-method
(zzz-DTTCompatibility-class),
14
- unitName,QuantityUnit-method
(zzz-DTTCompatibility-class),
14
- unitName,TemperatureUnit-method
(zzz-DTTCompatibility-class),
14
- unitName,TimeUnit-method
(zzz-DTTCompatibility-class),
14
- unitName,VolumeUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV
(zzz-DTTCompatibility-class),
14
- unitNameCV,ConcentrationUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV,DistanceUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV,MassUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV,QuantityUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV,TemperatureUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV,TimeUnit-method
(zzz-DTTCompatibility-class),
14
- unitNameCV,VolumeUnit-method
(zzz-DTTCompatibility-class),
14
- upperLeftX
(zzz-DTTCompatibility-class),
14
- upperLeftX,Zone-method
(zzz-DTTCompatibility-class),
14
- upperLeftY
(zzz-DTTCompatibility-class),
14
- upperLeftY,Zone-method
(zzz-DTTCompatibility-class),
14
- uri (zzz-DTTCompatibility-class), 14
- uri,BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- uri,Database-method
(zzz-DTTCompatibility-class),
14
- uri,DatabaseEntry-method
(zzz-DTTCompatibility-class),
14
- uri,Description-method
(zzz-DTTCompatibility-class),
14
- uri,Hardware-method
(zzz-DTTCompatibility-class),
14
- uri,Image-method
(zzz-DTTCompatibility-class),
14
- uri,Organization-method
(zzz-DTTCompatibility-class),

- 14
- uri, Person-method
(zzz-DTTCompatibility-class),
14
- uri, Protocol-method
(zzz-DTTCompatibility-class),
14
- uri, Software-method
(zzz-DTTCompatibility-class),
14
- user (accessors), 3
- user, NetAffxResource-method
(NetAffxResource-class), 6
- userAttribute
(zzz-ARRCompatibility-class),
13
- userAttribute, UserAttributes-method
(zzz-ARRCompatibility-class),
13
- UserAttribute-class
(zzz-ARRCompatibility-class),
13
- userAttributes
(zzz-ARRCompatibility-class),
13
- userAttributes, ArraySetFile-method
(zzz-ARRCompatibility-class),
13
- userAttributes, TemplateFile-method
(zzz-ARRCompatibility-class),
13
- UserAttributes-class
(zzz-ARRCompatibility-class),
13
- userAttributeValue
(zzz-ARRCompatibility-class),
13
- userAttributeValue, UserAttribute-method
(zzz-ARRCompatibility-class),
13
- UserAttributeValue-class
(zzz-ARRCompatibility-class),
13
- value (zzz-ARRCompatibility-class), 13
- value, Control-method
(zzz-ARRCompatibility-class),
13
- value, Datum-method
(zzz-DTTCompatibility-class),
14
- value, Measurement-method
(zzz-DTTCompatibility-class),
14
- value, NameValueType-method
(zzz-DTTCompatibility-class),
14
- value, NodeValue-method
(zzz-DTTCompatibility-class),
14
- value, OntologyEntry-method
(zzz-DTTCompatibility-class),
14
- value, ParameterValue-method
(zzz-DTTCompatibility-class),
14
- value_assn
(zzz-DTTCompatibility-class),
14
- value_assn, FactorValue-method
(zzz-DTTCompatibility-class),
14
- Value_assn-class
(zzz-DTTCompatibility-class),
14
- vector, 4
- version (zzz-ARRCompatibility-class), 13
- version, ArrayDesign-method
(zzz-DTTCompatibility-class),
14
- version, ArraySetFile-method
(zzz-ARRCompatibility-class),
13
- version, Database-method
(zzz-DTTCompatibility-class),
14
- version, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- version, SoftwareApplication-method
(zzz-DTTCompatibility-class),
14
- version, TemplateFile-method
(zzz-ARRCompatibility-class),
13
- volume (zzz-DTTCompatibility-class), 14
- volume, BibliographicReference-method

- (zzz-DTTCompatibility-class),
14
- volumeUnit
(zzz-DTTCompatibility-class),
14
- volumeUnit,Unit_assn-method
(zzz-DTTCompatibility-class),
14
- VolumeUnit-class
(zzz-DTTCompatibility-class),
14
- warningType_assn
(zzz-DTTCompatibility-class),
14
- warningType_assn,Reporter-method
(zzz-DTTCompatibility-class),
14
- WarningType_assn-class
(zzz-DTTCompatibility-class),
14
- width (zzz-DTTCompatibility-class), 14
- width,ArrayGroup-method
(zzz-DTTCompatibility-class),
14
- x (zzz-DTTCompatibility-class), 14
- x,Position-method
(zzz-DTTCompatibility-class),
14
- xclass, 9, 11, 12
- xmlValue.XMLAttributeValue, 13
- y (zzz-DTTCompatibility-class), 14
- y,Position-method
(zzz-DTTCompatibility-class),
14
- year (zzz-DTTCompatibility-class), 14
- year,BibliographicReference-method
(zzz-DTTCompatibility-class),
14
- zone (zzz-DTTCompatibility-class), 14
- zone,ZoneLocations_assnlist-method
(zzz-DTTCompatibility-class),
14
- Zone-class
(zzz-DTTCompatibility-class),
14
- zone_assnref
(zzz-DTTCompatibility-class),
14
- zone_assnref,Feature-method
(zzz-DTTCompatibility-class),
14
- zone_assnref,ZoneDefect-method
(zzz-DTTCompatibility-class),
14
- Zone_assnref-class
(zzz-DTTCompatibility-class),
14
- zone_ref (zzz-DTTCompatibility-class),
14
- zone_ref,Zone_assnref-method
(zzz-DTTCompatibility-class),
14
- Zone_ref-class
(zzz-DTTCompatibility-class),
14
- zoneDefect
(zzz-DTTCompatibility-class),
14
- zoneDefect,Adjustments_assnlist-method
(zzz-DTTCompatibility-class),
14
- ZoneDefect-class
(zzz-DTTCompatibility-class),
14
- zoneGroup (zzz-DTTCompatibility-class),
14
- zoneGroup,ZoneGroups_assnlist-method
(zzz-DTTCompatibility-class),
14
- ZoneGroup-class
(zzz-DTTCompatibility-class),
14
- zoneGroups_assnlist
(zzz-DTTCompatibility-class),
14
- zoneGroups_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
14
- ZoneGroups_assnlist-class
(zzz-DTTCompatibility-class),
14
- zoneLayout
(zzz-DTTCompatibility-class),

14
zoneLayout, ZoneLayout_assn-method
(zzz-DTTCompatibility-class),
14
ZoneLayout-class
(zzz-DTTCompatibility-class),
14
zoneLayout_assn
(zzz-DTTCompatibility-class),
14
zoneLayout_assn, ZoneGroup-method
(zzz-DTTCompatibility-class),
14
ZoneLayout_assn-class
(zzz-DTTCompatibility-class),
14
zoneLocations_assnlist
(zzz-DTTCompatibility-class),
14
zoneLocations_assnlist, ZoneGroup-method
(zzz-DTTCompatibility-class),
14
ZoneLocations_assnlist-class
(zzz-DTTCompatibility-class),
14
zonesPerX (zzz-DTTCompatibility-class),
14
zonesPerX, ZoneGroup-method
(zzz-DTTCompatibility-class),
14
zonesPerY (zzz-DTTCompatibility-class),
14
zonesPerY, ZoneGroup-method
(zzz-DTTCompatibility-class),
14
zzz-.TypedAssociation-class
(Internal-class), 4
zzz-AffyCompatible-class
(AffyCompatible-class), 3
zzz-ARRCompatibility-class, 13
zzz-DTTCompatibility-class, 14
zzz-NetAffxCompatibility-class, 16