Package 'seqPattern'

May 25, 2024

Title Visualising oligonucleotide patterns and motif occurrences across a set of sorted sequences

Version 1.37.0

Date 08-05-2015

Author Vanja Haberle <vanja.haberle@gmail.com>

Maintainer Vanja Haberle <vanja.haberle@gmail.com>

Imports Biostrings, GenomicRanges, IRanges, KernSmooth, plotrix

Depends methods, R (>= 2.15.0)

Suggests BSgenome.Drerio.UCSC.danRer7, CAGEr, RUnit, BiocGenerics, BiocStyle

Enhances parallel

Description Visualising oligonucleotide patterns and sequence motifs occurrences across a large set of sequences centred at a common reference point and sorted by a user defined feature.

License GPL-3

biocViews Visualization, SequenceMatching

Collate MotifScanningFunctions.R MotifScanningMethods.R PatternOccurrenceFunctions.R PatternOccurrenceMethods.R PlottingFunctions.R PlottingMethods.R

git_url https://git.bioconductor.org/packages/seqPattern

git_branch devel

git_last_commit 03617d6

git_last_commit_date 2024-04-30

Repository Bioconductor 3.20

Date/Publication 2024-05-24

Contents

seqPattern-package		 2
getPatternOccurrenceList	•••	 3

seqPattern-package

motifScanHits	4
motifScanScores	6
plotMotifDensityMap	7
plotMotifOccurrenceAverage	10
plotMotifScanScores	12
plotPatternDensityMap	14
plotPatternOccurrenceAverage	17
TBPpwm	19
zebrafishPromoters	19
zebrafishPromoters24h	20
	21

Index

seqPattern-package

Visualising oligonucleotide patterns and motif occurences across a set of ordered sequences

Description

Visualising oligonucleotide patterns and sequence motifs occurrences across a large set of sequences centred at a common reference point and sorted by a user defined feature.

Details

Package:	SeqPattern
Type:	Package
Version:	1.0
Date:	2014-12-05
License:	GPL-3
Depends:	R (>= 3.0.1), methods

Author(s)

Vanja Haberle

Maintainer: Vanja Haberle <vanja.haberle@gmail.com>

getPatternOccurrenceList

Occurrence of sequence patterns in a set of ordered sequences

Description

Finds positions of specified sequence patterns in a list of sequences of the same length ordered by a provided index. Sequence patterns can be consensus sequences of variable length and can contain IUPAC ambiguity code. Position of each pattern occurrence is specified in two-dimensional matrix, *i.e.* the first coordinate provides the ordinal number of the sequence and the second coordinate gives the position within the sequence where the pattern occurs.

Usage

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length in which to search for the patterns.
patterns	Character vector specifying one or more DNA sequence patterns (oligonucleotides). IUPAC ambiguity codes can be used and will match any letter in the subject that is associated with the code.
seqOrder	Integer vector specifying the order of the provided input sequences. Must have the same length as the number of sequences in the regionSeq. The default value will order the sequences as they are ordered in the regionSeq object.
useMulticore	Logical, should multicore be used. useMulticore = TRUE is supported only on Unix-like platforms.
nrCores	Number of cores to use when useMulticore = TRUE. Default value NULL uses all detected cores.

Details

This function uses the matchPattern function to find occurrences of given sequence patterns in a set of input sequences. Input sequences must all be of the same length and are ordered according to the index provided in the seqOrder argument, creating a n * m matrix, where n is the number of sequences and m is the length of the sequences. Positions of pattern matches in the resulting matrix are returned as two-dimensional coordinates.

Value

The function returns a named list with one element for each sequence pattern specified in the patterns argument. Each element of the list is a data.frame with positions of the corresponding pattern in the set of input sequences. The input sequences of the same length are sorted according

to the index in seqOrder argument and the positions of pattern matches in the resulting n * m matrix (where n is the number of sequences and m is the length of the sequence) are provided. The sequence column in the data.frame provides the ordinal number of the sequence in the ordered list of sequences and the position column provides the start position of the pattern match within that sequence.

Author(s)

Vanja Haberle

See Also

```
plotPatternDensityMap
motifScanHits
```

Examples

```
library(GenomicRanges)
load(system.file("data", "zebrafishPromoters.RData", package="seqPattern"))
promoterWidth <- elementMetadata(zebrafishPromoters)$interquantileWidth</pre>
```

motifScanHits Occurrence of motifs in a set of ordered sequences

Description

Finds positions of sequence motif hits above a specified threshold in a list of sequences of the same length ordered by a provided index. Motif is specified by a position weight matrix (PWM) that contains estimated probability of base b at position i and is usually constructed via call to PWM function. Position of each motif hit is specified in two-dimensional matrix, *i.e.* the first coordinate provides the ordinal number of the sequence and the second coordinate gives the position within the sequence where the motif occurs.

Usage

```
motifScanHits(regionsSeq, motifPWM, minScore = "80%",
      seqOrder = c(1:length(regionsSeq)))
```

4

motifScanHits

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length in which to search for the motif hits.
motifPWM	A numeric matrix representing the Position Weight Matrix (PWM), such as re- turned by PWM function. Can contain either probabilities or log2 probability ratio of base b at position i.
minScore	The minimum score for counting a motif hit. Can be given as a character string containing a percentage (<i>e.g.</i> "85%") of the PWM score or a single number specifying score threshold. If a percentage is given, it is converted to a score value taking into account both minimal and maximal possible PWM scores as follows: minPWMscore + percThreshold/100 * (maxPWMscore - minPWMscore) This differs from the formula in the matchPWM function from the Biostrings package which takes into account only the maximal possible PWM score and considers the given percentage as the percentage of that maximal score: percThreshold/100 * maxPWMscore
seqOrder	Integer vector specifying the order of the provided input sequences. Must have the same length as the number of sequences in the regionSeq. The default value will order the sequences as they are ordered in the input regionSeq object.

Details

This function uses the matchPWM function to find matches to given motif in a set of input sequences. Only matches above specified minScore are considered as hits. Input sequences must all be of the same length and are ordered according to the index provided in the seqOrder argument, creating a n * m matrix, where n is the number of sequences and m is the length of the sequences. Positions of motif hits in the resulting matrix are returned as two-dimensional coordinates.

Value

The function returns a data.frame with positions of the motif hits in the set of input sequences. The input sequences of the same length are sorted according to the index in seqOrder argument and the positions of motif hits in the resulting n * m matrix (where n is the number of sequences and m is the length of the sequence) are provided. The sequence column in the data.frame provides the ordinal number of the sequence in the ordered list of sequences and the position column provides the start position of the motif hit within that sequence.

Author(s)

Vanja Haberle

See Also

plotMotifDensityMap
getPatternOccurrenceList

Examples

motifScanScores

Motif scanning scores for a set of ordered sequences

Description

Provides motif scanning scores along the full length of a sequence for a list of sequences of the same length ordered by a provided index. Motif is specified by a position weight matrix (PWM) that contains estimated probability of base b at position i and is usually constructed via call to PWM function. Scanning scores are returned in the form of a two-dimensional matrix, where the rows are sequences ordered by the specified index and the columns are relative positions within the sequence. Each cell in the matrix contains the score of the specified motif in the given sequence starting at the given position. The resulting matrix can be used to visualise motif occurrences and their strength in an ordered set of sequences centered at a common reference point.

Usage

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length to be scanned with the motif.
motifPWM	A numeric matrix representing the Position Weight Matrix (PWM), such as re- turned by PWM function. Can contain either probabilities or log2 probability ratio of base b at position i.
seqOrder	Integer vector specifying the order of the provided input sequences. Must have the same length as the number of sequences in the regionSeq. The default value will order the sequences as they are ordered in the input regionSeq object.
asPercentage	Logical, should the scores represent percentage of the maximal motif PWM score (TRUE) or raw scores (FALSE).

6

Details

This function uses the PWMscoreStartingAt function to get scores for a given motif starting at each position (nucleotide) in a set of input sequences. Input sequences must all be of the same length and are ordered according to the index provided in the seqOrder argument, creating an $n \star m$ matrix, where n is the number of sequences and m is the length of the sequences. Each cell in the matrix contains the score of the specified motif in the given sequence starting at the given position.

Value

The function returns a matrix with motif scanning scores for each position in the set of input sequences.

Author(s)

Vanja Haberle

See Also

plotMotifScanScores
motifScanHits

Examples

plotMotifDensityMap Plotting density maps of motif occurrence

Description

Plots density of motif occurrences in an ordered set of sequences of the same length in the form of a two dimensional map centered at a common reference position. Motif is specified by a position weight matrix (PWM) that contains estimated probability of base b at position i, and only motif hits above specified threshold are taken into account and plotted.

Usage

```
plotMotifDensityMap(regionsSeq, motifPWM, minScore = "80%",
    seqOrder = c(1:length(regionsSeq)), flankUp = NULL, flankDown = NULL,
    nBin = NULL, bandWidth = NULL, color = "blue", transf = NULL, xTicks = NULL,
    xTicksAt = NULL, xLabel = "", yTicks = NULL, yTicksAt = NULL, yLabel = "",
    cexAxis = 8, plotScale = TRUE, scaleLength = NULL, scaleWidth = 15,
    addReferenceLine = TRUE, plotColorLegend = TRUE, outFile = "DensityMap",
    plotWidth = 2000, plotHeight = 2000)
```

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length for which the motif occurrence density should be visualised.
motifPWM	A numeric matrix representing the Position Weight Matrix (PWM), such as re- turned by PWM function. Can contain either probabilities or log2 probability ratio of base b at position i.
minScore	The minimum score for counting a motif hit. Can be given as a character string containing a percentage (<i>e.g.</i> "85%") of the PWM score or a single number specifying score threshold. If a percentage is given, it is converted to a score value taking into account both minimal and maximal possible PWM scores as follows: minPWMscore + percThreshold/100 * (maxPWMscore - minPWMscore) This differs from the formula in the matchPWM function from the Biostrings package which takes into account only the maximal possible PWM score and considers the given percentage as the percentage of that maximal score: percThreshold/100 * maxPWMscore
seqOrder	Integer vector specifying the order of the provided input sequences. Must have the same length as the number of sequences in the regionSeq. Input sequences will be sorted according to this index in an ascending order form top to the bottom of the plot, <i>i.e.</i> the sequence labeled with the lowest number will appear at the top of the plot. The default value will order the sequences as they are ordered in the input regionSeq object.
flankUp, flankDo	
	The number of base-pairs upstream and downstream of the reference position in the provided sequences, respectively. flankUp + flankDown must sum up to the length of the sequences. If no values are provided both flankUp and flankDown are set to be half of the length of the input sequences, <i>i.e.</i> the reference position is assumed to be in the middle of the sequences.
nBin	Numeric vector with two values containing the number of equally spaced points in each direction over which the density is to be estimated. The first value spec- ifies number of bins along x-axis, <i>i.e.</i> along the nucleotides in the sequence, and the second value specifies the number of bins along y-axis, <i>i.e.</i> across ordered input sequences. The values are passed on to the gridsize argument of the bkde2D function to compute a 2D binned kernel density estimate. If nBin is not specified it will default to $c(n, m)$, where n is the number of input sequences and m is the length of sequences.
bandWidth	Numeric vector of length 2, containing the bandwidth to be used in each coor- dinate direction. The first value specifies the bandwidth along the x-axis, <i>i.e.</i>

	along the nucleotides in the sequence, and the second value specifies the band- width along u axis, i.e. agrees ordered input sequences. The values are passed	
	width along y-axis, <i>i.e.</i> across ordered input sequences. The values are passed on to the bandwidth argument of the bkde2D function to compute a 2D binned	
	kernel density estimate and are used as standard deviation of the bivariate Gaus-	
	sian kernel. If bandWidth is not specified it will default to $c(3,3)$.	
color	Character specifying the color palette for the density plot. One of the following	
	color palettes can be specified: "blue", "brown", "cyan", "gold", "gray",	
	"green", "pink", "purple", "red". Please refer to the vignette for the ap-	
transf	pearance of these palettes. The function mapping the density scale to the color scale. See Details.	
xTicks		
XTICKS	Character vector of labels to be placed at the tick-marks on x-axis. The default NULL value produces five tick-marks: one at the reference point and two equally	
	spaced tick-marks both upstream and downstream of the reference point.	
xTicksAt	Numeric vector of positions of the tick-marks on the x-axis. The values can	
	range from 1 (the position of the first base-pair in the sequence) to input se-	
	quence length. The default NULL value produces five tick-marks: one at the ref-	
	erence point and two equally spaced tick-marks both upstream and downstream of the reference point.	
xLabel	The label for the x-axis. The default is no label, <i>i.e.</i> empty string.	
yTicks	Character vector of labels to be placed at the tick-marks on y-axis. The default	
YTICKS	NULL value produces no tick-marks and labels.	
yTicksAt	Numeric vector of positions of the tick-marks on the y-axis. The values can	
	range from 1 (the position of the last sequence on the bottom of the plot) to	
	input sequence length (the position of the first sequence on the top of the plot). The default NULL value produces no tick-marks.	
yLabel	The label for the y-axis. The default is no label, <i>i.e.</i> empty string.	
cexAxis	The magnification to be used for axis annotation.	
plotScale	Logical, should the scale bar be plotted in the lower left corner of the plot.	
scaleLength	The length of the scale bar to be plotted. Used only when plotScale = TRUE. If	
	no value is provided, it defaults to one fifth of the input sequence length.	
scaleWidth	The width of the line for the scale bar. Used only when plotScale = TRUE.	
addReferenceLi		
Logical, should the vertical dashed line be drawn at the reference point. plotColorLegend		
p=0000101 =080	Logical, should the color legend for the pattern density be plotted. If TRUE a	
	separate .png file named outFile."ColorLegend.png" will be created, showing	
	mapping of pattern density values to colours.	
outFile	Character vector specifying the base name of the output plot file. The final name	
plotWidth,plot	of the plot file for each pattern will be outFile."pattern.jpg".	
μιστωτατή, μιστ	Width and height of the density plot in pixels.	
	here and the grant of the constraint processing procesing processing processing processing processing processi	

Value

The function produces a PNG file in the working directory, visualising density of the motif occurrence above specified threshold in the set of ordered input sequences.

Author(s)

Vanja Haberle

References

Haberle *et al.* (2014) Two independent transcription initiation codes overlap on vertebrate core promoters, *Nature* **507**:381-385.

See Also

motifScanHits
plotPatternDensityMap

Examples

plotMotifOccurrenceAverage

Plotting average profile of motif occurrence

Description

Plots average profile of motif occurrence for a set of input sequences of the same length. Motif is specified by a position weight matrix (PWM) that contains estimated probability of base b at position i, and only motif hits above specified threshold are taken into account.

Usage

```
plotMotifOccurrenceAverage(regionsSeq, motifPWM, minScore = "80%", flankUp =
    NULL, flankDown = NULL, smoothingWindow = 1, color = "black", xLabel =
    "Distance to reference point (bp)", yLabel = "Relative frequency", cexAxis =
    1, addReferenceLine = TRUE, plotLegend = FALSE, cexLegend = 1, add = FALSE,
    ...)
```

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length for which the patterns occurrence profile should be visualised.	
motifPWM	A numeric matrix representing the Position Weight Matrix (PWM), such as re- turned by PWM function. Can contain either probabilities or log2 probability ratio of base b at position i.	
minScore	The minimum score for counting a motif hit. Can be given as a character string containing a percentage (<i>e.g.</i> "85%") of the PWM score or a single number spec- ifying score threshold. If a percentage is given, it is converted to a score value taking into account both minimal and maximal possible PWM scores as fol- lows: minPWMscore + percThreshold/100 * (maxPWMscore - minPWMscore) This differs from the formula in the matchPWM function from the Biostrings package which takes into account only the maximal possible PWM score and considers the given percentage as the percentage of that maximal score: percThreshold/100 * maxPWMscore	
flankUp,flankD		
	The number of base-pairs upstream and downstream of the reference position in the provided sequences, respectively. flankUp + flankDown must sum up to the length of the sequences. If no values are provided both flankUp and flankDown are set to be half of the length of the input sequences, <i>i.e.</i> the reference position is assumed to be in the middle of the sequences.	
smoothingWindo	W	
	Integer specifying the size of a window (in base-pairs) to be used for smoothing the signal. Default value of 1 corresponds to no smoothing.	
color	Value specifying the color for plotting.	
xLabel, yLabel	Character strings for x and y axis labels, respectively.	
cexAxis	The magnification to be used for axis annotation relative to the current setting of cex.	
addReferenceLine		
	Logical, should the vertical dashed line be drawn at the reference point.	
plotLegend	Logical, should the legend be plotted at the top.	
cexLegend	The magnification to be used for legend relative to the current setting of cex.	
add	Logical, should the pattern occurrence profiles be added to the existing plot.	
	Further arguments to be passed to plot() or lines() methods, such as lty, <i>etc</i> .	

Value

The function finds all hits matching the motif above the specified score threshold in the set of input sequences and plots an average profile reflecting the occurrence of the motif across input sequences.

Author(s)

Vanja Haberle

See Also

motifScanHits
plotMotifDensityMap

Examples

```
smoothingWindow = 3)
```

plotMotifScanScores Plotting heatmap of motif scanning scores

Description

Plots heatmap of motif scanning scores for a set of sorted sequences of the same length in the form of a two dimensional map centered at a common reference position. Motif is specified by a position weight matrix (PWM) that contains estimated probability of base b at position i, and the percentage of the maximal PWM score is plotted for each position (nucleotide) in each sequence.

Usage

```
plotMotifScanScores(regionsSeq, motifPWM, seqOrder = c(1:length(regionsSeq)),
  flankUp = NULL, flankDown = NULL, xTicks = NULL, xTicksAt = NULL,
  xLabel = "", yTicks = NULL, yTicksAt = NULL, yLabel = "", cexAxis = 8,
  plotScale = TRUE, scaleLength = NULL, scaleWidth = 15,
  addReferenceLine = TRUE, plotColorLegend = TRUE, outFile =
  "MotifScanningScores.png", plotWidth = 2000, plotHeight = 2000)
```

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length for which the motif occurrence density should be visualised.
motifPWM	A numeric matrix representing the Position Weight Matrix (PWM), such as re- turned by PWM function. Can contain either probabilities or log2 probability ratio of base b at position i.
seqOrder	Integer vector specifying the order of the provided input sequences. Must have the same length as the number of sequences in the regionSeq. Input sequences will be sorted according to this index in an ascending order form top to the bottom of the plot, <i>i.e.</i> the sequence labeled with the lowest number will appear at the top of the plot. The default value will order the sequences as they are ordered in the input regionSeq object.

12

flankUp, flankDown

The number of base-pairs upstream and downstream of the reference position in the provided sequences, respectively. flankUp + flankDown must sum up to the length of the sequences.

- xTicks Character vector of labels to be placed at the tick-marks on x-axis. The default NULL value produces five tick-marks: one at the reference point and two equally spaced tick-marks both upstream and downstream of the reference point.
- xTicksAt Numeric vector of positions of the tick-marks on the x-axis. The values can range from 1 (the position of the first base-pair in the sequence) to input sequence length. The default NULL value produces five tick-marks: one at the reference point and two equally spaced tick-marks both upstream and downstream of the reference point.
- xLabel The label for the x-axis. The default is no label, *i.e.* empty string.
- yTicks Character vector of labels to be placed at the tick-marks on y-axis. The default NULL value produces no tick-marks and labels.
- yTicksAt Numeric vector of positions of the tick-marks on the y-axis. The values can range from 1 (the position of the last sequence on the bottom of the plot) to input sequence length (the position of the first sequence on the top of the plot). The default NULL value produces no tick-marks.
- yLabel The label for the y-axis. The default is no label, *i.e.* empty string.
- cexAxis The magnification to be used for axis annotation.
- plotScale Logical, should the scale bar be plotted in the lower left corner of the plot.
- scaleLength The length of the scale bar to be plotted. Used only when plotScale = TRUE. If no value is provided, it defaults to one fifth of the input sequence length.
- scaleWidth The width of the line for the scale bar. Used only when plotScale = TRUE. addReferenceLine

Logical, should the vertical dashed line be drawn at the reference point.

plotColorLegend

Logical, should the color legend for the scanning score be plotted on the right side of the plot.

outFile Character vector specifying the base name of the output plot file. The final name of the plot file for each pattern will be outFile."pattern.jpg".

Width and height of the density plot in pixels.

Value

The function produces a PNG file in the working directory, visualising motif scanning scores in the set of ordered input sequences.

Author(s)

Vanja Haberle

plotWidth,plotHeight

See Also

motifScanScores
plotPatternDensityMap

Examples

plotPatternDensityMap Plotting density maps of sequence pattern occurrence

Description

Plots density of sequence pattern occurrences in an ordered set of sequences of the same length in the form of a two dimensional map centered at a common reference position. Multiple sequence patterns can be processed at once and one plot per pattern will be created with the same color scale across all plots, allowing visual density comparison across different patterns.

Usage

```
plotPatternDensityMap(regionsSeq, patterns, seqOrder = c(1:length(regionsSeq)),
  flankUp = NULL, flankDown = NULL, nBin = NULL, bandWidth = NULL,
  color = "blue", transf = NULL, xTicks = NULL, xTicksAt = NULL, xLabel = "",
  yTicks = NULL, yTicksAt = NULL, yLabel = "", cexAxis = 8, plotScale = TRUE,
  scaleLength = NULL, scaleWidth = 15, addPatternLabel = TRUE, cexLabel = 8,
  labelCol = "black", addReferenceLine = TRUE, plotColorLegend = TRUE,
  outFile = "PatternDensityMap", plotWidth = 2000, plotHeight = 2000,
  useMulticore = FALSE, nrCores = NULL)
```

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length for which the patterns occurrence density should be visualised.
patterns	Character vector specifying one or more DNA sequence patterns (oligonucleotides). IUPAC ambiguity codes can be used and will match any letter in the subject that is associated with the code.
seqOrder	Integer vector specifying the order of the provided input sequences. Must have the same length as the number of sequences in the regionSeq. Input sequences will be sorted according to this index in an ascending order form top to the bottom of the plot, <i>i.e.</i> the sequence labeled with the lowest number will appear

14

at the top of the plot. The default value will order the sequences as they are ordered in the input regionSeq object.

flankUp,flankDown

The number of base-pairs upstream and downstream of the reference position in the provided sequences, respectively. flankUp + flankDown must sum up to the length of the sequences. If no values are provided both flankUp and flankDown are set to be half of the length of the input sequences, *i.e.* the reference position is assumed to be in the middle of the sequences.

- nBin Numeric vector with two values containing the number of equally spaced points in each direction over which the density is to be estimated. The first value specifies number of bins along x-axis, i.e. along the nucleotides in the sequence, and the second value specifies the number of bins along y-axis, i.e. across ordered input sequences. The values are passed on to the gridsize argument of the bkde2D function to compute a 2D binned kernel density estimate. If nBin is not specified it will default to c(n, m), where n is the number of input sequences and m is the length of sequences.
- bandWidth Numeric vector of length 2, containing the bandwidth to be used in each coordinate direction. The first value specifies the bandwidth along the x-axis, *i.e.* along the nucleotides in the sequence, and the second value specifies the bandwidth along y-axis, *i.e.* across ordered input sequences. The values are passed on to the bandwidth argument of the bkde2D function to compute a 2D binned kernel density estimate and are used as standard deviation of the bivariate Gaussian kernel. If bandWidth is not specified it will default to c(3,3).
- color Character specifying the color palette for the density plot. One of the following color palettes can be specified: "blue", "brown", "cyan", "gold", "gray", "green", "pink", "purple", "red". Please refer to the vignette for the appearance of these palettes.

transf The function mapping the density scale to the color scale. See Details.

- xTicks Character vector of labels to be placed at the tick-marks on x-axis. The default NULL value produces five tick-marks: one at the reference point and two equally spaced tick-marks both upstream and downstream of the reference point.
- xTicksAt Numeric vector of positions of the tick-marks on the x-axis. The values can range from 1 (the position of the first base-pair in the sequence) to input sequence length. The default NULL value produces five tick-marks: one at the reference point and two equally spaced tick-marks both upstream and downstream of the reference point.
- xLabel The label for the x-axis. The default is no label, *i.e.* empty string.
- yTicks Character vector of labels to be placed at the tick-marks on y-axis. The default NULL value produces no tick-marks and labels.
- yTicksAt Numeric vector of positions of the tick-marks on the y-axis. The values can range from 1 (the position of the last sequence on the bottom of the plot) to input sequence length (the position of the first sequence on the top of the plot). The default NULL value produces no tick-marks.
- yLabel The label for the y-axis. The default is no label, *i.e.* empty string.
- cexAxis The magnification to be used for axis annotation.

plotScale	Logical, should the scale bar be plotted in the lower left corner of the plot.	
scaleLength	The length of the scale bar to be plotted. Used only when plotScale = TRUE. If no value is provided, it defaults to one fifth of the input sequence length.	
scaleWidth	The width of the line for the scale bar. Used only when plotScale = TRUE.	
addPatternLabe	1	
	Logical, should the pattern label be written in the upper left corner of the plot.	
cexLabel	The magnification to be used for pattern label.	
labelCol	The color to be used for pattern label and scale bar.	
addReferenceLi	ne	
	Logical, should the vertical dashed line be drawn at the reference point.	
plotColorLegen	d	
	Logical, should the color legend for the pattern density be plotted. If TRUE a separate .png file named outFile."ColorLegend.png" will be created, showing mapping of pattern density values to colours.	
outFile	Character vector specifying the base name of the output plot file. The final name of the plot file for each pattern will be outFile."pattern.png".	
plotWidth, plotHeight		
	Width and height of the density plot(s) in pixels.	
useMulticore	Logical, should multicore be used. useMulticore = TRUE is supported only on Unix-like platforms.	
nrCores	Number of cores to use when useMulticore = TRUE. Default value NULL uses all detected cores.	

Value

The function produces PNG files in the working directory, visualising density of patterns occurrence in the set of ordered input sequences. One file/plot per specified pattern is created.

Author(s)

Vanja Haberle

References

Haberle *et al.* (2014) Two independent transcription initiation codes overlap on vertebrate core promoters, *Nature* **507**:381-385.

See Also

getPatternOccurrenceList
plotMotifDensityMap

plotPatternOccurrenceAverage

Examples

```
library(GenomicRanges)
load(system.file("data", "zebrafishPromoters.RData", package="seqPattern"))
promoterWidth <- elementMetadata(zebrafishPromoters)$interquantileWidth
# dinucleotide patterns
plotPatternDensityMap(regionsSeq = zebrafishPromoters, patterns = c("TA", "GC"),
            seqOrder = order(promoterWidth), flankUp = 400, flankDown = 600,
            color = "blue")
# motif consensus sequence
plotPatternDensityMap(regionsSeq = zebrafishPromoters, patterns = "TATAWAWR",
            seqOrder = order(promoterWidth), flankUp = 400, flankDown = 600,
            color = "cyan")</pre>
```

plotPatternOccurrenceAverage

Plotting average profile of sequence pattern occurrence

Description

Plots average profile of sequence pattern occurrence for a set of input sequences of the same length. Multiple sequence patterns can be processed at once and visualised in the same plot, allowing comparison across different patterns.

Usage

```
plotPatternOccurrenceAverage(regionsSeq, patterns, flankUp = NULL,
    flankDown = NULL, smoothingWindow = 1, color = rainbow(length(patterns)),
    xLabel = "Distance to reference point (bp)", yLabel = "Relative frequency",
    cexAxis = 1, addReferenceLine = TRUE, plotLegend = TRUE, cexLegend = 1,
    useMulticore = FALSE, nrCores = NULL, add = FALSE, ...)
```

Arguments

regionsSeq	A DNAStringSet object. Set of sequences of the same length for which the patterns occurrence profile should be visualised.	
patterns	Character vector specifying one or more DNA sequence patterns (oligonucleotides). IUPAC ambiguity codes can be used and will match any letter in the subject that is associated with the code.	
flankUp, flankDown		
	The number of base-pairs upstream and downstream of the reference position in the provided sequences, respectively. flankUp + flankDown must sum up to the length of the sequences. If no values are provided both flankUp and flankDown are set to be half of the length of the input sequences, <i>i.e.</i> the reference position	

is assumed to be in the middle of the sequences.

smoothingWindow	
	Integer specifying the size of a window (in base-pairs) to be used for smoothing the signal. Default value of 1 corresponds to no smoothing.
color	A vector of values specifying the colors for plotting. Number of colors must match the number of patterns that should be plotted.
xLabel, yLabel	Character strings for x and y axis labels, respectively.
cexAxis	The magnification to be used for axis annotation relative to the current setting of cex.
addReferenceLine	
	Logical, should the vertical dashed line be drawn at the reference point.
plotLegend	Logical, should the legend be plotted at the top.
cexLegend	The magnification to be used for legend relative to the current setting of cex.
useMulticore	Logical, should multicore be used. useMulticore = TRUE is supported only on Unix-like platforms.
nrCores	Number of cores to use when useMulticore = TRUE. Default value NULL uses all detected cores.
add	Logical, should the pattern occurrence profiles be added to the existing plot.
	Further arguments to be passed to plot() or lines() methods, such as lty, <i>etc.</i>

Value

The function finds all hits matching the specified patterns in the set of input sequences and plots one average profile per pattern reflecting the occurrence of patterns across sequences.

Author(s)

Vanja Haberle

See Also

getPatternOccurrenceList
plotPatternDensityMap

Examples

```
load(system.file("data", "zebrafishPromoters.RData", package="seqPattern"))
```

TBPpwm

Description

This is a matrix object representing a position-weight matrix (PWM) for TATA-box binding protein motif. The scores in the matrix correspond to log2 ratio between probability of observing base b (rows) at position i (columns) and background probability of base b. This PWM has been derived using position frequency matrix from Jaspar database. It is intended to be used as an input data in running examples from seqPattern package help pages.

Usage

data(TBPpwm)

Format

A matrix object

Value

A matrix object

zebrafishPromoters Zebrafish promoters sequences

Description

This is a DNAStringSet object that contains sequence of 1000 randomly selected promoters active in 24 hpf zebrafish *Danio rerio* embryos. The data is taken from Nepal *et al.* Genome Research 2013, and represents regions flanking 400 bp upstream and 600 bp downstream of the dominant TSS detected by Cap analysis of gene expression (CAGE). It is intended to be used as an input data in running examples from seqPattern package help pages.

Usage

```
data(zebrafishPromoters)
```

Format

A DNAStringSet object

Value

A DNAStringSet object

References

Nepal *et al.* (2013) Dynamic regulation of the transcription initiation landscape at single nucleotide resolution during vertebrate embryogenesis, *Genome Research* **23**(11):1938-1950.

zebrafishPromoters24h Zebrafish promoters genomic coordinates

Description

This is a data.frame object that contains genomic coordinates of 12078 promoters active in zebrafish *Danio rerio* embryos at 24 hours past fertilisation. For each promoter additional information is provided, including position of the dominant (most frequently used) TSS position, number of CAGE tags per million supporting that promoter and the interquantile width of the promoter (width of the central region containing >= 80% of the CAGE tags). The data is taken from Nepal *et al.* Genome Research 2013, and it is intended to be used for running examples from seqPattern package vignette.

Usage

```
data(zebrafishPromoters24h)
```

Format

A data.frame object

Value

A data.frame object

References

Nepal *et al.* (2013) Dynamic regulation of the transcription initiation landscape at single nucleotide resolution during vertebrate embryogenesis, *Genome Research* **23**(11):1938-1950.

Index

* datasets plotPatternOccurrenceAverage,DNAStringSet-method TBPpwm, 19 (plotPatternOccurrenceAverage), zebrafishPromoters, 19 17 PWM, 4–6, 8, 11, 12 zebrafishPromoters24h, 20 PWMscoreStartingAt, 7 * package seqPattern-package, 2 seqPattern, 19, 20 seqPattern (seqPattern-package), 2 bkde2D, 8, 9, 15 seqPattern-package, 2 data.frame, 20 TBPpwm, 19 DNAStringSet, 3, 5, 6, 8, 11, 12, 14, 17, 19 zebrafishPromoters, 19 getPatternOccurrenceList, 3, 5, 16, 18 ${\tt zebrafish Promoters 24h, } 20$ getPatternOccurrenceList,DNAStringSet-method (getPatternOccurrenceList), 3 matchPattern. 3 matchPWM, 5, 8, 11 matrix. 19 motifScanHits, 4, 4, 7, 10, 12 motifScanHits,DNAStringSet,matrix-method (motifScanHits), 4 motifScanScores, 6, 14 motifScanScores, DNAStringSet, matrix-method (motifScanScores), 6 plotMotifDensityMap, 5, 7, 12, 16 plotMotifDensityMap,DNAStringSet,matrix-method (plotMotifDensityMap), 7 plotMotifOccurrenceAverage, 10 plotMotifOccurrenceAverage,DNAStringSet,matrix-method (plotMotifOccurrenceAverage), 10 plotMotifScanScores, 7, 12 plotMotifScanScores, DNAStringSet, matrix-method (plotMotifScanScores), 12 plotPatternDensityMap, 4, 10, 14, 14, 18 plotPatternDensityMap,DNAStringSet-method (plotPatternDensityMap), 14 plotPatternOccurrenceAverage, 17